



JAN 30, 2024

🌐 Immunophenotyping of immune cells by high dimensional flow cytometry

📁 In 1 collection

mariangela.massarocenere^{1,2,3}, Valerio Chiurchiù^{4,5}

¹Department of Experimental Neuroscience, Santa Lucia Foundation IRCCS, Rome, Italy;

²Department of Systems Medicine, University of Rome Tor Vergata, Rome, Italy;

³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, United States;

⁴Laboratory of Resolution of Neuroinflammation, Santa Lucia Foundation IRCCS, Rome, Italy;

⁵Institute of Translational Pharmacology, National Research Council, Rome, Italy

OPEN  ACCESS



n.berretta

ABSTRACT

Protocol to assess cellular phenotypes by multiparametric flow cytometry panels containing markers to identify cell types and activation states

DOI:

dx.doi.org/10.17504/protocols.io.n2bvj3w2plk5/v1

Protocol Citation: mariangela.massarocenere, Valerio Chiurchiù 2024. Immunophenotyping of immune cells by high dimensional flow cytometry. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n2bvj3w2plk5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Jan 25, 2024

Tissue dissociation

- 1 After animal sacrifice (see protocol “intracardiac perfusion”), the dissected and collected cerebral region of interest (SNpc and striatum) is immersed in D-PBS with high glucose buffer at 4°C
- 2 Immediately after, the brain regions are dissociated to single-cell suspension by enzymatic degradation and myelin removal using the adult brain dissociation kit and GentleMACS dissociator (Miltenyi Biotec.)
- 3 The isolated tissue are placed in C-tubes and a mix containing Enzyme P and Buffer Z is added
- 4 Add a second mix containing Enzyme A and Buffer Y, tightly close C-tubes and attach upside down onto the sleeve of the gentleMACS Octo Dissociator with Heaters (Miltenyi Biotec.)
- 5 Run the appropriate program 37C_ABDK_02 for 30 min
- 6 Detach the C-tube and centrifugate to collect sample
- 7 Resuspend and filter to a 70 um MACS SmartStrainer with 10 ml of cold D-PBS with high glucose

- 8 After centrifugation at 300g for 10 minutes, resuspend the sample in the debris removal solution, gently overlay in cold D-PBS with high glucose and centrifuge at 3000g for 10 minutes with full acceleration brake.
- 9 Discard the two phases and resuspend the sample in cold D-PBS with high glucose
- 10 Centrifuge at 1000g for 10 minutes
- 11 Resuspend in the staining buffer (sample is ready for cell count and flow cytometry staining)

Immunophenotyping of brain tissues by flow cytometry

- 12 Stain single suspension cells with PerCP5.5-conjugated anti-CD45 (1:100, Biolegend), APC-Vio770-conjugated anti-CD11b/c (1:100, REAffinity Miltenyi Biotec.), APC-conjugated anti-CD161 (1:100, Biolegend) or CD86 (1:100, Miltenyi Biotec.), PE-conjugated anti-CD45RA (1:100, Biolegend), PE-Cy7-conjugated anti-CD3 (1:100, Miltenyi Biotec.) or anti-MHC-II (1:100, REAffinity Miltenyi Biotec.), and FITC-conjugated anti-granulocytes (1:100 Miltenyi Biotec.) or VioBright-conjugated anti-ACSA-2/O4/CD31 (1:100 REAffinity Miltenyi Biotec.) for 15 minutes at 4°C in the dark
- 13 Wash and resuspend in PBS
- 14 Acquire cells on an appropriate flow cytometer

- 15** Exclude astrocytes, oligodendrocytes, and endothelial cells based on their expression of ACSA-2, O4, and CD31 (Lineage). The remaining cells are identified as CD45+ total leucocytes/resident immune cells
- 16** Inside this gate, microglial cells are identified as CD45lowCD11b/c+ cells and infiltrated myeloid cells as CD45highCD11bhigh, subsequently identified as macrophages according to the expression of F4/80. The remaining of CD45highCD11blow cells are further gated to identify T-lymphocytes as CD3+, B-lymphocytes as CD45RA+, and NK cells as CD161+. The expression of CD68 and MHC-II is further assessed inside the microglial cell population
- 17** For each analysis, acquire at least 0.2×10^6 live by gating on aqua Live/Dead negative cells and then analyze using Flowjo software (Tree Star)