


JAN 15, 2024

🌐 Preparation of mouse tissue homogenates for RT-QulC assay

arpine.sokratian¹, andrew.west¹

¹Duke University



andrew.west

ABSTRACT

This protocol is designed for a standardized and efficient procedure to homogenize mouse tissue, conducive for RT-QulC analysis. The process involves treating samples with PBS mixed with Triton-X100, followed by homogenization using a prob-tip sonicator. Centrifugation is then employed to remove detergent-insoluble fraction, optimizing the tissue for subsequent analysis. The use of Triton-X100, a mild detergent, ensures the release of alpha-synuclein aggregates, enhancing signal specificity without compromising structural integrity. The protocol's significance lies in addressing the critical need for a reliable method in preparing mouse tissue for RT-QulC analysis, with the distinct signal generated in assays validating its efficacy.

PROTOCOL MATERIALS

✕ Low Protein Binding Collection Tubes (1.5 mL) Thermo
Fisher Catalog #90411

Step 1

✕ Roche PhosSTOP EasyPack Merck MilliporeSigma (Sigma-Aldrich) Catalog #4906837001

Step 4

✕ Roche cOmplete Tablets EDTA-free, EasyPack Merck MilliporeSigma (Sigma-Aldrich) Catalog #4693132001

Step 4

✕ Protein LoBind Tubes, 1.5 mL Eppendorf Catalog #0030108116

In 2 steps

BEFORE START INSTRUCTIONS

The tissues need to be perfused with fresh PBS and flash-frozen in liquid nitrogen

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.e6nvwd2w9lmk/v1

Protocol Citation: arpine.sokratian, andrew.west 2024. Preparation of mouse tissue homogenates for RT-QulC assay. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvwd2w9lmk/v1>

MANUSCRIPT CITATION: Reference information: JCI Insight. 2023;8(23):e172192. <https://doi.org/10.1172/jci.insight.172192>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Nov 28, 2023

Last Modified: Jan 15, 2024

PROTOCOL integer ID:
91538

Keywords: ASAPCRN, alpha-synuclein, tissue homogenation, rt-quick

Funders

Acknowledgement:

Aligning Science Across
Parkinson's
Grant ID: ASAP-020527

- 1 Measure the weight of each eppendorf tube using analytical scale and label the tube with assigned number

⊗ Protein LoBind Tubes, 1.5 mL VWR
International Catalog #0030108116

⊗ Low Protein Binding Collection Tubes (1.5 mL) VWR
International Catalog #90411

- 2 Cut a tissue slice (around 100-200 mg) with a disposable blade on a plastic plate lid wrapped in a foil (all on ice), transfer to the Eppendorf protein low-binding tube (with known weight)

Safety information

biohazard; place disposable blades into biohazard needle container; material should be thrown away into biohazard bags previously placed into tube or a 15 mL falcon (secondary container)

- 3 Measure the weight of the tissue piece in a tube using analytical scales
Subtract weight of the tube to obtain the final weight of the brain slice (indicate the weight in the table)



- 4 Add 20 vol weight/volume of PBS buffer containing 1%Tx100 and protease/phosphatase inhibitors (

⊗ Roche PhosSTOP EasyPack Merck MilliporeSigma (Sigma-Aldrich) Catalog #4906837001



Roche cOmplete Tablets EDTA-free, EasyPack Merck MilliporeSigma (Sigma-Aldrich) Catalog #4693132001

(for example, 100 mg – 5 mL) - transfer to 15 mL falcon tube;

- 5 Sonicate tissue samples for  00:01:00, 5 sec ON and 15 sec OFF ( On ice); 1m
- Use 0.16 inch microtip for homogenization

Equipment

Fisherbrand Model 120 Sonic Dismembrator

NAME

Sonicator

TYPE

Fisherbrand

BRAND

FB120110

SKU

<https://www.fishersci.com/shop/products/fisher-scientific-model-120-sonic-dismembrator-4/p-3974654> LINK

Equipped with 1/8" probe

SPECIFICATIONS

- 6 Spin down sonicated samples at  10.000 x g,  4 °C,  00:10:00 10m


- 7 Transfer supernatant into a new 15 mL falcon tube;

- 8 Aliquot the samples in lobind eppendorf



Protein LoBind Tubes, 1.5 mL VWR
International Catalog #0030108116

tubes and store at

 -70 °C