



DEC 04, 2022

WORKS FOR ME

1

Whole blood T cell assay for NHPs, containment protocol

Forked from [Immunophenotyping for NHPs, containment protocol](#)

DOI

dx.doi.org/10.17504/protocols.io.bp2l692prlqe/v1Jonathan Audet¹, Courtney Meilleur¹¹Public Health Agency of Canada

Jonathan Audet

Public Health Agency of Canada

COMMENTS 0

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

This is a protocol used to perform T cell assays using whole blood (collected in EDTA tubes) or cells from bronchoalveolar lavage (BAL). We have successfully used this protocol for rhesus macaques, cynomolgus macaques, and African green monkeys (although antibody mix provided was titrated on macaques). It allows the characterization of T cells into Th1 and Th2 (or Tc1 and Tc2) phenotypes. It technically offers 64 different activation profiles (4 cytokines + CD107a + CD69).

The nature and concentration of the stimulant is left unspecified as we have used different ones depending on reagent availability.

This protocol was designed to deal with samples coming from containment labs (> CL2) at our facility. If you are using this protocol at a different facility please ensure that proper testing and approvals are in place. The protocol can be used for experiments completed entirely in CL2, simply go from step 17 directly to 31.

DOI

dx.doi.org/10.17504/protocols.io.bp2l692prlqe/v1

PROTOCOL CITATION

Jonathan Audet, Courtney Meilleur 2022. Whole blood T cell assay for NHPs, containment protocol.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.bp2l692prlqe/v1>

FORK NOTE

FORK FROM

[Forked from Immunophenotyping for NHPs, containment protocol, Jonathan Audet](#)

KEYWORDS

flow cytometry, nonhuman primate, whole blood, broncho-alveolar lavage, T cell assay

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 02, 2022

LAST MODIFIED

Dec 04, 2022






PROTOCOL INTEGER ID

73500

GUIDELINES

For BAL fluid, make sure to set FSC voltage much lower than for whole blood.

MATERIALS TEXT

-  Human TruStain FcX **BioLegend Catalog #422302** Step 8
-  BD FACS Lysing Solution (10X) **BD Biosciences Catalog #349202** Step 7
-  BD Cytofix/Cytoperm **BD Biosciences Catalog #554722** In 2 steps
-  BD Perm/Wash buffer **BD Biosciences Catalog #554723** In 3 steps
-  Ghost Dye Red 780 **Tonbo Biosciences Catalog #13-0865-T500** Step 10

Equipment			
FACS Tube		NAME	
Tube		TYPE	
Falcon		BRAND	
14-959-2A		SKU	

Equipment			
2 ml screw-cap tubes		NAME	
Microtubes		TYPE	
Sarstedt		BRAND	
72.694.006		SKU	

-  PBS **Contributed by users** Step 18

SAFETY WARNINGS

-

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK








The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://dx.doi.org/10.17504/protocols.io.bp2l692prlqe/v1) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not

substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

BEFORE STARTING

The antibody mix described in the methods was tested on Cynomolgus macaque whole blood and BAL fluid. All antibodies should cross-react with rhesus macaques and African green monkeys but the panel might need to be re-titrated.

Stimulation

- 1 Aliquot  100 µL whole blood into a tube or well.
- 2 Add  200 µL cRPMI containing 1.5 X the desired concentration of stimulant.
(Stimulant(s) and concentration depend on the experiment and can include anti-CD28 and anti-CD49d)
- 3 Mix and incubate at  37 °C 5% CO₂  Overnight .
- 4 Add monensin and brefeldin A, assuming 1e6 cells.
Anti-CD107a BV421 can be added here as well ( 5 µL)
- 5 Mix and incubate at  37 °C 5% CO₂  06:00:00 .

6h

Preparations

1h 10m

- 6 Prepare the staining mix. (for 1 sample:)

A	B	C	D	E
---	---	---	---	---

A	B	C	D	E
Supplier	Antibody	Clone	Channel	Volume per test
BD Biosciences	CD3	SP34-2	Alexa Fluor 700	5
BD Biosciences	CD8	RPA-T8	BV786	5
BD Biosciences	CD45RA	5H9	PE-CF594	2.5
BD Biosciences	CD4	L200	PerCP-Cy5.5	20
BD Biosciences	CCR7	G043H7	BV711	5
BD Biosciences	CD69	FN50	BV650	5
	Total			42.5
	BD Brilliant Stain			57.5



- 7 Prepare the 1X FACS Lysing solution by diluting the 10X stock with Milli-Q water. You will need 2 ml of 1X solution for each sample.

 BD FACS Lysing Solution (10X) **BD Biosciences Catalog #349202**


Surface Staining


30m

- 8 Put 5 µl of TruStain FcX in each tube/well.  Human TruStain FcX **BioLegend Catalog #422302**

- 9 Incubate 10 min at Room Temperature (RT).  00:10:00  Room temperature

10m

- 10 Add 5 µl of the  Ghost Dye Red 780 **Tonbo Biosciences Catalog #13-0865-T500**
If staining BAL: Add 5 µl of 1:100 diluted Ghost Dye.

- 11 Incubate 20 min at RT in the dark.  Room temperature

20m

🕒 00:20:00

12 Add stain mix. Mix. 🧪 100 µL

13 Incubate 20 min at RT in the dark. 🌡 Room temperature

🕒 00:20:00

20m

RBC Lysis

14 Add 2 ml of 1X FACS Lysing solution. Vortex immediately, but gently. 🧪 2 mL

15 Incubate no more than 12 min at RT in the dark. 🕒 00:10:00 No more than 12 minutes

10m

16 Spin at 300 x g for 5 min. 🌀 300 x g, 20°C, 00:05:00

17 Decant supernatant.

18 Add 2 ml of PBS. Vortex. 🧪 2 mL

🔗 PBS Contributed by users

19 Spin at 300 x g for 5 min. 🌀 300 x g, 20°C, 00:05:00

20 Decant supernatant.

Sample Inactivation


21 Resuspend in Cytofix/Cytoperm.  BD Cytofix/Cytoperm **BD Biosciences Catalog #554722**

21.1 100 µl per 5×10^5 cells.

21.2 Use at least 400 µl for easy decanting.

22 Incubate at least 30 min at RT in the dark.  00:30:00  Room temperature In the dark.

30m

23 Spin at 500 x g for 8-10 min.  500 x g, 20°C, 00:08:00

24 On a clean bench, decant supernatant.

25 Use same volume of Cytofix/Cytoperm as before to resuspend the cells.

 BD Cytofix/Cytoperm **BD Biosciences Catalog #554722**

26 Transfer in a 2 ml screwcap tube. Shake the tube to cover all surfaces with Cytofix/Cytoperm.

Equipment

2 ml screw-cap tubes

NAME

Microtubes

TYPE

Sarstedt

BRAND

72.694.006

SKU

27 Transfer tubes from containment space to CL2 space according to the facility's approved protocols/SOPs.


28 Tubes can be opened in CL2 (in a BSC) no less than 30 min after the resuspension (step 17).
(Samples are generally processed the next day; keep at 4 C overnight, in the dark)


30m

Intracellular staining

29 Give tubes a quick spin in a tabletop centrifuge.


30 Ensure all tubes have at least a few hundred microliters of Cytofix/Cytoperm.

31 Transfer the samples into 1 ml of BD  BD Perm/Wash buffer **BD Biosciences Catalog #554723** in FACS tubes.

32 Spin at 500 x g for 8-10 min.  500 x g, 20°C, 00:08:00

33 Decant supernatant.

34 Repeat steps 31-33.



35 Prepare the staining mix. (for 1 sample:)
(uses  BD Perm/Wash buffer **BD Biosciences Catalog #554723**)

A	B	C	D	E
Supplier	Antibody	Clone	Channel	Volume per test
BioLegend	IFNgamma	B27	APC	2.5
BioLegend	TNF	MAb11	PE-Cy7	5
BioLegend	IL-2	MQ1-17H12	Alexa Fluor 488	0.625
BioLegend	IL-4	8D4-8	PE	2.5
	Total			10.625
	PermWash			89.375

36 Resuspend in 100 µl of intracellular staining mix.

37 Incubate  4 °C  00:30:00

30m

38 Add  1 mL  BD Perm/Wash buffer **BD Biosciences Catalog #554723** .

39  500 x g, 4°C, 00:08:00

8m

40 Decant supernatant and blot tubes/plate.

41 Repeat steps 38-40.

42 Resuspend in PBS 1% formaldehyde.

43 Run on flow cytometer.

44 Gating strategy for whole blood: