

OPEN ACCESS



Protocol Citation: Annan SI Cook, Xuefeng Ren 2023. ADP Glo Protocol. protocols.io https://protocols.io/view/adp-

glo-protocol-c2n6ydhe

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: Sep 25, 2023

Last Modified: Sep 30,

2023

PROTOCOL integer ID:

88510

ADP Glo Protocol

Annan SI Cook^{1,2}, Xuefeng Ren¹

¹University of California, Berkeley; ²Aligning Sciences Across Parkinson's



Annan SI Cook

ABSTRACT

Protocol for ADP Glo (Promega) assay to measure the catalytic activity of the lipid kinase VPS34 on small unilammelar vesicles.

ATTACHMENTS

852-2202.pdf

Materials

Small unilamellar vesicles (SUVs) composed of

A
40% dioleoyl phosphatidylcholine (DOPC)
20% dioleoyl phosphatidylethanolamine
20% dioleoyl phosphatidylserine
20% liver phosphatidylinositol

All lipids from Avanti Polar Lipids.

- HEPES pH 7.5
- NaCl
- MgCl₂
- MnCl₂
- TCEP
- ATP
- PI3KC3-C1
- ADP-Glo ATP Depletion Reagent (Promega)

Equipment

- Glass tube
- N2 stream source
- Vacuum dessicator
- Benchtop centrifuge
- Sonicator with a small probe tip
- Water bath at 37°C
- Liquid nitrogen

Preparation of Small Unilamellar Vesicles (SUVs)

15m

- Dehydrate 40.4 mg of the lipid mixture (40% DOPC, 20% DOPE, 20% DOPS, 20% liver phosphatidylinositol, Avanti Polar Lipids) in a glass tube using a gentle N2 stream.
- 2 Allow the lipids to dry Overnight in a vacuum desiccator to remove excess chloroform.



Rehydrate the dried lipids by adding [M] 25 millimolar (mM) HEPES (pH 7.5) and



[м] 200 millimolar (mM) NaCl.

4 Vortex the mixture to resolubilize the lipids.



- Perform 10 cycles of freeze-thaw by immersing the lipids in liquid nitrogen for a few seconds and then in a \$\circ\$ 37 °C water bath.

15m

Preparation of Buffers

7 Prepare a kinase dilution buffer.

Kinase dilution buffer

A	В
HEPES pH 7.5	25 mM
NaCl	200 mM
MgCl ₂	10 mM
MnCl ₂	1 mM
TCEP	2 mM

8 Prepare a 2x lipid dilution buffer.

2x lipid dilution buffer

Α	В
HEPES pH 7.5	25 mM
NaCl	200 mM

A	В
MgCl ₂	20 mM
MnCl ₂	2 mM
TCEP	4 mM

ADP-Glo assay

1h 40m

9 In a reaction tube, combine the following components:

A	В
50 nM PI3KC3-C1	9 μL
0.1 mg/mL SUVs in the 2x lipid dilution buffer	9 μL
500 μM ATP in the same buffer	2 μL

10 Mix the components gently.



Allow the reaction to proceed for 01:00:00 at 8 Room temperature

1h





Incubate the reaction tubes with the ADP-Glo ATP depletion reagent for 00:40:00 at

40m



! Room temperature

Measure the luminescent output for each condition tested using a luminometer or a plate reader.