



Apr 15, 2021

MMH1 Solution

Erçağ Pinçe¹¹AMOLF

1

Works for me

dx.doi.org/10.17504/protocols.io.mknc4ve

Erçağ Pinçe

ABSTRACT

Minimal medium H1 protocol for bacterial competition experiments

DOI

dx.doi.org/10.17504/protocols.io.mknc4ve

PROTOCOL CITATION

Erçağ Pinçe 2021. MMH1 Solution. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.mknc4ve>



LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jan 12, 2018

LAST MODIFIED

Apr 15, 2021

PROTOCOL INTEGER ID

9582

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Prepare stock solutions of the salts

1 Prepare Potassium Phosphate:

Name of the chemical	Formula	Brand	Molecular Weight(gr/mol)	Final Concentration	Stock	Factor	Stock in 100mL in gr
Potassium phosphate(Monobasic)	KH ₂ PO ₄	Sigma-Aldrich	139.06	50 mM	1M	20x	9.5
Potassium Phosphate(Dibasic)	K ₂ HPO ₄	Sigma-Aldrich	174.18	50 mM	1M	20x	5.25

- Dissolve the reagents in 80 mL of ddH₂O
- Adjust pH to 7.0
- Make the volume up to 100 mL with ddH₂O
- Filter sterilize and store at room temperature.

Prepare stock solutions of other salts:

Name of the chemical	Formula	Brand	Molecular Weight(gr/mol)	Final Concentration	Stock Con.	Factor	Stock in 200 mL in gr
Ammonium Sulfate	(NH ₄) ₂ SO ₄	Sigma-Aldrich	132.14	7.6mM	760mM	100x	20.09
Magnesium Sulfate heptahydrate	MgSO ₄	Sigma-Aldrich	246.47	0.5mM	50mM	100x	2.46
Sodium Chloride	NaCl	Sigma-Aldrich	58.44	67mM	1.34M	20x	15.66
Iron Sulfate (Cover with foil)	Fe ₂ (SO ₄) ₃	Sigma-Aldrich	399.88	1.25μM	1.25mM	1000x	0.1

- Mix each compound in 200 mL in separate beakers and stir for ~20mins at room temperature.
- Filter sterilize each solution and store at room temperature.
- **Important Note:** You can either use the stock solution of NaCl [1.34M] or you can directly add it to the final MMH1 mix in its solid form.

The protocol for preparing 1M potassium phosphate was taken from Cold Spring Harbor Protocols

http://cshprotocols.cshlp.org/content/2017/7/pdb.rec090290.full?text_only=true

Prepare MMH1 mix

- MMH1:** 1) Add the salts and Thiamine
 - 435 mL ddH₂O
 - 5 mL Ammonium Sulfate [100x]
 - 5 mL Magnesium Sulfate [100x]
 - 5 mL Iron(III) Sulfate [100x]
 - 25 mL Sodium Chloride [20x]
 - 25 mL Potassium Phosphate [20x]
 - 50 mg Thiamine Hydrochloride (Sigma-Aldrich)

MMH1[2x]:

- 370 mL ddH₂O
- 10 mL Ammonium Sulfate [100x]
- 10 mL Magnesium Sulfate [100x]
- 10 mL Iron(III) Sulfate [100x]
- 50 mL Sodium Chloride [20x]
- 50 mL Potassium Phosphate [20x]
- 100 mg Thiamine Hydrochloride (Sigma-Aldrich)

2-) Stir mixes at room temperature for ~30mins.

3-) Measure pH and adjust it to 7.0 if necessary.

4-) Filter sterilize and store in the cold room(+4C)