



Aug 13, 2020

Universal sandwich ELISA for investigating the binding of Protein L to avian immunoglobulins using anti-IgY-peroxidase as conjugate.

Angel A Justiz-Vaillant¹, Monica F. Smikle²¹University of the West Indies St. Augustine; ²University of the West Indies. Mona Campus**1** Works for me dx.doi.org/10.17504/protocols.io.bjqykmxw[University of the West Indies](#) angel.vaillant@sta.uwi.eduAngel Justiz-Vaillant
University of the West Indies St. Augustine

DOI

dx.doi.org/10.17504/protocols.io.bjqykmxw

PROTOCOL CITATION

Angel A Justiz-Vaillant, Monica F. Smikle 2020. Universal sandwich ELISA for investigating the binding of Protein L to avian immunoglobulins using anti-IgY-peroxidase as conjugate.. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bjqykmxw>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 13, 2020

LAST MODIFIED

Aug 13, 2020

PROTOCOL INTEGER ID

40440

MATERIALS

NAME	CATALOG #	VENDOR
Anti-Chicken IgY, HRP Conjugate, 300ul	G1351	Promega
Nunc™ 96-Well Polystyrene Round Bottom Microwell Plates, V 96 well plate, Non-Treated, clear, without lid, Sterile	260210	Thermo Fisher
Protein-L from P. Magnus		

- 1 This ELISA is used to study the interaction of recombinant protein-L (SpL) with diverse avian immunoglobulins.
- 2 The 96 well microtitre plate is coated overnight at 4°C with 2 µg/µl per well of SpL in carbonate-bicarbonate buffer pH 9.6.

- 3 Then plate is treated with bovine serum albumin solution and washed 4X with PBS-Tween.
- 4 50 µl of avian egg yolk or egg white (1 mg/ml) is added and incubated for 1.30h at room temperature and the microplate is then rewashed 4X with PBS-Tween.
- 5 Then 50 µl of peroxidase-labeled-anti-IgY conjugate diluted 1:15000 in PBS-non-fat milk is added to each well and incubated for 1.30h at RT. The plate is washed 4X with PBS-Tween.
- 6 Pipette 50 µl of 3,3',5,5' - tetramethylbenzidine (TMB; Sigma-Aldrich) to each well.
- 7 The reaction is stopped with 50 µl of 3M H₂SO₄ solution.
- 8 The plate is visually assessed for the development of colour and read in a microplate reader at 450 nm.
- 9 A cut-off point should be calculated as the mean of the optical density of negative controls x 3. The higher the OD value the higher will be the affinity of SpL to avian immunoglobulins.