



## HBSS/HEPES for ROS measurement in PMN

Marco Cosentino<sup>1</sup>, Elisa Storelli<sup>1</sup>, Alessandra Luini<sup>1</sup>, Massimiliano LM Legnaro<sup>1</sup>, Emanuela Rasini<sup>1</sup>, Marco Ferrari<sup>1</sup>, Franca Marino<sup>1</sup>

<sup>1</sup>Center for Research in Medical Pharmacology, University of Insubria (Varese, Italy)

**1** Works for me [dx.doi.org/10.17504/protocols.io.bhdcj22w](https://dx.doi.org/10.17504/protocols.io.bhdcj22w)

Mattia Di Rocco

DOI

[dx.doi.org/10.17504/protocols.io.bhdcj22w](https://dx.doi.org/10.17504/protocols.io.bhdcj22w)

DOCUMENT CITATION

Marco Cosentino, Elisa Storelli, Alessandra Luini, Massimiliano LM Legnaro, Emanuela Rasini, Marco Ferrari, Franca Marino 2020. HBSS/HEPES for ROS measurement in PMN. **protocols.io**  
[dx.doi.org/10.17504/protocols.io.bhdcj22w](https://dx.doi.org/10.17504/protocols.io.bhdcj22w)

LICENSE

This is an open access document distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 10, 2020

LAST MODIFIED

Jul 07, 2020

DOCUMENT INTEGER ID

38020

### HBSS/HEPES x ROS composition (g/l):

**NaCl:** 4.2300 g

**KCl:** 0.1864 g

**MgSO<sub>4</sub>:** 0.1232 g

**CaCl<sub>2</sub>:** 0.735 g

**Glucose:** 0.900 g

**Hepes:** 1.300 g

To prepare 0.5 L of HBSS/HEPES, dissolve the reagents listed above in 450 mL of ultrapure H<sub>2</sub>O.

Adjust the pH to 7.4 with HCl or NaOH, and then add H<sub>2</sub>O to 0.5 L.

Storage: 4°C Fridge 1- (Room TS08)

**NaCl code:** S9625, Sigma

**KCl code:** P9541, Sigma

**MgSO<sub>4</sub> code:** 1.05886, Sigma

**CaCl<sub>2</sub> code:** 1.02382 Sigma

**Glucose code:** 1.08337, Sigma

**Hepes code:** H7006 Sigma