



APR 24, 2023

## IHC\_cFOS+Parvalbumin

Bilge Buyukdemirtas<sup>1</sup>

<sup>1</sup>Wesleyan University



Bilge Buyukdemirtas

### ABSTRACT

Immunohistochemistry protocol for c-FOS and Parvalbumin co-staining of mouse brain tissue slices (30 microns thick).

### GUIDELINES

All shaker incubations are done at 200 rpm speed.

### OPEN ACCESS

**DOI:**  
[dx.doi.org/10.17504/protocols.io.81wgbyb9yvpk/v1](https://dx.doi.org/10.17504/protocols.io.81wgbyb9yvpk/v1)

**Protocol Citation:** Bilge Buyukdemirtas 2023. IHC\_cFOS+Parvalbumin. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.81wgbyb9yvpk/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
 We use this protocol and it's working


**Created:** Apr 24, 2023

**Last Modified:** Apr 24, 2023


**PROTOCOL integer ID:**  
 80981

## Staining

- 1** Antigen Retrieval  
 Remove PBS from wells. Add 500ul 0.3% citrate buffer into each well.

- 1.1 Incubate at RT on shaker for 30 min  
Incubate at  65 °C for 45 min  
Cool down on the bench for 20 min

## 2 3X PBS Washes

- 2.1 Remove citrate buffer from wells and replace with 500ul 1X PBS (or use a dropper pipette and fill the well halfway)
- 2.2 Incubate at RT on shaker for 15 min  
 remove PBS and replace with new 1X PBS, repeat two more times

## 3 Block Non-Specific Signal Make enough: 10% NGS (Vector Labs, #S-1000-20) in 0.3% triton-X PBS for all samples (500ul per well)

- 3.1 Remove PBS and replace with blocking reagent 500ul per well.
- 3.2 Incubate on shaker RT for 60min.

## 4 Primary Antibody Make enough: mouse anti-parv (1:1000, Sigma-Aldrich, #P3088) and rabbi anti-cfos (1:1000, cell Signaling, #2250S) in 1% NGS (Vector Labs, #S-1000-20) in 0.3% triton-X PBS


**4.1** Replace the blocking reagent with primary antibody mix, 500 ul per well, and incubate at

 4 °C

 Overnight

**5** DAY 2  
3X PBS Washes

**5.1** Remove primary Ab mix from wells and replace with 500ul 1X PBS (or use a dropper pipette and fill the well halfway)

**5.2** Incubate at RT on shaker for 15 min  repeat two more times

**6** Secondary Antibody  
Make enough: goat anti-mouse AF-555 (1:500, Invitrogen, #A-21422) and biotin goat anti-rabbit (1:1000, Vector Labs, #PK-6101) in 1% NGS (Vector Labs, #S-1000-20) in 0.3% triton-X PBS

**6.1** Remove PBS from wells and replace with 500ul Secondary antibody mix. Cover with aluminum foil. Note: secondary antibodies are light sensitive, so from this point on the plate MUST be covered from light as much as possible.

Incubate at RT on shaker for 120 min.

**7** 3X PBS washes as described in Step 5 and 5.1 with the plate covered 

**8** Tertiary Antibody (for cFOS only)  
make enough: AF-488 streptavidin conjugated (1:1000, Thermo Fisher, #S11223) in 1% NGS (Vector Labs, #S-1000-20) in 0.3% triton-X PBS

**8.1** Replace PBS in wells with 500uL of tertiary antibody solution, cover the plate with aluminum foil, and incubate at RT on shaker for 120 min.

**9** 3X PBS washes as described in step 5 and 5.1 with the plated covered 

## Mounting

**10** Mounting with medium with and without DAPI  
Step 10 includes a Step case.  
**If the mounting medium has DAPI included**  
**If the mounting medium doesn't have DAPI included**