



2 ▼

Mar 23, 2022

DNA metabarcoding protocol for siphonophore gut contents V.2

Alejandro Damian Serrano¹¹Yale University

1

dx.doi.org/10.17504/protocols.io.5qpvo57o7l4o/v2

Dunn Lab

Alejandro Damian Serrano
Yale University

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Protocol for the SiphWeb DNA metabarcoding of siphonophore gut content.

Starting with DNA extractions from pooled frozen gastrozooids, this protocol takes the user through the steps of PCR amplification, PCR cleanup and pooling of amplicons from six complementary 18S barcodes that can be submitted into an Illumina MiSeq Lane.

DOI

dx.doi.org/10.17504/protocols.io.5qpvo57o7l4o/v2https://docs.google.com/spreadsheets/d/1x71z9YLqxo9XszNcjAPB_p0cnHHYyMTEHNL4wb84j4/edit#gid=656913706

Alejandro Damian Serrano 2022. DNA metabarcoding protocol for siphonophore gut contents.

protocols.io<https://dx.doi.org/10.17504/protocols.io.5qpvo57o7l4o/v2>

Alejandro Damian Serrano



Corrected the barcode names and start-end positions on the 18S gene to reflect the hypervariable regions as reference. Corrected mistyped annealing time values in the PCR thermocycler programs.

protocol ,

Mar 23, 2022

Mar 23, 2022

59815

MATERIALS

GoTaq(R) Long PCR Master Mix, 10

Reactions Promega Catalog #M4020

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

DNA Extraction

4h

4h

- 1 Mainly identical to:
<http://www.bea.ki.se/documents/EN-DNeasy%20handbook.pdf>
With the following modifications:

-Digestion at δ 56 °C 1-2h

-Elution using 2 rounds of incubation and centrifuge with \square 50 μ L of AE buffer , incubating \odot 00:10:00 at δ 56 °C , to a total of \square 100 μ L

-Label \square 2 mL Eppendorf tubes to hold the extracted DNA. Include the extraction number, the siphonophore species, "GCDNA", the date of extraction, and your initials.

-Use Nanodrop to assess DNA yield (ng/ μ L) and 260:280 ratio.

NanoDrop™ 3300 Fluorospectrometer
Fluorospectrometer

NanoDrop™ 3300 nd-3300 [↗](#)



Store eluted DNA in δ -20 °C freezer.

PCR

2h 52m

- 2 Based off Promega GoTaq protocol.

(Plan for a 96-well final plate [90 wells used for 13 templates, a positive, and a negative control, across 6 primer pairs (detailed in T1)], can be re-scaled as needed).

See reference planning spreadsheet here:

https://docs.google.com/spreadsheets/d/1x71z9YLqxo9XszNcjAPB_p0cnHHYyMTEHNL4wb84j4/edit#gid=656913706

Only select templates that have a DNA yield > 10ng/ μ L and a 260:280 > 1.7

3

A	B	C	D	E	F	G	H	I	J	K	L	M	N
Original Name	Type (F/R)	Barcode Region	Primer sequence	Length (bp)	TM © min	TM © max	%GC	Hairpin TMC	Self Dimer TMC	Pair Dimer TMC	Insert size (bp)	Start position	End position
152F	F	V5-V7S	TGACGGAAGGGCACCACCAG	20	62.7	64.6	63.2	63	0.1	None	152	1187	1207
152R	R	V5-V7S	TCCACCAACTAAGAACGGCC	20	60	60	55	None	None	None	152	1319	1339
166F	F	V3	AACGGCTACCACATCCAAGG	20	60	60	55	None	None	None	146	420	440
166R	R	V3	CACCAGACTTGCCCTCCAAT	20	60	60	55	31.9	None	None	146	546	566
272F	F	V5-V7L	AAACGATGCCGACTAGCGAT	20	59.9	59.9	50	44.6	7.9	6.7	272	1067	1087
272R	R	V5-V7L	TCCACCAACTAAGAACGGCC	20	60	60	55	None	None	6.7	272	1319	1339
179F	F	V7	GGCCGTTCTTAGTTGGTGA	20	60	60	55	None	None	11	170	1319	1339
179R	R	V7	TGCGGCCAGAACATCTAAG	20	60.1	60.1	55	None	None	11	170	1469	1489
261F	F	V7p+V8	AACAGGTCTGTGATGCCCTT	20	59.2	59.2	50	44.1	4.2	14	215	1472	1492
261R	R	V7p+V8	TGTGTACAAAGGGCAGGGAC	20	59.9	59.9	55	None	None	14	215	1667	1687
134F	F	V9	CTTTGTACACACCGCCGTC	20	61.6	61.6	60	None	None	None	115	1675	1695
134R	R	V9	CCTTGTTACGACTTTTACTTCTCT	25	58.8	58.8	40	None	None	None	115	1765	1790

T1. Primer sequences and properties. Positions based on the 18S Gene of *Lymnaea diaphana* (GenBank JF909497.1).

4 -Pull all reagents in T2 from the freezer and let them thaw.

1m 45s

-Separate and label two PCR plates,

one for V9 (annealing at **48 °C** for **00:00:45** per cycle)

one for V3, V5-V7S, V5-V7L, V7, and V7p+V8 (all annealing at **54 °C** for **00:01:00** x30 cycles). Divide and label each pair of rows that will contain each primer pair mix.

5

10m

Reagent	Volume (µm)
Template (variable)	2
Forward primer (10µM)	0.5
Reverse primer (10µM)	0.5
MgCl2 (25mM)	2.5
dNTPs (10mM)	1
GoTaq	0.2
Buffer 5X (green)	2.5
BSA	1
H2O	14.8
TOTAL	25

T2. Reagent breakdown goal per well in the plate.

-Make master-masternmix for all primer pairs as in T3. Always vortex each reagent before pipetting, and vortex the mixture at the end.

Reagent	Volume (µm)
MgCl2 (25mM)	350
dNTPs (10mM)	140
GoTaq	28
Buffer 5X	350
BSA	140
H2O	2072
TOTAL	3080

T3. Master-masternmix volumes for 6 primer pair masternmixes downstream.

- 6 -Serve **440 µL** of this master-mastermix into 6 **2 mL** Eppendorf tubes. Lable the tubes with the name of each primer pair. ^{10m}
 -Add **10 µL** of **10 micromolar (µM)** Forward and **10 µL** of **10 micromolar (µM)** Reverse primer into each primer-specific mastermix to end up with T4.
 -Vortex .

A	B
Reagent	Volume (µl)
Forward primer (10µM)	10
Reverse primer (10µM)	10
MgCl2 (25mM)	50
dNTPs (10mM)	20
GoTaq	4
Buffer 5X	50
BSA	20
H2O	296
TOTAL	460

T4. Volumes in primer pair mastermix.

- 7 In a set of 14 PCR tubes (can use two rows in either PCR plate too), serve an excess amount of template to seed each well with. I usually go for **20 µL** . Serve the templates in the same spatial arrangement as they will be seeded in each set of rows for each primer pair. ^{10m}

Include the positive and negative controls. I tend to pick the elution AE buffer stock I used for extractions as negative control, and a previously successful template as positive control.

- 8 Using a single **200 µL** pipette tip on a multi-step pipette set to **23 µL** doses in 4 steps, seed each primer-specific pair of rows. Remember to get rid of the excess (before or after serving the steps, depending on the model of the pipette). Change tips in between primer-pair-specific mastermix. ^{30m}

Using a multichannel 10-100µl, seed the **2 µL** of template in each well. Be careful to not accidentally aspire air at the bottom of the tubes/wells and end up with un-seeded wells. Always visually inspect the tips before serving into the wells. Change tips every time to avoid contamination.

- 9 Seal the plates (**00:00:06** under the plate press) with cellophane cover. ^{1m}
 Vortex plates.
1000 rpm, 00:00:10 , Benchtop manual plate centrifuge and spin it down.

- 10 Insert plates (V9) and (V3, V5-V7S, V5-V7L, V7, V7p+V8) in two separate thermocyclers. Run PCRs following the programs specified in T5A and T5B respectively. ^{1h 45m}

SimpliAmp Thermal Cycler

PCR

Applied Biosystems A24811 [↗](#)

Any standard PCR thermocycler will suffice



A	B	C	D	E	F
Temperature (C)	95	95	48	72	72
Time	2m	30s	45s	1m	5m
Cycles		x30	x30	x30	

T5A. PCR program for barcode V9.

A	B	C	D	E	F
Temperature (C)	95	95	54	72	72
Time	2m	30s	1m	1m	5m
Cycles		x30	x30	x30	

T5B. PCR program for barcodes V3, V5-V7S, V5-V7L, V7, and V7p+V8.

11 Vortex plates.

1m

[🌀1000 rpm, 00:00:10](#) , [Benchtop manual plate centrifuge](#) and spin it down.

Un-seal the plates ([🕒00:00:06](#) under the plate press)

Gel Electrophoresis

1h 25m

12 [📦2 g](#) of Agarose for every [📦100 mL](#) of TAE buffer.

30m

[📦TAE Buffer, 10X,](#)

[1000ml Promega Catalog #V4271](#)

[📦Agarose Thermo](#)

[Fisher Catalog #75510019](#)

For 90 wells, we will need 3 gels with 2 16-well combs each.

For 3 gels, we need to make [📦200 mL](#) of agar in one batch, and [📦100 mL](#) in another batch.

In the double batch, add [📦4 g](#) of agarose and [📦200 mL](#) of TAE buffer. Microwave for [🕒00:03:00](#) . Stir and microwave for another [🕒00:01:00](#) if solids are still not dissolved. When cool, add [📦4 µL](#) of SYBR-Safe.

In the single batch, add [📦2 g](#) of agarose and [📦100 mL](#) of TAE buffer. Microwave for [🕒00:02:30](#) . Stir and microwave for another [🕒00:01:00](#) if solids are still not dissolved. When cool, add [📦2 µL](#) of SYBR-Safe.

[📦SYBR SAFE DNA stain Life](#)

[Technologies Catalog #S33102](#)

Pour gels and let them cool until solid.

Annotate distribution of samples on wells. I usually arrange a row per primer pair ordered by amplicon size: V9:V5-V7S, V3:V7, V7p+V8:V5-V7L in 3 gels, and then the well sequence would go: Ladder, A1-B6 (letters vary with position in the PCR plate), -ve, +ve.

- 13 Set up power source to 125V, 500mA, ⌚ 00:45:00 55m
- Load 📄 6 µL of 100bp ladder on the left side of each row of wells (2 per gel).
- Use adjustable-spacer multichannel pipette to load 📄 5 µL of PCR products onto the wells.
- Well distance is ➦ 6.25 mm , PCR plate distance is ➦ 9 mm .
- Run gels.

- 14 Place gels under UV light, take photograph, print. 2m



You should see a band over each well except the negative control. The bands should land at approximately the height (in bp) of the ladder corresponding to the name (expected amplicon length) of each primer pair.

PCR Cleanup

47m 6s

- 15 Put some nuclease-free water to warm up to ⚡ 56 °C . 10m

🔗 [ExcelaPure™ 96-Well UF PCR Purification Plates \(No Receivers\) Edge](#)

Bio Catalog #36181

Add 📄 55 µL of water into each PCR well (except controls and failed reactions) using a multichannel 20-200µl and filter tips. Mix by pipetting up and down ~5 times. Transfer mix to ExcelaPure purification plate. Change tips every time.

- 16 ⚙️ 5000 rpm, ⌚ 00:07:00 Centrifuge the purification plate with a collecting plate underneath. 7m
- Make sure to balance the **plate centrifuge** with a couple of adequately weighted collecting plates.

- 17 Add 📄 50 µL of warm (⚡ 56 °C) water to each well in the purification plate. Mix by pipetting up and down ~15 10m
- times with the multichannel pipette and filter tips. Change tips every time.

- 18 Incubate at ⚡ 56 °C for ⌚ 00:10:00 . 10m
- Prepare a final clean product PCR plate: Label and annotate adequately to preserve sample location. I recommend drawing a line between each pair of rows, segregating primer pairs for barcodes V3 to V9.

- 19 Mix the warm eluted clean products in the purification plates by pipetting up and down ~15 times with the 10m
- multichannel pipette and filter tips.
- Transfer clean products to the final plate.
- Change tips every time.

- 20 Seal the plates (⌚ 00:00:06 under the hot plate press) with cellophane cover. Store in ⚡ -20 °C freezer. 6s

Amplicon isomolar pooling

2h 50m

- 21 The goal here is to obtain an equal concentration of each barcode for each sample, so they can be pooled with equal 40m
- representation for MiSeq Illumina sequencing.

First step is to quantify the DNA content in each cleaned PCR product well using Qubit HS.

Qubit 2.0 Fluorometer instrument Q33226
with Qubit RNA HS Assays

If you have the pre-mixed mastermix, just run the standards at **190 µL** of mastermix and **10 µL** of standard.
Vortex mix and incubate for 1min.

Prepare the assay tubes with **198 µL** MM and **2 µL** of template.

Vortex mix and incubate for 1min.

Pick the ng/µl units.

22 -Calculate amount of **template** to seed in each pooling well.

-Pool target concentrations are equal to the lowest concentration across barcodes for a sample excluding those that are <2ng/µl, which is the minimum required by YCGA for amplicon sequencing.

A	B	C	D	E	F
V9	V5-V7S	V3	V7	V7p+V8	V5-V7L
0.5	0.5	0.6	0.6	0.9	1

T6. Molar correction factors for each barcode based on the molecular weight of each expected amplicon length.

Template volume to add = (**50 µL** / (Qubit ng/µl value of original PCR product/Pool target concentration)) * Molar_correction_factor

~ See T6 for Molar correction factors ~

23 -Calculate amount of **water** to add in each pooling well.

Water volume to add = (**50 µL** * Molar_correction_factor) - **Template volume to add**

24 Prepare a pooling PCR plate, labeled.

2h

Add the water and template quantities for each well. Change the pipette tip every time when handling PCR products.

25 Label Eppendorf **2 mL** tubes to hold to final amplicon pools. Include the extraction number, the genus of the siphonophore, "GC DNA pooled", and your initials.

10m

Using a multichannel pipette of 20-200µl, proceed to pool the amplicons by sample. If the volumes do not fit in a row on the pooling PCR plate, you can make 2 semi-pools, one combining V9:V5-V7S:V3, and one combining V7:V7p+V8:V5-V7L.



Be very careful on this step, since mistakes during pooling will result in having to re-do the whole protocol for the samples affected.

Using a single-channel 20-200µl pipette, transfer the pooled amplicons to the Eppendorf tubes. Vortex and close caps.

Store the amplicon pools in a **-20 °C** freezer.