

Jun 25, 2024

LYSOSOMAL ISOLATION PROTOCOL

DOI

dx.doi.org/10.17504/protocols.io.kqdg32kdqv25/v1

Scott Vermilyea¹

¹Univ. of Minnesota, Team Lee

ASAP Collaborative Rese...

Team Lee



Jane Balster

ASAP - Team Lee

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.kqdg32kdqv25/v1

Protocol Citation: Scott Vermilyea 2024. LYSOSOMAL ISOLATION PROTOCOL. protocols.io

<https://dx.doi.org/10.17504/protocols.io.kqdg32kdqv25/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 07, 2024

Last Modified: June 25, 2024

Protocol Integer ID: 102392

Keywords: Crude lysosomal fraction, Homogenisation , Precipitation, Lysis buffer

Funders Acknowledgement:

ASAP

Grant ID: 000592

Abstract

This protocol details the isolation of lysosomes.

Materials

Homogenisation Buffer:

A	B
Sucrose	250 mM
EDTA	2 mM
Magnesium chloride	1.5 mM
Potassium Chloride	10 mM
HEPES	20 mM

TNE lysis buffer:


A	B
Tris	50 mM
NaCl	150 mM
EDTA, pH 7.4	5 mM
SDS	1%
NP-40	0.5%
DOC	0.5%
protease/phosphatase inhibitors	



Subcellular Lysosome Isolation Protocol


30m

1 In Vitro Harvest

- 1.1 Following designated treatment, wash cells three times with DMEM followed by the application of  0.5 mL of homogenization buffer supplemented with proteinase and phosphatase buffer.

Homogenisation Buffer:

A	B
Sucrose	250 mM
EDTA	2 mM
Magnesium chloride	1.5 mM
Potassium Chloride	10 mM
HEPES	20 mM

- 1.2 Gently detach cells with cell scraper.
- 1.3 Homogenize using Teflon homogenizer (12 strokes).
- 1.4 Separate  50 μ L of total homogenate (TH) and lyse with *TNE lysis buffer.


TNE lysis buffer:

A	B
Tris	50 mM
NaCl	150 mM
EDTA, pH 7.4	5 mM
SDS	1%
NP-40	0.5%
DOC	0.5%
protease/phosphatase inhibitors	

1.5 Centrifuge remaining volume of TH at  1000 x g, 4°C, 00:10:00 , collect the supernatant.

10m



1.6 Further centrifuge the supernatant for  20000 x g, 4°C, 00:20:00 to collect the precipitate as crude lysosomal fraction (CLF). Lyse the CLF with *TNE lysis buffer, and prepare lysates for western blot.

20m



TNE lysis buffer:

A	B
Tris	50 mM
NaCl	150 mM
EDTA, pH 7.4	5 mM
SDS	1%
NP-40	0.5%
DOC	0.5%
protease/phosphatase inhibitors	