

Rab8a expression
and purification

May 23, 2022

Rab8a expression and purification

David M. Snead^{1,2}, Mariusz Matyszewski¹¹Department of Cellular and Molecular Medicine, University of California, San Diego, La Jolla, CA 92093;²Department of Biochemistry and Molecular Biology, Johns Hopkins University Bloomberg School of Public Health, Baltimore, MD, 21205

1

dx.doi.org/10.17504/protocols.io.6qpvr63mzvmk/v1

Mariusz Matyszewski

Recombinant Rab8a expression and purification protocol as used by the Leschziner and Reck-Peterson Labs.

Original protocol by David Snead. Adapted for protocols.io by Mariusz Matyszewski.

Current version as used in Snead, Matyszewski, Dickey et al. 2022.

DOI

dx.doi.org/10.17504/protocols.io.6qpvr63mzvmk/v1

David M. Snead, Mariusz Matyszewski 2022. Rab8a expression and purification.
protocols.io

<https://dx.doi.org/10.17504/protocols.io.6qpvr63mzvmk/v1>

ASAP

Grant ID: ASAP-000519

MJFF

Grant ID: 18321

Rab8a, protein purification, protein expression

_____ protocol ,

Mar 03, 2022

May 23, 2022

59034

Materials:

Ni-NTA beads

IgG beads

Purified TEV protease

Equipment:


FPLC

Ultracentrifuge

Sonicator


Buffers:

Lysis Buffer (make  150 mL):

- [M]50 millimolar (mM) HEPES pH 7.4
- [M]200 millimolar (mM) NaCl
- [M]10 % volume Glycerol
- [M]5 millimolar (mM) MgCl₂
- [M]2 millimolar (mM) DTT
- [M]0.5 millimolar (mM) Pefabloc
- 1 cOmplete Protease Inhibitor Cocktail tablet per  50 mL (3 for suggested volume in this prep)

Wash Buffer (make  750 mL):

- [M]50 millimolar (mM) HEPES pH 7.4
- [M]150 millimolar (mM) NaCl
- [M]5 % volume Glycerol
- [M]5 millimolar (mM) MgCl₂
- [M]1 millimolar (mM) DTT

Elution Buffer (make  100 mL):

- [M]50 millimolar (mM) HEPES pH 7.4
- [M]150 millimolar (mM) NaCl
- [M]5 % volume Glycerol
- [M]5 millimolar (mM) MgCl₂
- [M]1 millimolar (mM) DTT
- [M]300 millimolar (mM) Imidazole pH 8.0

TEV Buffer (make **200 mL** ; will also need a little bit of extra TEV buffer with **10 % volume Glycerol** for storage):

- **50 millimolar (mM) HEPES pH 7.4**
- **200 millimolar (mM) NaCl**
- **5 % volume Glycerol**
- **5 millimolar (mM) MgCl₂**
- **1 millimolar (mM) DTT**

S200 Buffer (make **250 mL**):

- **50 millimolar (mM) HEPES pH 7.4**
- **200 millimolar (mM) NaCl**
- **1 % volume Glycerol**
- **5 millimolar (mM) MgCl₂**
- **1 millimolar (mM) DTT**

Expression

1h 6m

- 1 Transform n-terminally His6-ZZ tagged Rab8a (pET28a backbone) into BL21(DE3) E. coli cells.

BL21-CodonPlus (DE3)-RIPL have also been used successfully. Please adjust the antibiotics used based on the cells used.
pET28a backbone comes with Kanamycin resistance

- 2 Make fresh LB media. This protocol assumes **4 L (2L per flask)** are made. Scale accordingly.
- 3 Grow an overnight culture in LB media with **50 µg/mL Kanamycin**
Make **50 mL** . This is enough for the main 4L growth.
- 4 Add the overnight culture into main LB flasks with antibiotic present (**50 µg/mL Kanamycin**). Use 100x dilution (**10 mL per 1L**).

- 5 Let it grow in a shaker at 🧪 **200 rpm, 37°C** until OD600 of 0.4-0.5 is reached.
- 6 Chill the cells and set the shaker to 🌡 **18 °C** . We chill the cells in the cold room at 🌡 **4 °C** ^{1h} for ⌚ **01:00:00**
- 7 Add IPTG at [M] **0.5 millimolar (mM) final concentration** into the chilled culture.

Higher IPTG concentrations can be used if needed.

- 8 Let it grow in a shaker at 🧪 **200 rpm, 18°C, 18:00:00**
- 9 Harvest the cells. We spun them down using a JLA 9.1000 rotor at ^{6m} ⚙ **6000 rpm, 4°C, 00:06:00** in 1L Beckman bottles. We then resuspended the cells in 🧴 **15-30 mL LB** and transferred them to 50 mL Falcon tubes (1 tube per 2L harvested). Flash frozen with liquid nitrogen and stored in 🌡 **-80 °C** freezer.

Purification: Day 1 2h

- 10 Resuspend the cell pellet in **lysis buffer**. Add about 🧴 **40 mg Lysozyme** and incubate ^{30m} ⌚ **00:30:00** at 🌡 **On ice**
- 11 Lyse the cells by sonication 🌡 **On ice** .
- 12 Ultracentrifugation in Ti-70 rotor ⚙ **30000 rpm, 4°C, 00:30:00** ^{30m}

Done using ultracentrifugation in our lab, but regular centrifugation should work too.
Adjust accordingly.

12.1 While centrifuging, prepare **2 mL Ni-NTA beads** by equilibrating with **lysis buffer**.

13 Once centrifugation is done, add the beads to the supernatant. Let it nutate for **01:00:00**^{1h} **4 °C**.

13.1 Prepare **Wash and Elution buffers** while waiting.

14 Apply beads to a gravity column and wash with **250 mL Wash buffer**






15 Elute with **40 mL Elution buffer**. Best to resuspend in **20 mL**, elute, resuspend again, and finish elution.

16 Equilibrate IgG beads with **Wash Buffer**.


17 Dilute the eluted protein to **90 mL Wash Buffer**. Split in 2 (**45 mL each**). Add **2 mL IgG beads** to each.

18 Nutate at **4 °C** for 2-3 hours.


19 Apply to gravity column and wash with **250 mL Wash buffer**

- 20 Make **TEV buffer** and transfer beads with  **10 mL TEV Buffer** into 15 mL tubes.
- 21 Add  **400 µL**  **1.3 mg/mL** TEV protease and incubate  **Overnight** at  **4 °C**

Purification: Day 2

- 22 Equilibrate Ni-NTA beads with **TEV buffer** in gravity column.
- 23 Pour cleaved protein onto the beads. Collect flowthrough.
- 24 Wash with leftover **TEV buffer**. Collect flowthrough as well.
- 25 Combine all flowthroughs and concentrate to  **1 mL**
- 26 Equilibrate S200I column with **S200 buffer**.

Other SEC columns can be used instead. S75 should have good separation too.

- 27 Apply concentrated protein to the S200I column and run S200I program.
- 28 Concentrate protein and do buffer exchange with **TEV buffer modified with 10% glycerol**.
- 29 Concentrate protein, quantify, and flash freeze for  **-80 °C** storage.

