

Aug 28, 2024 Version 1

Automated immunodetection of proteinase-K-resistant Lewy Pathology in FFPE human colon tissue V.1

DOI

dx.doi.org/10.17504/protocols.io.kxygxy8ekl8j/v1

Elizabeth Videlock¹


¹UCLA

Neurodegeneration Method Development Community
Tech. support email: ndcn-help@chanzuckerberg.com



Elizabeth Videlock

UCLA



OPEN ACCESS



DOI: dx.doi.org/10.17504/protocols.io.kxygxy8ekl8j/v1

Protocol Citation: Elizabeth Videlock 2024. Automated immunodetection of proteinase-K-resistant Lewy Pathology in FFPE human colon tissue. [protocols.io https://dx.doi.org/10.17504/protocols.io.kxygxy8ekl8j/v1](https://dx.doi.org/10.17504/protocols.io.kxygxy8ekl8j/v1)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 28, 2024

Last Modified: August 28, 2024

Protocol Integer ID: 106608

Funders Acknowledgement:

Chan Zuckerberg Initiative

Grant ID: 2021-235154

Abstract


Automated detection of gut Lewy Pathology using Leica Bond RX processor based on Protocol F using Bond Polymer Refine Detection kit (Leica Biosystems, Cat#: DS9800). Adapted from "Method 5" as described in Beach TG et al., J Parkinsons Dis. 2016.

CITATION

Beach TG, Corbillé AG, Letournel F, Kordower JH, Kremer T, Munoz DG, Intorcia A, Hentz J, Adler CH, Sue LI, Walker J, Serrano G, Derkinderen P (2016). Multicenter Assessment of Immunohistochemical Methods for Pathological Alpha-Synuclein in Sigmoid Colon of Autopsied Parkinson's Disease and Control Subjects..

LINK

<https://doi.org/>



a. Adjacent sections from a healthy control stained for phosphorylated α -Syn (MJFR-13, abcam ab168381) are shown with (left) and without (right) proteinase K treatment demonstrating robust detection of insoluble LP. b. Higher power view of Lewy Pathology (solid arrow) in a PD myenteric plexus ganglion. Samples obtained from Banner Health Brain and Body Donation Program. Scale bars are 50 μ m.


Acknowledgments

We are grateful to the Banner Sun Health Research Institute Brain and Body Donation Program of Sun City, Arizona, for the provision of human colon tissue and data. The Brain and Body Donation Program is supported by the National Institute of Neurological Disorders and Stroke (U24 NS072026 National Brain and Tissue Resource for Parkinsons Disease and Related Disorders), the National Institute on Aging (P30 AG19610 Arizona Alzheimers Disease Core Center), the Arizona Department of Health Services (contract 211002, Arizona Alzheimers Research Center), the Arizona Biomedical Research Commission (contracts 4001, 0011, 05-901 and 1001 to the Arizona Parkinsons Disease Consortium) and the Michael J. Fox Foundation for Parkinsons Research.

We wish to acknowledge Yunfeng Li and the Translational Pathology Core Laboratory (TPCL) at UCLA for assistance with tissue processing, IHC staining, and scanning.

Materials

Note

-  BOND Polymer Refine Detection Kit **Leica Biosystems Catalog #DS9800** contains
- Peroxide Block (30 mL) 3–4% (v/v) Hydrogen peroxide.
 - Polymer (30 mL) Anti-rabbit Poly-HRP-IgG (<25 µg/mL) containing 10% (v/v) animal serum in tris-buffered saline/0.1% ProClin™ 950.
 - DAB Part 1 (2.4 mL) 66 mM 3,3'-Diaminobenzidine tetrahydrochloride hydrate, in a stabilizer solution.
 - DAB Part B (30 mL) ≤0.1% (v/v) Hydrogen Peroxide in a stabilizer solution.
 - Hematoxylin (30 mL) <0.1% Hematoxylin

Equipment

new equipment	NAME
Leica BOND RX	BRAND
3342171	SKU
https://www.leicabiosystems.com/ihc-ish-fish/ihc-ish-instruments/products/leica-bond-rx/ ^{SPECIFICATIONS}	




Protocol materials

- ⊗ Bond™ Wash Solution 10X Concentrate **Leica Biosystems Catalog #AR9590** In [2 steps](#)
- ⊗ Anti-Alpha-synuclein (phospho S129) antibody [MJF-R13 (8-8)] **Abcam Catalog #AB168381** Step 5
- ⊗ Histo-Clear **National Diagnostics Catalog #HS2001GLL** Step 17
- ⊗ Chemical Permunt Mounting Medium **Fisher Scientific Catalog #SP15-100** Step 17
- ⊗ Epitope Retrieval Solution 1 BOND **Leica Biosystems Catalog #AR9961** Step 2
- ⊗ EnVision System-HRP, Labelled Polymer (Rabbit) **Agilent Technologies Catalog #K4003** Step 7
- ⊗ BOND Polymer Refine Detection Kit **Leica Biosystems Catalog #DS9800** Materials
- ⊗ Dewax Solution **Leica Biosystems Catalog #AR9222** Step 1
- ⊗ BOND Enzyme Pretreatment Kit **Leica Biosystems Catalog #AR9551** Step 2




- 1 Leica Bond RX routine factory based “Bake and Dewax” protocol.

 Dewax Solution **Leica Biosystems Catalog #AR9222**


- 2 **Enzyme Antigen Retrieval:** Incubate in Enzyme 1 solution for  00:20:00 at

20m


 Room temperature

 BOND Enzyme Pretreatment Kit **Leica Biosystems Catalog #AR9551**

Note


For comparison, this step can be replaced with Heat Induced Antigen Retrieval using ER1 buffer for  00:20:00

 Epitope Retrieval Solution 1 BOND **Leica Biosystems Catalog #AR9961**

- 3 Incubate in Peroxide Block for  00:05:00


5m



- 4 Incubate in Bond Wash buffer three times  00:00:00 each.

 Bond™ Wash Solution 10X Concentrate **Leica Biosystems Catalog #AR9590**


- 5 **Primary antibody:** Incubate in 1:10,000

1h

 Anti-Alpha-synuclein (phospho S129) antibody [MJF-R13 (8-8)] **Abcam Catalog #AB168381**


for  01:00:00 at  Room temperature

- 6 Incubate in Bond Wash buffer three times  00:00:00 each.

 Bond™ Wash Solution 10X Concentrate **Leica Biosystems Catalog #AR9590**

- 7 Incubate with labeled polymer  00:10:00

10m

 EnVision+ System-HRP, Labelled Polymer (Rabbit) **Agilent Technologies Catalog #K4003**











- 8 Incubate in Bond Wash buffer three times  00:02:00 each.

2m

- 9 Incubate with polymer for  00:08:00 .

8m



- 10 Incubate in Bond Wash buffer five times  00:00:00 each.
- 11 Deionized water wash  00:00:00 .
- 12 Incubate with Mixed DAB Refine for  00:00:00 followed by Mixed DAB Refine for  00:10:00 . 10m
- 13 Deionized water wash for three times  00:00:00 each.
- 14 Incubate with Hematoxylin for  00:10:00 . 10m
- 15 Incubate in Bond Wash buffer three times  00:00:00 each.
- 16 Deionized water wash  00:00:00 .
- 17 Slides are dehydrated in series of alcohols, cleared with  Histo-Clear **National Diagnostics Catalog #HS2001GLL** and mounted with  Chemical Permunt Mounting Medium **Fisher Scientific Catalog #SP15-100**

Protocol references

Beach TG, Corbillé AG, Letournel F, Kordower JH, Kremer T, Munoz DG, Intorcía A, Hentz J, Adler CH, Sue LI, Walker J, Serrano G, Derkinderen P. Multicenter Assessment of Immunohistochemical Methods for Pathological Alpha-Synuclein in Sigmoid Colon of Autopsied Parkinson's Disease and Control Subjects. *J Parkinsons Dis*. 2016 Oct 19;6(4):761-770. doi: 10.3233/JPD-160888. PMID: 27589538; PMCID: PMC5501392.

Citations

Beach TG, Corbillé AG, Letournel F, Kordower JH, Kremer T, Munoz DG, Intorcía A, Hentz J, Adler CH, Sue LI, Walker J, Serrano G, Derkinderen P. Multicenter Assessment of Immunohistochemical Methods for Pathological Alpha-Synuclein in Sigmoid Colon of Autopsied Parkinson's Disease and Control Subjects.

<https://doi.org/>