

Jun 28, 2024

## Crude Subcellular fractionation of FAM177A1-GFP expressing cells

DOI

**dx.doi.org/10.17504/protocols.io.n2bvjn2ypgk5/v1**

Berrak Ugur<sup>1,2,3,4</sup>, Michael G. Hanna<sup>1,2,3,4</sup>, Pietro De Camilli<sup>1,2,3,4</sup>

<sup>1</sup>Departments of Cell Biology and Neuroscience, Yale University School of Medicine, New Haven, CT, USA;

<sup>2</sup>Program in Cellular Neuroscience, Neurodegeneration, and Repair, Yale University School of Medicine, New Haven, CT, USA;

<sup>3</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, USA;

<sup>4</sup>HHMI, Yale University School of Medicine, New Haven, CT, USA

ASAP Collaborative Rese...



Berrak Ugur

Yale University

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.n2bvjn2ypgk5/v1](https://dx.doi.org/10.17504/protocols.io.n2bvjn2ypgk5/v1)

**Protocol Citation:** Berrak Ugur, Michael G. Hanna, Pietro De Camilli 2024. Crude Subcellular fractionation of FAM177A1-GFP expressing cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n2bvjn2ypgk5/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** June 26, 2024

**Last Modified:** June 28, 2024

**Protocol Integer ID:** 102587

**Keywords:** ASAPCRN

**Funders Acknowledgement:****ASAP**

Grant ID: 000580

## Abstract

This protocol details the crude subcellular fractionation of FAM177A1-GFP expressing cells.

## Materials

**Fractionation buffer:**

	A	B
	Tris, pH 7.4	25 mM
	NaCl	150 mM
	Protease inhibitor	




## Crude Subcellular fractionation

2d 1h 5m

1

### Note

All preparations were performed  On ice .

Culture and transfect HeLa cells as described in [dx.doi.org/10.17504/protocols.io.eq2lyp55mlx9/v1](https://dx.doi.org/10.17504/protocols.io.eq2lyp55mlx9/v1)


2

 24:00:00 -  48:00:00 after transfection, wash cells in 6 well plates with PBS.

2d




2.1

Add  200  $\mu$ L PBS (or fractionation buffer) to each well.



2.2

Scrape cells to release from well and transfer to a 1.7 mL eppendorph tube with additional  100  $\mu$ L fractionation buffer wash.


3

Spin the lysate at  1500 rpm, 00:05:00 in a benchtop centrifuge.

5m




4

Remove the supernatant and resuspend cell pellet in  1 mL of cold fractionation buffer.

5

Homogenize resuspended cells with cell cracker (Isobiotec; 8-12 strokes),  1 mL at a time.

5.1


Wash with  1 mL fractionation buffer between samples.



5.2

Always use new syringes for each sample.



- 6 Ultracentrifuge lysates in a Beckman-Coulter table-top ultracentrifuge (TLA100 rotor) at  50000 rpm, 4°C, 01:00:00 to pellet membrane.
- 7 After the centrifugation, transfer the supernatant to a new 1.7 mL Eppendorf tube and save it for western blot analysis.
- 8 Solubilize the membrane fractions from the bottom of the tubes using 4X Laemni buffer for western blot analysis.

1h

