



Oct 30, 2020

© CITE-seq for PBMCs

Nicola Wilson¹, Fernando Calero-Nieto¹, Javier Rodríguez-Ubreva², Bertie Gottgens¹, Esteban Ballestar², Roser Vento-Tormo³

¹Cambridge Institute For Medical Research, United Kingdom;

²Epigenetics and ImmuneDisease Group, Josep Carreras Research Institute (IJC), 08916 Badalona, Barcelona, Spain;

³Wellcome Sanger Institute, Cambridge, UK

1 Works for me

dx.doi.org/10.17504/protocols.io.bjhmkj46

Human Cell Atlas Method Development Community | Coronavirus Method Development Community

Alaina Spivey

ABSTRACT

This is a protocol for CITE-seq for PBMCs. It includes the steps for Thawing samples and Antibody Staining.

DOI

dx.doi.org/10.17504/protocols.io.bjhmkj46

PROTOCOL CITATION

Nicola Wilson, Fernando Calero-Nieto, Javier Rodríguez-Ubreva, Bertie Gottgens, Esteban Ballestar, Roser Vento-Tormo 2020. CITE-seq for PBMCs. **protocols.io**

https://dx.doi.org/10.17504/protocols.io.bjhmkj46

KEYWORDS

CITE-sequencing, PBMCs, COVID

LICENSE

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 09, 2020

LAST MODIFIED

Oct 30, 2020

PROTOCOL INTEGER ID

40205

GUIDELINES

For antibody staining: The aim is to collect a total of \sim 7,000 cells per patient. For this, **25,000 cells of the pool** of 2 samples will be loaded into one channel of the 10x Chromium.

MATERIALS TEXT

MATERIALS

users Catalog #H13680-0040

Solution) BioLegend Catalog #422301, 422302

mprotocols.io

10/30/2020

Citation: Nicola Wilson, Fernando Calero-Nieto, Javier RodrÃÂguez-Ubreva, Bertie Gottgens, Esteban Ballestar, Roser Vento-Tormo (10/30/2020). CITE-seq for PBMCs. https://dx.doi.org/10.17504/protocols.io.bjhmkj46

- RPMI
- 10% FBS
- water bath
- 15 ml tube
- 70% EtOH
- Centrifuge
- Class II MSC
- RBC lysis buffer
- FACS buffer
- Ficoll
- Trypan blue
- CD3+
- Low binding tubes
- MS magnetic column
- CITE-seq Antibody mix
- PBS
- 0.04 % BSA

SAFETY WARNINGS

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

ABSTRACT

This is a protocol for CITE-seq for PBMCs. It includes the steps for Thawing samples and Antibody Staining.

Thawing Samples

- 1 Warm up 10x the volume of thawing medium (RPMI + 10% FBS).
- 2 Thaw frozen PBMCs in § 37 °C water bath.
- 3

Transfer sample into a 15 ml tube.

4

Add 10x volume of thawing medium dropwise into tube.

- 5 Spray down with 70% EtOH.
- 6

Centrifuge **\$500 rcf, 00:05:00**.

7 Remove the sealed bucket from the centrifuge and place into the Class II MSC.

8	Discard supernatant.	
9	Resuspend in □1 mL RBC lysis buffer for ७00:03:00 to remove erytrocytes.	
10	Add 9 ml FACs buffer .	
11	Centrifuge at 3400 rcf , 00:05:00 .	
12	Remove supernatant.	
13	Resuspend the cell pellet in $\ \Box 500\ \mu l$ FACS buffer .	
Equal volu		e for all patients.
14	Perform ficoll to eliminate dead cells.	
	14.1	Add 300 μl Ficoll and 500 μl cell suspension to a centrifuge conical tube.
	14.2	Centrifuge at 10000 rcf , 00:00:10 .
	14.3	Swap the tube and centrifuge again at \$\circ{10000}{10000} \text{rcf, 00:00:30}\$.

15

Recover live cells.

16 Wash with 11 mL FACS buffer . 17 Centrifuge and discard supernatant. 18 Resuspend in 250 mL FACS buffer . 19 Count cells using trypan blue. One fraction of these PBMCs will be used directly and another fraction will be used for CD3+ cells depletion. With the first fraction of PBMCs: pool 125,000 of live PBMCs from each patient (2 patients) in low 19.1 binding tubes; 250,000 cells in total. IMPORTANT: Keep the cells & On ice all times!! 20 Cells should be pooled in a final volume of exactly 100μ . 20.1 Centrifuge and remove $\blacksquare 87.5~\mu l$ supernatant . Final volume $\blacksquare 12.5~\mu l$. Proceed with CD3 negative cells isolation using the **second fraction** of PBMCs: 21 21.1 Use $2x10^6$ of the PBMCs ($200 \mu l$ cell suspension + $20 \mu l$ anti-CD3 magnetic beads). mprotocols.io 10/30/2020

21.2 Incubate **© 00:15:00** at **§ 4 °C** (not ice). 21.3 Resuspend in 1 mL FACS buffer . 21.4 Centrifuge and remove supernatant. 21.5 Resuspend in \$\boxed{\subseteq} 500 \mu I FACS buffer . 21.6 Pass the cells through a MS magnetic column with 3 washes of $\Box 500 \mu I$. 21.7 22 Count CD3 negative cells using trypan blue. 23 Pool 125,000 of CD3 negative cells from each patient (2 patients) in low binding tubes = 250,000 cells in total IMPORTANT: Keep the cells on ice at all times!! 24 Cells should be pooled in a final volume of exactly $\Box 100 \ \mu I$. 25 Centrifuge and remove $\;\; {\color{red}\square}87.5\;\mu l\; supernatant$. Final volume $12.5 \mu l$.

⋈ protocols.io 5 10/30/2020

26 Prepare CITE-seq Antibody mix:

26.1

Use one aliquot of AB pool per sample; centrifuge at @14000 rcf, 4°C, 00:01:00.

26.2

Add **50** µl FACS buffer to reconstitute.

26.3 Vortex.

26.4

Incubate in the dark for **© 00:05:00** to ensure the antibodies are completely resuspended.

26.5

Centrifuge **14000 rcf, 4°C, 00:10:00** .

27

During the centrifugation, add **□1.25** µl Human TruStain FcX™ Fc Blocking Reagent per tube.

28


Incubate © 00:10:00 at & 4 °C.

29

Add 12.5 µl CITE-seq AB mix to FC blocked samples, taking care to avoid any pellet (if one exists).

30

Incubate for **© 00:30:00** at **§ 4 °C** in dark.



⋈ protocols.io 7 10/30/2020

41

Take $\blacksquare 5 \mu l$ cell suspension + $\blacksquare 10 \mu l$ Trypan blue (final dilution 1:3).

42 Proceed immediately to the 10x Genomics Single Cell protocol.

NB: The aim to is collect a total of \sim 7,000 cells per patient. For this, **25,000 cells of the pool** of 2 samples will be loaded into one channel of the 10x Chromium.