



MAR 29, 2024

🌐 Proteolytic Peptides from Conditioned Media Desalting with C18 Hydrophilic–Lipophilic Balance (HLB) Cartridges

🔗 Forked from [Proteolytic Peptide Desalting with C18 Hydrophilic–Lipophilic Balance \(HLB\) Cartridges](#)

Joanna Bons¹, J P Rose¹, M A Watson¹, B Schilling¹

¹Buck Institute for Research on Aging

Cellular Senescence Network (SenNet) Method Development Community

Schilling Lab

OPEN  ACCESS



M A Watson
Buck Institute

ABSTRACT

Desalting of proteolytic peptides from conditioned media elution using C₁₈ Hydrophilic–Lipophilic Balance (HLB) Cartridges in preparation for downstream proteomic profiling.

For desalting of proteolytic peptides from conditioned media Step 4 is modified as we typically reconstitute in an appropriate volume of 0.2% FA rather than by weight/volume due to the input material for digestion also being based on a proportion of the volume of concentrated conditioned media.

DOI:

dx.doi.org/10.17504/protocols.io.eq2lywdzpvx9/v1

Protocol Citation: Joanna Bons, J P Rose, M A Watson, B Schilling 2024. Proteolytic Peptides from Conditioned Media Desalting with C18 Hydrophilic–Lipophilic Balance (HLB) Cartridges . [protocols.io](#) <https://dx.doi.org/10.17504/protocols.io.eq2lywdzpvx9/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 25, 2024

Last Modified: Mar 29, 2024

PROTOCOL integer ID: 97347

Keywords: Desalting, HLB, Peptides, Proteomics, Mass Spectrometry, Conditioned Media, Secretome, SASP

MATERIALS

- Oasis HLB 1cc cartridges, 10 mg sorbent (Waters, Cat. #WAT058882)
- Vacuum manifold
- Centrifuge
- Centrifugal vacuum concentrator
- 1.5-mL microcentrifuge tubes
- Condition buffer/Elution buffer: 50% acetonitrile (ACN), 0.2% formic acid (FA) in water
- Equilibration buffer/Wash buffer: 0.2% formic acid (FA) in water
- HPLC-grade water

- 1 Centrifuge the samples at 1,850 x *g* for 5 minutes at room temperature to pellet insoluble material.
- 2 Desalt the samples using Oasis HLB solid-phase extraction cartridges placed on top of a vacuum manifold as follows:
 - 2.1 Condition each cartridge two times with 800 µL of the condition buffer (50% ACN, 0.2% FA).
 - 2.2 Equilibrate each cartridge three times with 800 µL of the equilibration buffer (0.2% FA).
 - 2.3 Load the peptide samples.

2.4 Wash each cartridge three times with 800 μ L of the wash buffer (0.2% FA).

2.5 Elute peptides with 800 μ L of the elution buffer (50% ACN, 0.2% FA) followed by a further 400 μ L of the same elution buffer.

3 Vacuum dry the eluted peptide solution in a centrifugal vacuum concentrator.

4 Reconstitute the dried peptides in an appropriate volume of 0.2% FA and thoroughly mix the solution.

5 Centrifuge at 12,000 x *g* at room temperature for 2 min, and store at -20°C.