



Jun 14, 2022

DNA extraction (Zymo kit)

Tsu-Chun Hung¹¹KMUdx.doi.org/10.17504/protocols.io.kxygxz11ov8j/v1

Tsu-Chun Hung

DNA extraction (Zymo kit)

DOI

dx.doi.org/10.17504/protocols.io.kxygxz11ov8j/v1

Tsu-Chun Hung 2022. DNA extraction (Zymo kit). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kxygxz11ov8j/v1>






_____ protocol ,

May 24, 2022

Jun 14, 2022

63108


Sample Collection 2m 30s

- 1 Add  **200 μ L** of **1mm beads** to 1.5ml enppendorf tube 30s
- 2 Add  **200 μ L** of **0.5mm beads** to 1.5ml enppendorf tube 30s
- 3 Add  **600 μ L** of **lysis solution** to 1.5ml enppendorf tube 30s


4 Collect  **9-10 µg** of **sample** to 1.5ml enppendorf tube

1m

Sample crush

5 Put 1.5ml eppendorf tube in the mixmill for sample crush; at condition: 30 rpm/s, for 4mins^{4m}
 **00:04:00**

Centrifugation

6 Put 1.5ml eppendorf tube in centrifuge for centrifugation, at this condition:
 **17.0 x g, 25°C, 00:03:00**

3m

DNA purification 37m 30s


7 Add  **600 µL** of **DNA binding buffer** to 1st well of 96 deep well plate

30s

7.1 Add  **25 µL** of **magnetic beads** to 1st well of 96 deep well plate

30s

Shake the bottle and pipetting before using magnetic beads

7.2 Add  **200 µL** of **sample (lysate)** from 15ml centrifuge microtube to 1st well of 96 deep well plate

30s



Usally add last

8 Add  **500 µL** of **Wash 1** to 2nd well of 96 deep well plate

30s

9 Add  **500 µL** of **Wash 2** to 3rd well of 96 deep well plate

30s

- | | | |
|----|--|-----|
| 10 | Add  500 µL of Wash 2 to 4th well of 96 deep well plate | 30s |
| 11 | Add  100 µL of water to 5th well of 96 deep well plate | 30s |
| 12 | Put in the automated DNA extraction machine | 34m |