




MAR 23, 2023

SuperSoil - Soil DNA Extraction

Jie Hao Ou¹, Yin-Tse Huang²

¹National Chung Hsing University, Taichung; ²Kaohsiung Medical University



Jie Hao Ou

ABSTRACT

SuperSoil - Soil DNA Extraction

GUIDELINES

Keep all solutions at room temperature

MATERIALS

1. Bead (ball mill)

A	B
Zirconia Balls	150 μ m : 500 μ m = 1 : 2

2. Solution SD1 (Chaotropic)

A	B	C
	(w/w)	Price
Sodium thiocyanate	5%	1400/250 g
Na ₂ HPO ₄	2.5%	1350/500 g
pH	~8.4	

3. Solution SD2

A	B	C
	(w/w)	Price
AlCl ₃ · 6H ₂ O	1.5%	900/500 g
Ammonium acetate	25%	820/ 500 g
pH (Adjust with acetic acid)	6.3	

4. Solution SD3 (Chaotropic)

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.4r3l27864g1y/v1

Protocol Citation: Jie Hao Ou, Yin-Tse Huang 2023. SuperSoil - Soil DNA Extraction. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4r3l27864g1y/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
 We use this protocol and it's working

Created: Sep 26, 2022

Last Modified: Mar 23, 2023

PROTOCOL integer ID:
 70491

Keywords: Soil DNA extraction

A	B	C
	(w/w)	Price
Guanidinium thiocyanate	65%	1800/250 g
pH	~7.0	

5. Solution SEA (Wash buffer 1)

A	B	C
	(w/w)	Price
Ethanol	50%	1000/4000 mL
Isopropanol	10%	150/ 500 g
Guanidine hydrochloride	30%	2900/1000 g

6.Solution S5 (Wash buffer 2)

A	B	C
	(w/w)	Price
Ethanol	70%	1000/4000 mL

7.Solution S6 (Elution buffer)

A	B	C
	(w/w)	Price
Tris Base	0.12%	3300/1000 g
pH (Adjust with HCl)	8	

SAFETY WARNINGS




Some of the chemicals in this protocol may produce hazardous gases if mixed with bleach.

Lysis

- 1 Add **following materials** to 1.8 ml centrifuge tube
 1. 250 mg soil

2. 1.5 g Bead
3. 800 µl SD1

2 **Vortex** for  00:10:00 at max speed

10m

Note

Vortex horizontally and ensure the cap is tightly sealed

3  14000 rpm, 00:01:00

1m

4 **Transfer the supernatant** to a clean 1.5 ml or 2 ml tube

Note

Expect 500–600 µl

Remove inhibitor

1m

5 Add **200 µl SD2**

Note

White precipitate forms.

If too much AlCl_3 (SD2) is added, it will cause all the DNA to be precipitated along with the inhibitory substances.

See Mustafa (2017) for more info.

CITATION

Irfan Mustafa, Hadiatullah and Sustiyah (2017). Removal of humic acid from peat soils by using AlCl_3 prior to DNA extraction. AIP Conference Proceedings.

LINK

doi.org/10.1063/1.4983434

6 **Vortex** for  00:00:05 at max speed

5s

7  14000 rpm, 00:01:00

1m

Note

Don't worry if the supernatant is yellow.

8 **Transfer the supernatant** to a clean 1.5 ml or 2 ml tube



Note

Expect 500–600 μl



Bind DNA to silica membrane

2m 5s


9 **Add 600 μl SD3**

- 10 **Vortex** for  00:00:05 at max speed 5s
- 11 **Transfer 650 ul lysate** to spin column
- 12  14000 rpm, 00:01:00, discard flow-through 1m
- 13 **Transfer remaining lysate** to spin column

Wash silica membrane 1m

- 14 Add **500 ul SEA** to spin column
- 15  14000 rpm, 00:01:00, discard flow-through 1m
- 16 Add **500 ul S5** to spin column
- 17  14000 rpm, 00:01:00, discard flow-through 1m

18

 14000 rpm, 00:01:00 , Centrifuge again to remove remaining ethanol (S5) , discard flow-through

1m

Elution

19

Add **50–100 ul S6**