



Dec 02, 2021

# 0.05 M Glycine

Jacquelina.Woods<sup>1</sup><sup>1</sup>FDA

1

[dx.doi.org/10.17504/protocols.io.bybapsie](https://dx.doi.org/10.17504/protocols.io.bybapsie)

## GenomeTrakr

Tech. support email: [genomeTrakr@fda.hhs.gov](mailto:genomeTrakr@fda.hhs.gov)Jessica Jones  
US Food and Drug Administration

This method was developed at the FDA's Center for Food Safety and Applied Nutrition for GenomeTrakr's pandemic response project, monitoring SARS-CoV-2 variants in wastewater. Protocols developed for this project cover wastewater collection, concentration, RNA extraction, RT-qPCR detection, library prep, genome sequencing, quality control checks, and data submission to NCBI. This method provides a reagent formula required in the rapid concentration of intact viruses from wastewater using a combination of PEG precipitation and ultracentrifugation.

DOI

[dx.doi.org/10.17504/protocols.io.bybapsie](https://dx.doi.org/10.17504/protocols.io.bybapsie)Jacquelina.Woods 2021. 0.05 M Glycine. **protocols.io**<https://dx.doi.org/10.17504/protocols.io.bybapsie>

protocol ,

Sep 16, 2021

Dec 02, 2021

Sep 16, 2021



Ruth Timme

US Food and Drug Administration

Sep 23, 2021




Jessica Jones

US Food and Drug Administration

53314

In steps of

[Virus Concentration from Wastewater Using PEG Precipitation and Ultracentrifugation](#)



 Glycine Sigma

**Aldrich Catalog #G7126** Step 1

Deionized or ultrapure water

1

 Glycine Sigma

Mix  **3.75 g** **Aldrich Catalog #G7126**  
to  **1 L**

in deionized or ultrapure water

2

Adjust  **9.5**

3

Autoclave  **121 °C**  **00:15:00**

15m

4

Store at  **2-8 °C**