



Sep 02, 2021

# Making electro-competent cells

Ashwinuday <sup>1</sup><sup>1</sup>Indian Institute of Science Education and Research

In Development



Share

[dx.doi.org/10.17504/protocols.io.bxxhppj6](https://dx.doi.org/10.17504/protocols.io.bxxhppj6)

iGEM IISER Pune India 2021

Ashwinuday

## ABSTRACT

This protocol can be used for obtaining of electro-competent bacterial cells.

## DOI

[dx.doi.org/10.17504/protocols.io.bxxhppj6](https://dx.doi.org/10.17504/protocols.io.bxxhppj6)

## PROTOCOL CITATION

Ashwinuday 2021. Making electro-competent cells . **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bxxhppj6>



## LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Sep 02, 2021

## LAST MODIFIED
















Sep 02, 2021

## PROTOCOL INTEGER ID

52937

- 1 Streak the cells on LB agar plates and incubate at  $\uparrow$  37 °C for  $\odot$  12:00:00 12h
- 2 Pick single colony and inoculate it in  $\square$  20 mL 2X LB broth and incubate at  $\uparrow$  37 °C for  $\odot$  12:00:00 1d 4h  
-  $\odot$  16:00:00 on shaker incubator
- 3 Take 2ml of primary inoculum and inoculate it in  $\square$  200 mL 2X LB broth and incubate it at  $\uparrow$  37 °C on shaker.

### 3.1 Grow till the OD reaches 0.6

- 4 Place the culture at  **4 °C** for  **01:30:00** 1h 30m
- 5 Aliquot in  **50 mL** falcon and spin at  **4000 rpm, 4°C, 00:25:00** 25m
- 6 Discard the supernatant and resuspend well in  **40 mL** 10% glycerol
- 7 Spin at  **4000 rpm, 4°C, 00:25:00** 25m
- 8 Discard the supernatant and resuspend in  **25 mL** 10% glycerol
- 9 Spin at  **4000 rpm, 4°C, 00:25:00** 25m
- 10 Discard the supernatant and resuspend in  **20 mL** 10% glycerol. Mix the contents of 2 falcons together such that only 2 falcons of  **40 mL** contents is present
- 11 Spin at  **4000 rpm, 4°C, 00:25:00** 25m
- 12 Discard the supernatant and resuspend in  **1 mL** 10% glycerol.
- 13 Make aliquots of  **50 µl** in autoclaved  **600 µl** eppendorfs, flash freeze and store at  **-80 °C**