



Preparation of electrocompetent *Escherichia coli* V.1

lewis.bingle¹

¹University of Sunderland

Version 1

Oct 09, 2020

1

Works for me

This protocol is published without a DOI.

SunMedMicro

lewis.bingle

ABSTRACT

Preparation of competent E.coli cells for electroporation

PROTOCOL CITATION

lewis.bingle 2020. Preparation of electrocompetent *Escherichia coli*. protocols.io
<https://protocols.io/view/preparation-of-electrocompetent-escherichia-coli-bm7zk9p6>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 09, 2020

LAST MODIFIED

Oct 09, 2020

PROTOCOL INTEGER ID

42969

GUIDELINES

1. After harvesting and chilling cells (steps 1-2) they should be kept as cold as possible during the rest of the procedure – e.g. keep all tubes and buffers on ice whenever possible; do not hold in the hand for long periods; handle tubes by top (lid) rather than bottom.
2. As cells are washed, the pellet sticks increasingly poorly to the wall of the tube after centrifugation – it may be necessary to increase speed or duration of centrifugation step.

MATERIALS TEXT

10 % glycerol in distilled / purified water.

ABSTRACT

Preparation of competent E.coli cells for electroporation

BEFORE STARTING

Prechill 10 % glycerol solution (or SDW) on ice or in fridge and keep on ice during use. Prechill centrifuge tubes on ice. Prechill centrifuge (if refrigerated) or put rotor / buckets in fridge.

- 1 Inoculate starter culture of the desired strain from single colony (e.g. 5 mL LB broth in universal bottle). Incubate (37°C with shaking at 200-220 rpm) Overnight .
- 2 Dilute overnight culture 1 / 100 into 25 mL fresh LB broth (supplemented with appropriate antibiotic if required) and incubate (37°C with shaking at 200-220 rpm) to mid-log phase growth (OD_{600nm} of between 0.5 and 0.7).

- 3 Decant bacteria into 50 ml Falcon (centrifuge) tube and chill on ice (approx. ⌚00:20:00).
- 4 Centrifuge 🌀4000 x g, 4°C for ⌚00:10:00 (use refrigerated centrifuge if possible, or pre-chill centrifuge rotor / bucket). * Discard supernatant.
- 5 Wash the cells 3x with ice-cold, sterile aqueous glycerol solution [m]10 % (w/v) (can use sterile distilled water) in decreasing volumes of 📏25 mL , 📏12 mL and 📏6 mL (centrifuge as in step 4 between each wash).
- 6 Resuspend final cell pellet in 📏300 µl of sterile 10% (w/v) glycerol (or water) 🧊 On ice .
- 7 Store electrocompetent cells on ice for immediate use or snap-freeze in liquid nitrogen for long-term storage at 🧊 -80 °C .