



Oct 03, 2022

## CompetentCell

Ana Belem García González<sup>1</sup>, Georgina Diego<sup>1</sup>,  
Irán Alessandra Chaparro Rodríguez<sup>1</sup>, Jair Alexis Gardea Sáenz<sup>1</sup>

<sup>1</sup>Tecnologico de Monterrey Campus Chihuahua

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.j8nlkkzd6l5r/v1](https://dx.doi.org/10.17504/protocols.io.j8nlkkzd6l5r/v1)

Tec-Chihuahua\_2022\_Agrocapsi

**D** Denisse Chacón Ramírez

### ABSTRACT

A protocol for making chemically competent cells is shown.

### DOI

[dx.doi.org/10.17504/protocols.io.j8nlkkzd6l5r/v1](https://dx.doi.org/10.17504/protocols.io.j8nlkkzd6l5r/v1)

### PROTOCOL CITATION

Ana Belem García González, Georgina Diego, Irán Alessandra Chaparro Rodríguez, Jair Alexis Gardea Sáenz 2022. CompetentCell. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.j8nlkkzd6l5r/v1>



### LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

### CREATED

Aug 04, 2022

### LAST MODIFIED

Oct 03, 2022

### PROTOCOL INTEGER ID

68197



















#### MATERIALS TEXT











- LB Broth
- CaCl<sub>2</sub> 50 mM Cold
- CaCl<sub>2</sub> 100 mM Cold
- 50 mL Falcon Tubes
- *E. coli* cells

#### Equipment:

- Spectrophotometer
- Centrifuge
- Incubator at 37 °C
- Shaker
- Water Bath at 37 °C
- Cooler with ice

#### Competent cells protocol with CaCl<sub>2</sub> 1h 20m

- 1 Inoculate  100 µL of *E. coli* in  10 mL of antibiotic-free starter LB culture medium.  
Incubate during  Overnight at  37 °C and  225 rpm 
- 2 Heat the LB media in the water bath at  37 °C
- 3 Inoculate  40 mL of LB media with  1 mL of cell culture from step 1.
- 4 Incubate at  37 °C with constant agitation (  225 rpm during  00:30:00 ). 30m
- 5 Measure de optical density at 600 nm every  00:30:00 until an OD600 of 0.4-0.6 is reached. 30m
- 6 Chill for  00:20:00  On ice 20m
- 7 Centrifuge during  00:15:00 at  4 °C and  3000 rpm 15m

- 8 Discard the supernatant, add  **20 mL** of  **CaCl<sub>2</sub> Contributed by users** solution 50 mM and resuspend carefully.
- 9 Chill for  **00:20:00**  **On ice** 20m
- 10 Centrifuge during  **00:15:00** at  **4 °C** and  **3000 rpm** 15m
- 11 Discard the supernatant, add 5 mL of CaCl<sub>2</sub> solution 100 mM and resuspend carefully. Chill <sup>20m</sup>  
 **On ice** for  **00:20:00**
- 12 Preserve at  **4 °C**