



# First-strand synthesis and touchdown PCR for SARS-CoV-2 V.2

Leonardo Caserta<sup>1</sup>

<sup>1</sup>Cornell University

Version 2

Feb 05, 2021

1 Works for me dx.doi.org/10.17504/protocols.io.br54m88w

Diel Lab



Leonardo Caserta

DOI

[dx.doi.org/10.17504/protocols.io.br54m88w](https://dx.doi.org/10.17504/protocols.io.br54m88w)

DOCUMENT CITATION

Leonardo Caserta 2021. First-strand synthesis and touchdown PCR for SARS-CoV-2. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.br54m88w>  
Version created by [Leonardo Caserta](#)

LICENSE

This is an open access document distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 05, 2021

LAST MODIFIED

Feb 05, 2021

DOCUMENT INTEGER ID

47004

## Annealing primer to RNA

Note: Pre-warm 5X SS IV Buffer at room temperature before using

Combine the following in a tube:

A	B
	<b>(13 µl)</b>
Random primers 6 (NEB)	1 µl (60ng/µl)
10 mM dNTP mix (10 mM each)	1 µl
Template RNA	11 µl

1. Mix by pipetting and briefly centrifuge the components.

2. Heat the RNA-primer mix at 65°C for 5 minutes, and then incubate on ice for at least 1 minute.

### Prepare RT reaction mix

1. Vortex and briefly centrifuge the 5× SSIV Buffer.

2. Combine the following into a reaction tube:

**(7 µl)**

5× First Strand buffer 4 µl

100 mM DTT 1 µl

RNAse out 1 µl

SS IV Reverse Transcriptase (200 U/µL) 1 µl

Cap the tube, mix, and then briefly centrifuge the contents.

### Reverse Transcriptase step

1. Add RT reaction mix to the annealed RNA, mix well

2. Incubate under the following conditions:

23 °C for 10 minutes

50°C for 10 minutes

55°C for 10 minutes

80°C for 10 minutes

4 °C ∞

### Multiplex PCR

Set up two multiplex PCR reaction mix as follows:

Q5 High-Fidelity 2X Master Mix	12.5 µl
H2O	6.75 µl
Primers pool (1 or 2) at 10 µM	1.75 µl
cDNA	4 µl
<b>TOTAL</b>	<b>25 µl</b>

Cycling conditions for **Q5 High-Fidelity 2X Master Mix (~1500bp)**

A	B	C
Temp.	Time	Number of cycles
98	30 sec	1
98	10 sec	4
65	2 min	
72	3 min	
98	10 sec	3
64	2 min	
72	3 min	
98	10 sec	2
63	2 min	
72	3 min	
98	10 sec	2
62	2 min	
72	3 min	
98	10 sec	2
61	2 min	
72	3 min	
98	10 sec	22
60	2 min	
72	3 min	
72	3 min	1
4	∞	