



Nov 06, 2020

Total protein in microalgae: Pierce BCA protein assay

Yingyu YY Hu¹, Zoe V. Finkel¹¹Dalhousie University**1** Works for me This protocol is published without a DOI.**Marine Microbial Macroecology Lab**
Tech. support email: ruby.hu@dal.caYingyu Hu
Dalhousie University

ABSTRACT

Here we describe a protocol for extracting total protein from microalgae and quantifying total protein by Pierce BCA protein assay.

https://assets.thermofisher.com/TFS-Assets/LSG/manuals/MAN0011430_Pierce_BCA_Protein_Asy_UG.pdf

PROTOCOL CITATION

Yingyu YY Hu, Zoe V. Finkel 2020. Total protein in microalgae: Pierce BCA protein assay. **protocols.io**
<https://protocols.io/view/total-protein-in-microalgae-pierce-bca-protein-ass-bnd4ma8w>

KEYWORDS

microalgae, total protein, Pierce BCA, protein solubilization buffer, bead mill cell disruption, microplate

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 13, 2020

LAST MODIFIED

Nov 06, 2020

PROTOCOL INTEGER ID

43164

GUIDELINES

1. Considering the working range of Pierce BCA assay is 20-2000 ug/ml protein, an ideal protein biomass for each sample is around 200 ug/filter extracted by 1 ml protein extraction buffer. The minimum requirement of protein is no less than 50 ug (extracted by 500 ul protein extraction buffer).
2. Under replete condition, assuming protein is 32% in dry mass, 50 ug of protein requires around 160 ug of dry mass, which is about 60 ug total carbon for each sample.

ABSTRACT

Here we describe a protocol for extracting total protein from microalgae and quantifying total protein by Pierce BCA protein assay.

https://assets.thermofisher.com/TFS-Assets/LSG/manuals/MAN0011430_Pierce_BCA_Protein_Asy_UG.pdf

Sample collection

- 1 Calculate the volume to obtain enough biomass for the assay (see guidelines)

2 Use 5 inches Hg to filter microalgae samples onto polycarbonate filters (pore size varies).

3 Fold filter in quarter and place into an MP biomedical bead tube.

- Use 2 ml bead tube for 25 mm filter; 15 ml bead tube for 47 mm filter
- Large ceramic beads are recommended for diatoms and larger phytoplankton cells (2 mL: MP116913500/ 15 mL: MP116933050)
- Small yttrium beads are recommended for picoplankton (2 mL: MP116960100/ 15 mL: MP116975050)

LYSING TUBES

MATRIX D 2 mL/15 mL

MP BIOMEDICALS 116913500/116933050

Lysing tubes

Matrix Y 2 mL/15 mL

MP MP116960100/MP11697505
Biomedicals 0

4 Flash-freeze tubes in thermo flask with liquid nitrogen

5 Store at  **-80 °C** until further processing

Prepare protein solubilization buffer (PSB)

6

Ni G, Zimbalatti G, Murphy CD, Barnett AB, Arsenault CM, Li G, Cockshutt AM, Campbell DA (2017). Arctic Micromonas uses protein pools and non-photochemical quenching to cope with temperature restrictions on Photosystem II protein turnover.. Photosynthesis research.
<https://doi.org/10.1007/s11120-016-0310-6>

7 In order to obtain compatible results, prepare sufficient PSB so that the same PSB can be used for sample extraction, blank filter extraction and standard solutions

- (1) Extract all samples: Each sample requires 0.25 mL PSB
- (2) Extract all blank filters: Each filter requires 0.25 mL PSB

(3) Each standard solution (500 ul) requires 0.125 mL PSB

8 For each 10 g PSB

8.1 Use antistatics weighing dish to weigh the following chemicals:

Antistatic weighing dish
Fisherbrand 08-732-112

(1) 0.136 g Tris base

[Tris](#)

[base Bioshop Catalog #TRS001.500](#)

(2) 0.133 g Tris HCl

[Tris](#)

[HCl Bioshop Catalog #TRS002.500](#)

(3) 0.8 g Lithium dodecyl sulphate

[Lithium Dodecyl](#)

[Sulphate Bioshop Catalog #LDS701.25](#)

8.2 Place a plastic beaker on the top of the scale surface

8.3 Remove the cap of a 15 mL tube and sit it in the beaker

Falcon® Centrifuge Tubes
Polypropylene, Sterile, 15 mL
Corning® 352096

8.4 Tare the total weight of beaker and tube

8.5 Transfer all chemicals weighed in [go to step #8.1](#) into the tube, rinse the dish with small amount of milliQ water to make certain all of the solute is transferred into the tube

8.6 Use a transfer pipet to add 4 g glycerol into the tube

[Glycerol Bioshop Catalog #GLY001.500](#)

8.7 Add 40 µl 0.5 M EDTA into the tube

[EDTA buffer solution \(0.5 M\) Sigma](#)

Aldrich Catalog #4055-100ml

8.8 Top to 10 g with milliQ water

8.9 Mix by vortex or shaker until all solute is completely dissolved.

Prepare Pefabloc solution

9

[4-\(2-2-aminoethyl\)-benzenesulfonyl fluoride HCL \(AEBSF,](#)

[Pefabloc\) Bioshop Catalog #AEB602.100](#)

Pefabloc is a protease inhibitor, and it loses activity over 24 hours.

10 Add 20.86 mL milliQ into 100 mg Pefabloc to obtain a final concentration of 20 mM.

11 Aliquot into 2.5 ml portions and keep frozen at **-20 °C**

12 The solution can be frozen~thawed multiple times.

Extract protein from microalgae

13 Prepare protein extraction buffer (PEB):

In each 1 mL PEB:

250 µl PSB

20 µl 20 mM Pefabloc

730 µl milliQ water

14 Prepare ice-bath, keep all samples in the ice-bath

15 If samples and blank filters are in bead tubes, add 1 mL PEB to each tube directly.

Volume of PEB varies due to the actual biomass collected (see guideline)

16 If samples are in cryo tubes instead of bead tubes

16.1 Rinse forceps with 70% ethanol and air dry

Filter forceps
blunt end, stainless steel
Millipore XX6200006P

16.2 Label bead tubes and add 1 mL PEB into each bead tube.

Volume of PEB varies due to the actual biomass collected (see guideline)

16.3 Use clean forceps to transfer samples and blank filters into its corresponding bead tube.

17 Turn on FastPrep








FastPrep-24 5G
Bead beater
MP Biomedicals 116005500 [↗](#)

18 Check the cap of each tube to make certain cap is tightly screwed. Organize the tubes in order, take notes of the position of each tube, in case the labels get rubbed out during extraction.

19 Run 1 min at 6.5 m/s

20 Keep tubes On ice for 00:01:00

21 Check labels. Put tubes back into FastPrep.

- 22 Run 1 min at 6.5 m/s
- 23 Keep tubes  **On ice** for  **00:01:00**
- 24 Check labels. Put tubes back into FastPrep.
- 25 Run 1 min at 6.5 m/s
- 26 Keep tubes  **On ice** for  **00:01:00**
- 27 Check labels. Put tubes back into FastPrep.
- 28 Run 1 min at 6.5 m/s
- 29 Turn on refrigerated centrifuge
- 30 Keep tubes  **On ice** for  **00:01:00**
- 31 Centrifuge tubes at  **14800 x g, 4°C, 00:05:00**
- 32 Transfer 100 ul supernatant to a 600 ul microtube

CENTRIFUGE 5430 R
Eppendorf MP2231000510

Microcentrifuge Tubes
1.7 mL/0.6 mL
Axygen Scientific MCT-175-C/MCT-060-L-C

Sample volume for Pierce BCA is 25 ul for each assay. Considering that sample is dispensed into microplate by reverse pipetting, it requires at least 100 ul to dispense 2X25 ul sample. If triplicate is preferred, transfer 150 ul supernatant.

- 33 Freeze both the supernatant and the rest at -80°C

Prepare Bovine serum albumin (BSA) standard solutions

34

[BSA 2 mg/mL standard](#) Thermo Fisher

Scientific Catalog #23209

- 35 Thaw 20 mM pefabloc and transfer 150 ul to a 600 ul microtube. Put the rest of the stock back into the freezer immediately.

- 36 Organize eight 1.7 mL microtubes in the tube rack, label the tubes from SD1 to SD8.

- 37 Use 0.5 mL Finntip Stepper Tip to dispense 10 ul pefabloc into each microtube

Finntip™ Stepper Pipette Tips
500 uL
Thermo Scientific™ 9404170

Wipe or dab the liquid drop on the outside of the tip, avoid wiping the tip open before dispensing the liquid.

- 38 Use 1.25 mL Finntip Stepper Tip (Thermo Fisher 9404180) to dispense 125 ul PSB into each microtube.

Finntip stepper tip
1.25 mL
Thermo Fisher 9404180

- 39 Add milliQ into each microtube according to the sheet below:

Standard	Pefabloc (20 mM) (uL)	PSB (uL)	MQ (uL)	BSA (2mg/mL) (uL)	Final Conc. (mg/mL)
----------	-----------------------	----------	---------	-------------------	---------------------

SD1	10	125	365	0	0
SD2	10	125	360	5	0.02
SD3	10	125	353	12	0.048
SD4	10	125	340	25	0.1
SD5	10	125	315	50	0.2
SD6	10	125	265	100	0.4
SD7	10	125	165	200	0.8
SD8	10	125	115	250	1

40 Break the ampule of BSA standard

SCIENCEWARE® Break-Safe™ Ampule
Opener
Bel-Art® 89217-378

41 Reverse pipet certain amount of BSA (2 mg/mL) into each tube according to the sheet [go to step #39](#)

Wipe or dab the liquid drop on the outside of the tip, avoid wiping the tip open before dispensing the liquid.

42 Vortex each tube.

43 Dispense 4 μ L of each standard solution onto microdrop plate.


μ Drop™ Plates
Thermo Scientific N12391

44 Read absorbance of eight standard solutions at 205 nm

Varioskan LUX Multimode Microplate
Reader
Thermo Fisher VL0L00D0

45 Subtract absorbance at 205 nm of blank standard from the 205 nm measurements of all other standard solutions

- 46 Plot the blank-corrected 205 nm measurement for each standard solution versus its concentration in mg/mL.



- 47 If the standard curve has good Coefficient of Determination, i.e. $R^2 > 0.99$, the standard solutions are in good quality; otherwise, prepare a new series of standard solutions until the quality of standard solutions passes the screen.

Prepare BCA working reagent (WR)

- 48 Use the following formula to determine the total volume of WR required. Consider several mL of extra volume since Finn timer is unable to expel the entire volume from the tip
(# standards + # samples + # blank filters) X (# replicates) X (200 μ l) = total volume WR required
- 49 Prepare WR by mixing 50 parts of BCA reagent A with 1 part of BCA Reagent B in a 50 mL falcon tube

 **Pierce BCA Protein Assay Kit Thermo Fisher**

Scientific Catalog #23225

Falcon® Centrifuge Tubes
Polypropylene, Sterile, 50 mL
Corning® 352070

Pierce BCA assay 2h

- 50 Turn on incubator and preheat to **37 °C**

SHAKING INCUBATOR
71L
Corning® LSE™ 6753

- 51 Prepare 96-well microplate and lid:
Each microplate can load eight standard solutions and forty samples+blanks, duplicated

96-Well Microplates
Polystyrene, Clear,
Greiner Bio-One 82050-760

Microplate Lids
Polystyrene
Greiner Bio-One 07000288

52

	<u>1</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>	<u>6</u>	<u>7</u>	<u>8</u>	<u>9</u>	<u>10</u>	<u>11</u>	<u>12</u>
<u>A</u>	S1	S1	Samples and sample blanks: 40 with duplicate									
<u>B</u>	S2	S2										
<u>C</u>	S3	S3										
<u>D</u>	S4	S4										
<u>E</u>	S5	S5										
<u>F</u>	S6	S6										
<u>G</u>	S7	S7										
<u>H</u>	S8	S8										

- 53 Reverse pipette 25 ul of each standard or sample/blank replicate into the microplate well.

- 54 Use 2.5 mL Finntip stepper tip to load 200 ul WR to each well.

Finntip™ Stepper Pipette Tips
2.5 mL
Thermo Fisher 9404190

Wipe or dab the liquid drop on the outside of the tip, avoid wiping the tip open before dispensing the liquid.

Dispense WR column by column.
Refill the stepper in between two columns of wells.

55 Cover plate with lid.

56 Shake and incubate at  37 °C for  00:30:00

57 Read absorbance at 562 nm

Varioskan LUX Multimode Microplate
Reader

Thermo Fisher VL0L00D0

58 Put plate back into incubator. Shake and incubate for another  01:30:00

59 Read absorbance at 562 nm

Calculate protein content per filter

60 Subtract the average 562 nm absorbance measurement of the blank standard replicates from the 562 nm measurements of all other individual **standard**.


61 Subtract the average 562 nm absorbance measurement of the blank sample (filter) replicates from the 562 nm measurements of all other individual **sample**.

62 Prepare a standard curve by plotting the average Blank-corrected 562 nm measurement for each BSA standard versus its concentration in mg/ml.

If using curve-fitting algorithms associated with a microplate reader, a four-parameter (quadratic) or best-fit curve will provide more accurate results than a purely linear fit.

63 Use the standard curve to determine the protein concentration of each unknown sample by using its blank-corrected 562 absorbance.

64



Increasing the incubation time can increase the net 562 nm measurement for each well and lower the minimum detection of the assay.

65 Protein_mg/filter = Protein_mg/mL X PEB_mL