

Dec 11, 2024

# spotPCR: A Rapid and Efficient Approach for Indexing Individual Template Molecules using Unique Molecular Identifiers

DOI

**[dx.doi.org/10.17504/protocols.io.261ged66wv47/v1](https://dx.doi.org/10.17504/protocols.io.261ged66wv47/v1)**

Jason D Limberis<sup>1</sup>

<sup>1</sup>University of California, San Francisco



Jason D Limberis

University of California, San Francisco

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.261ged66wv47/v1](https://dx.doi.org/10.17504/protocols.io.261ged66wv47/v1)**

**Protocol Citation:** Jason D Limberis 2024. spotPCR: A Rapid and Efficient Approach for Indexing Individual Template Molecules using Unique Molecular Identifiers. **protocols.io** **<https://dx.doi.org/10.17504/protocols.io.261ged66wv47/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** June 16, 2023

**Last Modified:** December 11, 2024

**Protocol Integer ID:** 83554

**Funders Acknowledgement:**

**John Metcalfe**

**Grant ID:** R01AI177637

## Abstract

Low-frequency mutations provide valuable insights in various fields, including drug resistance identification, cancer and infectious disease research. One promising strategy to enhance the sensitivity and specificity of mutation detection is the incorporation of unique molecular identifiers (UMIs) during polymerase chain reaction (PCR) amplification and before deep sequencing. However, conventional methods for UMI incorporation often necessitate multiple labor-intensive steps. spotPCR (Specific Primer Limited Unique Molecular Identifier Tagging PCR) overcomes these challenges, streamlining the UMI tagging process.

## Materials

### Required

✕ Q5 Hot Start High-Fidelity DNA Polymerase - 500 units **New England Biolabs Catalog #M0493L**

✕ Agencourt AmPure XP beads **Contributed by users Catalog #A63880**

✕ dNTPs **Contributed by users**

A thermocycler and a qPCR machine

A magnetic rack

### Optional

✕ NEBNext Library Quant Kit for Illumina - 500 rxns **New England Biolabs Catalog #E7630L**

A	B	C
Primer Set	Direction	Sequence
pncA	F	TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCGCGGCGTCATGGAC CCTAT
pncA	R	GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGNNNNNNNNNNNTTTC GAAGCCGCTGTACGCTCC
gyrA	F	TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGTCACCCGCAACGCCA AGGAT
gyrA	R	GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGNNNNNNNNNNNTTAT TGCCTGGCGAGCCGAAGT
Illumina adapter primer	F	CAAGCAGAAGACGGCATACGAGAT[i7]GTCTCGTGGGCTCGGAGATGTG TATAAGAGACAG
Illumina adapter primer	R	AATGATACGGCGACCAACGAGATCTACAC[i5]TCGTCGGCAGCGTCAGAT GTGTATAAGAGACAG

You can use any compatible Illumina adapters such as the IDT for Illumina DNA/RNA UD Indexes.



## Stage 1 PCR

1

A	B
COMPONENT	Volume (μl)
5X Q5 Reaction Buffer	5
5X Q5 High GC Buffer	5
10 mM dNTPs	0.5
Q5 High-Fidelity DNA Polymerase	0.25
1μM Forward primer 1	1
1μM Forward primer 2	1
0.0625μM Reverse primer 1*	1
0.625μM Reverse primer 2*	1
Template DNA	2
Nuclease-Free Water	6.25

\*We recommend doing a dilution series starting at 500fM to asses new primersets.

A	B	C	D
Step	Temp (C)	Time (s)	Cycles
Denaturation	98	120	1
Denaturation	98	10	3
Annealing	62	45	
Extension	72	120	
Extension	4	Forever	1

Cycle parameters. Three cycles are performed to ensure there are amplicon copies with Illumina tags on both ends.

## Stage 2 PCR

2 Add the below into the reaction while at  4 °C and proceed to the second PCR cycling

A	B
COMPONENT	Volume (μl)
Illumina adapter primer F (10uM)	1
Illumina adapter primer R (10uM)	1






A	B	C	D
Step	Temp (C)	Time (s)	Cycles
Denaturation	98	120	1
Denaturation	98	5	34
Annealing	65	15	
Extension	72	20	
Extension	4	Forever	1

Cycle parameters

## Bead cleanup



21m

3 Add  25  $\mu$ L (1X) of resuspended AMPure XP Beads to the sample  
Mix by pipetting 10x

4 Incubate  00:02:00 at  Room temperature

2m


5 Place on the magnet, allow the beads to aggregate, and remove and discard the supernatant

6 Add  200  $\mu$ L [M] 70 % (v/v) ethanol and incubate (still on the magnet) for  
 00:00:30


30s

6.1 Remove the supernatant



6.2 Repeat  go to step #6 for a total of 2 washes

7 Air dry for  00:00:30 , don't allow the beads to become cracked






30s

8 Immediately after the bead pellet becomes opaque, remove the tube from magnetic rack and resuspend in  20  $\mu$ L of **Low EDTA Tris Buffer**. Ensure all beads are in solution.



- 9 Incubate at room temperature for  00:05:00
- 10 Place on magnetic rack, wait for the solution to become clear ~  00:02:00 , and transfer the eluted DNA to a new well-labeled tube

### Optional: NEB Illumina Quantification

- 11 Thaw the NEBNext Library Quant Master Mix and NEBNext Library Quant Primer Mix. Ensure mixing of NEBNext Library Quant Primer Mix by vortexing. Place reagents on ice.
- 12 Thaw the NEBNext Library Quant DNA Standards, tubes 1–6.  
Mix by pulse vortexing on a low setting. Briefly spin to collect material from the sides of the tubes. Place on ice.
- 13 Thaw the NEBNext Library Quant Dilution Buffer (10X). Mix well by vortexing.  
Centrifuge briefly to collect material from the sides of the tube.  
Place on ice.
- 14 Add  100  $\mu$ L NEBNext Library Quant Primer Mix to the tube of NEBNext Library Quant Master Mix (  1.5 mL ). Mix by vortexing. Write the date on the master mix tube to indicate that primer mix has been added.
- 15 Dilute the NEBNext Library Quant Dilution Buffer (10X) 1:10 with nuclease-free water. Mix by vortexing.  
Prepare sufficient buffer for quantitating the desired number of libraries, allowing  1.2 mL for each library.
- 16 Prepare a 1:1,000 dilution of each library sample in NEBNext Library Quant Dilution Buffer (1X)
- 17 Aliquot  16  $\mu$ L NEBNext Library Quant Master Mix (with primers) to each well
- 18 Add-in  4  $\mu$ L of sample or standard per well
- 19



A	B	C	D
Step	Temp (C)	Time (s)	Cycles
Denaturation	95	60	1
Denaturation	95	15	35
Annealing	63	45	

Cycle parameters

A denaturation/melt curve can be included if desired, but is optional.

A	B
Sample	Conc. (pM)
DNA Standard 1	100
DNA Standard 2	10
DNA Standard 3	1
DNA Standard 4	0.1
DNA Standard 5	0.01
DNA Standard 6	0.001

20 Adjusted Conc. = Calculated Conc. × 399 / library size (bp)

## Analysis

21 **[Link to script \(click here\)](#)**

```
bash spotPCR_process.sh \  
  -R1 "Read1.fastq" \  
  -R2 "Read2.fastq" \  
  -Ref_name "Reference.fasta" \  
  -sample_name "SampleName" \  
  -threads "Threads" \  
  -umi "UMI_Barcode_Pattern"
```

The output will allow you to determine if there was sufficient UNI tagging of the amplicons or if a lower initial primer concentration is needed.

