

May 23, 2025

Non-enzymatic passaging of hPSC

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.5qpvo32k7v4o/v1



Valeria Fernandez Vallone¹, Lyn Healy², Nathalie Lefort³, Katarzyna Ludwik¹, Tamer Onder⁴, Lisa Pavinato^{5,6}, Fatma Visal FVO OKUR⁷, Harald Stachelscheid¹

¹Core Unit pluripotent Stem Cells and Organoids - Berlin Institute of Health @ Charite, Berlin, Germany;

²The Francis Crick Institute;

³Université de Paris, Imagine Institute, iPSC Core Facility, INSERM UMR U1163, F-75015 Paris, France.;

⁴Koc University;

⁵Institute of Oncology Research (IOR), Bellinzona Institutes of Science (BIOS+), Bellinzona, Switzerland;

⁶Faculty of Biomedical Sciences, Università della Svizzera Italiana, Lugano, Switzerland.;

⁷Hacettepe University, Center for Stem Cell Research and Development (PEDISTEM) and Hacettepe University Faculty of Medicine, Department of Pediatrics



Harald Stachelscheid

Core Unit pluripotent Stem Cells and Organoids - Berlin Inst...

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.5qpvo32k7v4o/v1

Protocol Citation: Valeria Fernandez Vallone, Lyn Healy, Nathalie Lefort, Katarzyna Ludwik, Tamer Onder, Lisa Pavinato, Fatma Visal FVO OKUR, Harald Stachelscheid 2025. Non-enzymatic passaging of hPSC. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.5qpvo32k7v4o/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 20, 2023



Last Modified: May 23, 2025

Protocol Integer ID: 91188

Keywords: passaging, EDTA, hPSC, iPSC, non-enzymatic, clumps, split

Funders Acknowledgements:

COST Action CorEuStem

Grant ID: CA20140

Abstract

This protocol describes the procedure to passage hPSC culture using enzyme-free reagents.

Guidelines


In this protocol, hPSC passaging is done in aggregate/clumps format. The protocol is adaptable to different enzyme-free reagents described in Materials section and variable plate formats.



Materials

LABORATORY EQUIPMENT AND CONSUMABLES

Use sterile material

- 1/5/10 mL pipettes
- 15/50 mL conical tubes
- Cell culture treated plastic vessels of choice e.g. 24, 12 or 6-well plates, T25, T75 flasks, 10cm dishes
- 10/200/1000µL tips and micropipettes (optional)
- Cell scraper
- Aspirator pump with disposable pipette
- Sterile filter unit with  0.22 µm filter
- Microscope, if available Stereo Microscope
- Incubator at 37°C and 5% CO₂ or under hypoxic conditions, 5% O₂/ 5% CO₂
- Class II Biosafety Cabinet

MEDIA AND REAGENTS

 EDTA (0.5 M), pH 8.0 Life Technologies Catalog #AM9260G

 ReLeSR™ 100 mL STEMCELL Technologies Inc. Catalog #5872

 Versene Solution Thermo Fisher Catalog #15040033

 Gentle Cell Dissociation Reagent 100 mL STEMCELL Technologies Inc. Catalog #7174

 DPBS no calcium, no magnesium Invitrogen - Thermo Fisher Catalog #14190136

 UltraPure Distilled Water Invitrogen - Thermo Fisher Catalog #10977-015

hPSC culture conditions and survival factors choice depend on hPSC line and individual lab practices. For options refer to protocols:




Maintenance of hPSC

Coating of tissue Culture Vessels for hPSC



Before start

If using EDTA as enzyme-free reagent, prepare EDTA [M] 0.5 millimolar (mM) as follows:

1. Dilute 1 to 1000 EDTA [M] 0.5 Molarity (M) in distilled water; e.g.  5 mL EDTA [M] 0.5 Molarity (M) in  995 mL distilled water.
2. Filter sterile using  0.22 μm filter.
3. Store at RT.



Preparation of destination vessel

5m

- 1 Prepare coated tissue culture vessels and culture media according to hPSC line requirements and desired format.
Refer to protocols: [Coating of tissue Culture Vessels for hPSC](#) and [Maintenance of hPSC](#).
- 2 Aspirate and discard coating solution from wells of destination vessel.
- 3 Add hPSC maintenance media per destination vessel, refer to **Table 1** for recommended volumes.

5m

A	B	C	D	E
Culture Vessel	DPBS (mL)	non-enzymatic reagent (mL)	Media for hPSC harvesting (mL)	Media in destination vessel (mL)
12 WP (per well)	0.5	0.5	1.5	1
6 WP (per well)	1	1	2.5	2
T25	2	2	5	7
10 cm	3	3	7.5	10
T75	5	5	10	15

Table 1. Recommended volumes according to vessel format

- 4 Keep prepared destination vessel at 37 °C until usage.

hPSC non-enzymatic passaging

51m

- 5 Prepare required volume of the reagents and media according to the **Table 1**. Equilibrate the media to Room temperature .

30m

Note

Culture media aliquot to be used can be warmed up at 37°C for 00:15:00 . However, to preserve recombinant proteins activity Room temperature equilibration is recommended.



- 6 Aspirate and discard media from selected source vessel. 1m
- 7 Wash the well once using DPBS (no Calcium/no Magnesium). 2m
- 8 Add enzyme-free reagent, refer to **Table 1** for recommended working volumes. 1m
- 9 Incubate the vessel at 37 °C monitoring hPSC detachment using a microscope. Refer to **Table 2** for recommended incubation times. 5m

Cells will be ready to harvest when colonies start to be loosen though still attached.

A	B	C	D	E
Reagent	EDTA 0.5mM	GDR	ReLeSR	Versene
Time (min)	2-3	2-3	5-7	2-5

Table 2. Recommended incubation times according to enzyme-free reagent.


Note

This incubation can also be done at Room temperature , recommendation at 37 °C is to avoid hPSC stress.

- 10 Aspirate and discard enzyme-free reagent. 1m
- 11 Add fresh hPSC maintenance media to the wells, refer to **Table 1** for recommended harvesting volumes. 2m
- 12 Pipette up and down using 5 mL serological pipette against the bottom of the well to dislodge cells from the culture surface. 3m

**Note**

Repeated execution of this step leads to progressively smaller hPSC aggregates, which can increase cellular stress and promote cell death. Aim to maintain a balance that results in a uniform aggregate suspension - ideally with sizes ranging from approximately

 50-200 μm while minimizing shear stress.

Note

The use of scrapers is generally not recommended. However, for hPSC cultures with strong adherence properties, employing a scraper may be preferable to prolonged incubation with enzyme-free dissociation reagents.

- 13 Collect hPSC aggregates suspension using 5 mL pipette and directly distribute drop-wise the approximate volume for splitting ratios: 1/10, 1/15, 1/20 (or as desired) to the destination vessel/s.

5m


Note

Alternatively, aggregates suspension can be transferred to a conical tube prior distribution to better assess suspension volume.

- 14 Gently move destination vessels with freshly passed aggregates in cross shape to distribute them evenly.

1m

- 15 Incubate the vessels at  37 °C and at  5 % volume CO₂.

- 16 After  24:00:00 perform full media change. For further hPSC culture refer to protocol **Maintenance of hPSC** and protocol **Reference pictures of hPSC cultured in defined conditions**