



# A Payed Patel 1 Path or

Payal Patel<sup>1</sup>, Bethany Rozeboom<sup>1</sup>, Tarah-Anne Abrigo<sup>1</sup>

<sup>1</sup>Mercer University, Biochemistry Senior Capstone

Western Blot V.1

2 Works for me dx.doi.org/10.17504/protocols.io.bfiijkce

Apr 29, 2020 Bethany Rozeboom

#### MATERIALS

NAME	CATALOG #	VENDOR
2x Laemmli Sample Buffer	1610737	BioRad Sciences
TE buffer		Thermo Fisher Scientific
Bio Rad Precision Plus protein Standard	161-0374	Bio-rad Laboratories
Tryptic Soy Broth, E - 10mL	R07226	Thermo Fisher
SignalFire™ ECL Reagent	#6883S	Cell Signaling Technology
Tris Buffered Saline (TBS-10X)	#12498	Cell Signaling Technology
Tris-Glycine Transfer Buffer (10X)	#12539	Cell Signaling Technology
Tris-Glycine SDS Running Buffer (10X)	#4050	Cell Signaling Technology
Tris Buffered Saline with Tween® 20 (TBST-10X)	#9997	Cell Signaling Technology
Nonfat Dry Milk	#9999	Cell Signaling Technology
BSA	#9998	Cell Signaling Technology
Anti-mouse IgG HRP-linked Antibody	#7076	Cell Signaling Technology

### Cell Preparation

- 1 Grow cell **□50 ml** of TSB, let incubate **⊙ Overnight**
- 2 Distributed □17 ml of culture into two □50 ml centrifuge tubes, and □10 ml into one □20 ml centrifuge tube
- 3 Centrifuge down for **© 00:20:00** at 4300rpm ( 4303 xg) **§ 4 °C**
- 4 Aspirate media from cultures; wash cells 3x with TE aspirate.
- 5 Lyse cells by resuspending in **1 ml** of the different lysis buffers
  - 5.1 [M]20 Milimolar (mM) Tris CI , [M]15 Milimolar (mM) NaCI , 0.9696% Triton pH7.5 (labelled A)

    Phosphatase Inhibitor: [M]50 Milimolar (mM) NaF ,

[M]1 Milimolar (mM) Sodium Orthovanadate ,

[M]10 Milimolar (mM) Sodium Pyrophosphate and

Citation: Payal Patel, Bethany Rozeboom, Tarah-Anne Abrigo (04/29/2020). Western BlotÃÂ. https://dx.doi.org/10.17504/protocols.io.bfiijkce

```
[M]10 Milimolar (mM) B-glycerophosphate
                    EDTA-free Protease tablet (1 tablet per 10 ml)
             5.2 1% SDS, \blacksquare10 ml Tris Cl , [M]1 Milimolar (mM) EDTA (labelled B)
                    Phosphatase Inhibitor: [M] 50 Milimolar (mM) NaF,
                    [M]1 Milimolar (mM) Sodium Orthovanadate,
                    [M] 10 Milimolar (mM) Sodium Pyrophosphate and
                    [M]10 Milimolar (mM) B-glycerophosphate
                    EDTA-free Protease tablet (1 tablet per 10 ml)
     Freeze thaw 3 times
     Sonicated for © 00:01:00, repeat 5 times
     Centrifuged at 13000 rpm for © 00:05:00
     Bradford Assay
     Combine 40 µl 4X Laemmli to cell lysate, 40 µl of Lysis A and Lysis B
     Load ■40 µl onto SDS-PAGE gel and ■40 µl of All Blue Precision Plus for ladder
     Soaked sponges and transfer paper in Transfer Buffer
12
13
     Soak PVDF membrane in Methanol
     Soak fiber pads thoroughly in a Transfer Buffer
14
15
     Make the blotting sandwich
```

- 15.1 Add 1 cm depth of Transfer Buffer to container and insert plastic cassette with black side down, lay a wet fiber pad on the black side of the cassette
- 15.2 Lay one wet blotting paper on the fiber pad and roll out air bubbles
- 15.3 Lay gel squarely on blotting paper and roll out air bubbles
- 15.4 Lay wet PVDF on the gel and roll out air bubbles
- 15.5 Lay one wet blotting paper on the membrane and roll out air bubbles
- 15.6 Lay a wet fiber pad on top of the blotting paper and close the cassette and clamp together with the white clip.
- 16 Electrotransfer to PVDF membrane at 30 mV **Overnight**

## Membrane Blocking

- 17 After transfer, wash PVDF membrane with **25 ml TBS** for **00:05:00** at **8 Room temperature**
- 18 Incubate membrane in **□25 ml Blocking Buffer** for **⊙01:00:00** at **§ Room temperature** 
  - 18.1 Combine **□7.5** g Nonfat dry Milk with **□150** ml TBST to make Blocking Buffer
- 19 Wash three times for **⊙ 00:05:00** each with **□15 ml TBST**
- 20 Incubate membrane in Primary Antibody (at 1:2000) in a **10 ml Primary Antibody Dilution Buffer** with gentle agitation **Overnight** 
  - 20.1 Combine **□1** g 5% BSA to **□20** ml 1X TBST and mix well

- 21 Wash three times for **⊘00:05:00** each with **□15 ml TBST**
- 22 Incubate membrane with Anti-mouse IgG, HRP-linked Antibody (1:1000) to detect biotinylated protein markers in 

  10 ml Blocking Buffer with gentle agitation for © 01:00:00 at & Room temperature
- 23 Wash three times for **⊙ 00:05:00** each with **□15 ml TBST**

#### Detection of Proteins

- 24 Prepare 1X SignalFire™ ECL Reagent (#6883) by diluting one part 2X Reagent A and one part 2X Reagent B (based off of the size of our membrane we mixed ☐1 ml of each together)
- 25 Incubate substrate with membrane for **© 00:07:00**, remove excess solution (membrane remains wet) expose to X-ray film using an imager