



APR 05, 2024


Ventral Midbrain Genomic PCR

madalynn.erb Erb¹

¹Van Andel Research Institute

ASAP Collaborative Research Network

Team Lee

 Jane Balster
ASAP - Team Lee

ABSTRACT

This protocol details ventral midbrain genomic PCR.

MATERIALS

- ⊗ Stainless Steel Brain Matrices, 1.0mm **Stoelting Catalog #51386**
- ⊗ QIAgen DNeasy Blood and Tissue Kit, 50 rxn **Qiagen Catalog #69504**
- ⊗ GeneJET PCR Purification Kit **Thermo Fisher Scientific Catalog #K0702**

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.kqdg324zzv25/v1

Protocol Citation: madalynn.erb Erb 2024. Ventral Midbrain Genomic PCR. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kqdg324zzv25/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Feb 12, 2024



A	B	C
Reagent	Final conc.	ul/sample
5X KAPA2G Buffer A	1.3X	6.5 µl
25 mM MgCl2	2.60 mM	2.6 µl
10 mM KAPA dNTP Mix	0.26 mM	0.65 µl
Forward primer (10 µM) CTGCAGCTTCGAGAGGAAAG	0.5 µM	0.5 µl
Flox reverse primer (10 µM) CACTCTGTCCTCAGGCTTTC	0.5 µM	0.5 µl
KO reverse primer (10 µM) AGGTGGGAATCGGGCTAGAG	0.5 µM	0.5 µl
50% Glycerol	6.50 %	3.25 µl
5 U/ µl KAPA2G Fast Hotstart DNA Polymerase	0.5 U/ul	0.1 µl
DNA (diluted to 25 ng/uL)	150 ng total	6.0 µl
H2O	to 25 ul	4.4 µl

Last Modified: Apr 05, 2024

PROTOCOL integer ID: 97844


Keywords: ASAPCRN

Brain tissue collection


- 1 Euthanize mouse via cervical dislocation.
- 2 Isolate  coronal midbrain section using stainless steel brain matrix (Stoelting 51386).
 - 2.1 Remove cortex and dorsal midbrain tissue.
 - 2.2 Separate ipsilateral and contralateral ventral midbrain regions.
 - 2.3 Immediately freeze on dry ice and store tissue at  -80 °C until DNA extraction.

DNA extraction from frozen brain tissue

- 3 Use Qiagen DNeasy Blood and Tissue Kit (Cat 69504).

4 Equilibrate ventral midbrain tissue to  Room temperature .

Note

Starting material amount ~  15 mg of tissue.

5 Cut tissue into small pieces.

5.1 Use a pipette tip to transfer tissue to a sterile 6cm dish.



5.2 Add  180 μ L of **Buffer ATL** to the tissue.



5.3 Dice up tissue with sterile scalpel.

5.4 Transfer tissue and buffer to 1.5mL etube.



6 Add  20 μ L **Proteinase K** to each sample.



6.1 Mix thoroughly by vortexing.



6.2 Incubate at 56°C until samples are completely lysed.



Note

Use thermomixer – check samples after $01:00:00$, may take up to $03:00:00$.

7 Add $4\ \mu\text{L}$ **RNase A** (100 undetermined) to each sample.



7.1 Mix by vortexing.



7.2 Incubate at Room temperature for $00:02:00$.



2m

8 Mix **Buffer AL** with ethanol 1:1.



8.1 Add $400\ \mu\text{L}$ Buffer AL / ethanol mix to each sample.



8.2 Vortex for 00:00:15 .

15s

9 Transfer samples to DNeasy Mini spin columns (placed in 2mL collection tubes).



9.1 Centrifuge at ≥ 6000 x g for 00:01:00 .

1m



9.2 Discard flow through.

10 Transfer column to new collection tube.



10.1 Add 500 µL Buffer AW1.



10.2 Centrifuge at ≥ 6000 x g for 00:01:00 .

1m



10.3 Discard flow through and collection tube.


11 Place the column in a new collection tube – add  500 μ L **Buffer AW2**.




11.1 Centrifuge  20000 x g for  00:03:00 to dry the column. 3m





11.2 Carefully remove the column to avoid contamination with residual ethanol.

11.3 If the column touches ethanol flow through, spin again in new collection tube for  00:01:00 1m




12 Place the column in clean 1.5mL etube – add  150 μ L of elution buffer (**Buffer AE**).



12.1 Incubate at  Room temperature for  00:01:00 . 1m



12.2 Centrifuge at \geq  6000 x g for  00:01:00 to elute DNA. 1m



Genomic PCR

13 DNA concentration measured using a NanoDrop One Spectrophotometer (Thermofisher Scientific).

13.1 All samples diluted to 25 undetermined with Buffer AE.



14 Use the Kapa2g Fast HotStart PCR Kit (Roche 07960930001) according to the manufacturer’s instructions.



A	B	C
Reagent	Final conc.	ul/sample
5X KAPA2G Buffer A	1.3X	6.5 µl
25 mM MgCl2	2.60 mM	2.6 µl
10 mM KAPA dNTP Mix	0.26 mM	0.65 µl
Forward primer (10µM)	0.5 µM	0.5 µl
CTGCAGCTTCGAGAGGAAAG		
Flox reverse primer (10µM)	0.5 µM	0.5 µl
CACTCTGTCCTCAGGCTTTC		
KO reverse primer (10µM)	0.5 µM	0.5 µl
AGGTGGGAATCGGGCTAGAG		
50% Glycerol	6.50%	3.25 µl
5 U/ µl KAPA2G Fast Hotstart DNA Polymerase	0.5 U/ul	0.1 µl
DNA (diluted to 25ng/uL)	150ng total	6.0 µl
H2O	to 25 ul	4.4 µl

PCR Cycle:

A	B	C	D
Step	Temp (°C)	Time	Note
1	94.0 °C	5 min.	
2	94.0 °C	30 sec.	
3	65.0 °C	15 sec.	-0.5 °C per cycle decrease
4	68.0 °C	1 sec.	
5			repeat steps 2-4 for 10 cycles (touchdown)

	A	B	C	D
	6	94.0 °C	30 sec.	
	7	60.0 °C	15 sec.	
	8	72.0 °C	1 sec.	
	9			repeat steps 6-8 for 20 cycles
	10	72.0 °C	5 min.	
	11	4.0 °C	hold	Hold

15 Run samples on 2% Agarose gel at 100V for 35-45 minutes.



A	B
Wild type allele	400 bp
Flox allele	500 bp
KO allele	270 bp