



Nov 20, 2021

Basic Rolling Circle Amplification (RCA) Protocol (Thermo phi29 polymerase)

argawikandito¹¹Universitas Gadjah Mada

1

dx.doi.org/10.17504/protocols.io.bz8np9ve

Laboratory of Virology, Faculty of Agriculture, Universitas Gadjah Mada

argawikandito

DOI

dx.doi.org/10.17504/protocols.io.bz8np9ve

argawikandito 2021. Basic Rolling Circle Amplification (RCA) Protocol (Thermo phi29 polymerase). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bz8np9ve>



document ,

Nov 20, 2021

Nov 20, 2021

55278

Basic Rolling Circle Amplification (RCA) Protocol (Thermo phi29 polymerase)

Aim: RCA is useful method to amplify circular DNA/RNA. By principle, RCA conducted using phi29 polymerase to obtain copies of DNA/RNA target in concatemeric form. phi29 works on isothermal reaction at 36-38°C.

Materials :

gDNA

10U phi29 polymerase (Thermo Scientific)

10X reaction buffer (Thermo Scientific)
 10mM dNTP mix (Thermo Scientific)
 100μM n-hexamer primer (or specific primers to obtain specific target) (Thermo Scientific)
 Nuclease-free water (Thermo Scientific)
 Thermal cycler
 Micropipettes
 Tubes 0.2μL

Steps :

DNA Denaturation

- Prep denaturation mixture on 0.2μL tube :

Reagent	Volume (μL)
10X reaction buffer	0.5
100μM n-hexamer primer	1
Nuclease-free water	2-3
gDNA	0.5-1
Total	5
NOTE : you can change 100μM n-hexamer primer with specific primer to obtain more specific result	

- Incubate 95°C for 3 minutes
- Put on ice immediately for 3-5 minutes



Rolling Circle Amplification

- Prep amplification mixture

Reagent	Volume (μL)
Denaturation mixture	5
10X reaction buffer	1.5
10mM dNTP mix	2
10U phi29 polymerase	1
Nuclease-free water	10.5
Total	20
NOTE : you can also add 100μM n-hexamer primer to increase RCA yield	

- Incubate 36°C for 18 hours
- Inactivate reaction at 65°C for 10 minutes

