

Apr 19, 2024

ImmPRESS Polymer Detection Following BaseScope In Situ Hybridization

DOI

dx.doi.org/10.17504/protocols.io.kqdg325rzv25/v1

madalynn.erb Erb¹

¹Van Andel Research Institute

ASAP Collaborative Rese...



madalynn.erb Erb
Van Andel Research Institute

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.kqdg325rzv25/v1

Protocol Citation: madalynn.erb Erb 2024. ImmPRESS Polymer Detection Following BaseScope In Situ Hybridization. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kqdg325rzv25/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 19, 2024

Last Modified: April 19, 2024

Protocol Integer ID: 98492

Keywords: ASAPCRN

Funders Acknowledgement:

Madalynn Erb

Grant ID: ASAP-000592



Abstract

Protocol for performing ImmPRESS immunostaining with alkaline phosphatase labeling after BaseScope In Situ Hybridization mouse brain sections.

Materials

Reagents: ImmPRESS®-AP Horse Anti-Rabbit IgG Polymer Detection Kit (MP-5401; Vector Labs) Vector® Blue

Alkaline Phosphatase Substrate (SK-5300; Vector Labs)



Day 1

- 1 Complete BaseScope in situ hybridization using 35um mouse brain sections. Do not dry slides. For BaseScope protocol see : (dx.doi.org/10.17504/protocols.io.5qpvo364zv4o/v1)
- 2 Wash slides 2 times in PBS at room temp 5 minutes each wash
- 3 Block sections in 0.4% Bovine Serum Albumin (BSA), 10% horse serum, 0.3% triton X-100 for 1 hour at room temp
- 4 Incubate sections in rabbit anti TH (1:500) in 0.4% BSA, 5% horse serum and 0.1% tritonX-100 overnight at 4 degrees C

Day 2

- 5 Wash slides in PBS at room temp (3 times 5 minutes each wash)
- 6 Incubate slides in secondary antibody (ImmPRESS-AP Polymer Reagent) at room temp for 2 hours.
- 7 Wash slides in PBS at room temp (3 times 5 minutes each wash)
- 8
 1. Prepare the Alkaline Phosphatase Substrate
 - a. Make 2.5mL 100mM Tris-HCl, pH 8.2
 - b. Add 1 drop of Reagent 1
 - c. Add 1 drop of Reagent 2
 - d. Add 1 drop of Reagent 3
 - e. Mix well – use immediately
- 9 Incubate sections in alkaline phosphatase substrate for 5 minutes
- 10 Wash slides in PBS (2 times 5 min each wash)



- 11 Dry slides at 60°C for ≥ 15 min until they are completely dry
- 12 Coverslip slides using VectaMount
- 13 Store slides at room temp

Protocol references

Preparation of Free Floating Coronal Mouse Brain Sections : dx.doi.org/10.17504/protocols.io.dm6gp3pb8vzp/v1

BaseScope In Situ Hybridization: (dx.doi.org/10.17504/protocols.io.5qpvo364zv4o/v1)