




MAR 07, 2024

Pulldowns

 In 1 collection

Dan Tudorica¹

¹University of California, Berkeley



Dan Tudorica
Hurley Lab, QB3, UC Berkeley

ABSTRACT

Traditional in vitro pulldown

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.5jyl8pw6dg2w/v1

Protocol Citation: Dan Tudorica 2024. Pulldowns . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5jyl8pw6dg2w/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Nov 09, 2023

Last Modified: Mar 07, 2024

Pulldown

- 1 In a final volume of 30 μ L, 20 μ g of purified MBP-Rubicon RH domain was mixed with 20 μ g of RAB7A in 50 mM HEPES 7.5, 150 mM NaCl, 2 mM $MgCl_2$, 10 mM TCEP buffer.

Include a negative control consisting of only soluble Rab7, with no MBP-Rubicon, in order to test your washing efficacy.
- 2 Incubate samples on rocker at room temperature for 1 H
- 3 Add 20 μ L of 50% v/v amylose resin to each sample
- 4 Allow to rock for an additional 30 min to bind
- 5 Collect resin at bottom of tube via a tabletop centrifuge, aspirate off and discard supernatant, and wash with 500 μ L of ice cold buffer for three total washes.
- 6 Elute sample by resuspending beads in 20 μ L of buffer + 20 mM maltose for 10 min on shaker.
- 7 Prepare samples for SDS-PAGE by adding 1x loading buffer, and load onto a 4-12% gel. Run gel at 120 V until dye front reaches the bottom of the gel

- 8 Stain with Coomassie blue G-250, and destain with water until bands are clearly visible. Image.