



Mar 13, 2021

# Peptostreptococcal protein L and a Chimeric Protein-LAG sandwich ELISA

Angel A Justiz-Vaillant<sup>1</sup><sup>1</sup>University of the West Indies St. Augustine**1** Works for me [dx.doi.org/10.17504/protocols.io.bta6nihe](https://dx.doi.org/10.17504/protocols.io.bta6nihe)

Carbon

Angel Justiz-Vaillant  
University of the West Indies St. Augustine

SUBMIT TO PLOS ONE

## ABSTRACT

This ELISA was used to study the interactions between Peptostreptococcal protein-L (SpL) and protein LAG (PLAG) with different immunoglobulin preparations from mammalian and avian species.

## DOI

[dx.doi.org/10.17504/protocols.io.bta6nihe](https://dx.doi.org/10.17504/protocols.io.bta6nihe)

## PROTOCOL CITATION

Angel A Justiz-Vaillant 2021 . Peptostreptococcal protein L and a Chimeric Protein-LAG sandwich ELISA.  
**protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bta6nihe>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Mar 13, 2021

## LAST MODIFIED

Mar 13, 2021

## PROTOCOL INTEGER ID

48190

- 1 This ELISA was used to study the interactions between Peptostreptococcal protein-L (SpL) and protein LAG (PLAG) with different immunoglobulin preparations from mammalian and avian species. The 96 well microtiter plate was coated overnight at 4°C with 2 µg/µl per well of SpL in carbonate-bicarbonate buffer pH 9.6.
- 2 The plate was then treated with bovine serum albumin solution and washed 4X with PBS-Tween. 50 µl of immunoglobulins (1 mg/ml) is added and incubated for 1h at room temperature, and then, the microplate is rewashed 4X with PBS-Tween.
- 3 Then, 50 µL of peroxidase-labeled PLAG conjugate diluted 1:5000 in PBS-non-fat milk was added to each well and

incubated for 1 h at RT. The plate was washed 4X with PBS-Tween.

- 4 Then, 50  $\mu$ L of o-phenylenediamine solution (4 mg/mL) was added, and the plate was incubated for 15 min at RT in the dark. The reaction was stopped with 50  $\mu$ L of a 3M H<sub>2</sub>SO<sub>4</sub> solution.
- 5 The plate was visually assessed for color development and read on a microplate reader at 492 nm. A cut-off point was calculated as the mean of the optical density of the negative controls multiplied by two. The cut-off point was set to 0.28.