



Apr 05, 2020

Pillowcase & rubber band DIY face mask @

Jernej Turnsek¹

¹J. Craig Venter Institute (JCVI)



dx.doi.org/10.17504/protocols.io.bemjjc4n

Coronavirus Method Development Community









ABSTRACT

Wearing a face mask helps stop the spread of the new coronavirus. This protocol describes a very simple way to make one at home using a cotton pillowcase, 2 rubber bands, and a minute of your time.

EXTERNAL LINK

https://www.cdc.gov/coronavirus/2019-ncov/prevent-getting-sick/diy-cloth-face-coverings.html

BEFORE STARTING

Cut a cotton pillowcase into 4 ~equally sized pieces.

Start with a \sim 14 inch x \sim 20 inch (\sim 36 cm x \sim 51 cm) pillowcase piece and 2 rubber bands.



2 Fold the fabric twice along the longer edge.



3 Insert both ends of the folded piece through a rubber band.



2 long twisted rubber bands as shown below were used above.



4 Fold both outer "flaps" back to the middle creating a total of 9 fabric protection layers.



Secure the mask behind your ears.

5 Wash the fabric by boiling it in salty water for 5 min after each use. Dry it thoroughly before using it again.

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited