



Feb 25, 2020

## 'Uniclear' water-based brain clearing for light sheet imaging of electrode tracks

Susu Chen<sup>1</sup>, Karel Svoboda<sup>1</sup><sup>1</sup>Janelia Research Campus

1 Works for me

[dx.doi.org/10.17504/protocols.io.zndf5a6](https://doi.org/10.17504/protocols.io.zndf5a6)
**Susu Chen**  
 Janelia Research Campus 

## ABSTRACT

UniClear procedure for whole mouse brain clearing and refractive index matching. The advantages of this method are that the method is water-based and produces a mechanically robust specimen. Cleared brains can be imaged with lightsheet microscopes (e.g. Zeiss Z1 lightsheet).

## MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
DMSO	472301	Sigma Aldrich
OptiPrep™ Density Gradient Medium	D1556)	Sigma Aldrich
2-Propanol (99.5 %)	278475	Millipore Sigma

## MATERIALS TEXT

2-methyl-2-butanol

<https://www.sigmaaldrich.com/catalog/product/sial/152463?lang=en@ion=US>

2-propanol

<https://www.sigmaaldrich.com/catalog/product/sial/278475?lang=en@ion=US>
Sodium dodecyl sulfate solution (SDS), 10% in H<sub>2</sub>O
<https://www.sigmaaldrich.com/catalog/product/sigma/71736?lang=en@ion=US>

OptiPrep™ Density Gradient Medium

<https://www.sigmaaldrich.com/catalog/product/sigma/d1556?lang=en@ion=US>

Dimethyl sulfoxide (DMSO)

<https://www.sigmaaldrich.com/catalog/product/sial/276855?lang=en@ion=US>

Nycodenz

<https://www.progen.com/nycodenz.html>

Preparation of SBiP buffer: 500ml

\*\*\* Mix SBiP at 4C till fully dissolved, kept at 4C.

At room temperature, SBiP will get activated and emulsified for delipidation. Use each batch within a month for best effect.

- Ice cold H<sub>2</sub>O 350 ml
- 50mM Na<sub>2</sub>HPO<sub>4</sub> 2 ml
- 4% SDS (in H<sub>2</sub>O, pH7.4) 10 ml
- 2-methyl-2-butanol 80 ml
- 2-propanol 40 ml

## Preparation of Opti-prep refractive index matching solution

- 60% (w%) Opti-prep solution
- 20% (w%) DMSO
- 20% (w%) PBS
- 20 mM Tris
- add overall 1.12g/ml Nycodenz (dissolve gradually by adding ~1/10th of Nycodenz powder at a time until fully mixed, heating up and stirring can help to dissolve faster)

## References:


Economo, Clark, et al., A platform for brain-wide imaging and reconstruction of individual neurons.


[Elife](#). 2016 Jan 20;5:e10566. doi: 10.7554/eLife.10566

Winnubst et al., Reconstruction of 1,000 Projection Neurons Reveals New Cell Types and Organization of Long-Range Connectivity in the Mouse Brain

[Cell](#). 2019 Sep 19;179(1):268-281.e13. doi: 10.1016/j.cell.2019.07.042. Epub 2019 Sep 5.

- 1 Perfuse mouse with cold 4% PFA. Post fix brain for at least ⌚ **24:00:00** at room temperature (RT). Wash brain in PBS to remove fixative. Fixed brain samples can be kept in PBS at 🌡 **4 °C** before further processing.
- 2 Delipidate brain with SBiP buffer: 📄 **20 ml** per adult mouse brain, shaking (70 RPM) at 🌡 **37 °C** . On the first day: change buffer at ⌚ **03:00:00** , ⌚ **06:00:00** , ⌚ **12:00:00** and leave over night. Starting on the second day: change SBiP buffer every ⌚ **24:00:00** for consecutive 14 days, shaking (70 RPM) at 🌡 **37 °C** .
- 3 After delipidation, wash brain in PBS for at least ⌚ **12:00:00** before further processing. Sample can be stored in PBS at RT before refractive index matching.
- 4 Refractive index matching: first day in 📄 **20 ml** 50% PBS + 50% opti-prep mixed solution, rotating at RT for ⌚ **24:00:00** . On second day transfer the brain to 📄 **40 ml** 100% opti-prep solution, rotating at RT for ⌚ **24:00:00** . Index-matched brain should appear transparent.
- 5 For light-sheet imaging (e.g. Zeiss Z1) mount brain vertically.





Example of a mounted brain (fluorescent yellow color is DNA dye YO-PRO, which is optional). For electrode localization we use autofluorescence.

To image tissue autofluorescence use 488 nm laser. To image CM-Dil tracks use 561 nm laser. Brain can be imaged with a voxel size of  $1.22\ \mu\text{m} \times 1.22\ \mu\text{m} \times 8\ \mu\text{m}$  on Z1 with a 5x objective.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited