

FEB 27, 2023

OPEN BACCESS

Protocol Citation: Santosh Sathe 2023. Evans blue assay to stain dead cells.. **protocols.io**

https://protocols.io/view/evans-blue-assay-to-stain-dead-cells-cp3wvqpe

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: Feb 27, 2023

Last Modified: Feb 27, 2023

PROTOCOL integer ID:

77654

Evans blue assay to stain dead cells.

Santosh Sathe¹

¹universitat konstanz



Santosh Sathe

ABSTRACT

Evans blue (Sigma, USA, catalog no: E2129) is a dye that can enter the compromised cell membranes of dead cells and stain them blue. I used this assay to differentiate dead cells in algal cultures of *Chlorella* and *Chlamydomonas* (the same protocol can be used for other green algae).

Evans blue assay to stain dead cells.

1 Contact: Santosh Sathe (santosh.sathe@uni-konstanz.de)

Introduction

Evans blue (Sigma, USA, catalog no: E2129) is a dye that can enter the compromised cell membranes of dead cells and stain them blue. I used this assay to differentiate dead cells in algal cultures of *Chlorella* and *Chlamydomonas* (the same protocol can be used for other green algae).

Preparation

Prepare 1% Evans blue (w/v in MilliQ water) by dissolving 0.1 gm Evans blue in 10 milliliter MilliQ water (I use Falcon tubes). Filter the dye through 0.2 Micron filters using a syringe and filter assembly (bigger pore size filters can also be used). I filter the dye to remove large undissolved particles of dye. The solution (cover it with aluminium foil) is stable at room temperate for several months.

Assay

- Take 900 μl algal culture in 1.5 milliliter Eppendorf tube (I avoid 2 milliliter Eppendorf tubes as in these tubes the cells do not pellet properly; see below). If the algal cultures are less dense (less than a million cells/ milliliter), then a few milliliter of culture can be centrifuged, and the cell pellet combined and suspended in 900 μl culture media.
- 5 Add 100 ul Evans blue (from 1% dye suspension prepared above).
- 6 Incubate the mixture at room temperature for 20 minutes, preferably in the dark.
- Rinse the cells with MilliQ water to remove excess dye. First, centrifuge the stained cells at 7000 RPM for 2 minutes. Remove about 900 µl dye with a 1 milliliter pipette (do not disturb the cell pellet). Since the dye is blue, the cell pellet cannot be seen during the first wash. Add 1 milliliter MilliQ water and centrifuge the cells. The cell pellet should be visible now as the leftover dye is no more concentrated. Repeat the above steps 2 to 3 times (remove 900 to 950 µl dye from the supernatant using a pipette, add MilliQ water, and centrifuge). Leave the cell pellet (after final wash) in a small volume (10 to 20 µl). Take a small amount of cell pellet (ca. 5 µl) and mount it on a clean glass slide. Put a coverslip, press gently so that cells are dispersed and observe the cells using a microscope. I use an oil immersion lens (100 X; Axioskop 40, Carl Zeiss, Germany). Count blue (dead) cells and estimate the dead cells (%) from the total cells.