

NOV 11, 2022



**WORKS FOR ME** 

Ethanol Quantification assay

In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.bnuemete

Marcos Valenzuela-Ortega<sup>1</sup>

<sup>1</sup>University of Edinburgh



Marcos Valenzuela-Ortega University of Edinburgh

COMMENTS 0

**ABSTRACT** 

Chemoenzymatic method for quantification of ethanol with a spectrophotometer at 500 nm (not UV!)

Based on **Lewicka 2014** (DOI: 10.1021/sb500020g)

Adapted for analysis in plate reader of multiple samples at the same time.

dx.doi.org/10.17504/protocols.io.bnuemete

PROTOCOL CITATION

Marcos Valenzuela-Ortega 2022. Ethanol Quantification assay. protocols.io https://dx.doi.org/10.17504/protocols.io.bnuemete

**COLLECTIONS** ①



untitled collection



Citation: Marcos Valenzuela-Ortega Ethanol Quantification assay <a href="https://dx.doi.org/10.17504/protocols.io.bnuemete">https://dx.doi.org/10.17504/protocols.io.bnuemete</a>

#### LICENSE

This is an open access protocol distributed under the terms of the <u>Creative Commons</u>

<u>Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**CREATED** 

Oct 23, 2020

LAST MODIFIED

Nov 11, 2022

PROTOCOL INTEGER ID

43622

PARENT PROTOCOLS

Part of collection

### untitled collection

MATERIALS TEXT

Yeast Alcohol Dehydrogenase - Sigma...

NAD+ - Sigma

Phenazine methosulfate (PMS) - Sigma

Iodonitro-tetrazolium violet (INTV) -Sigma

50 mM Tris-HCl buffer (pH 7.5)

50 mM tris-HCl pH= 7.5 (0.118g tris base + 0.635g tris.HCl in 100 ml of H20)

(or buy pre-made tris pH 7.5 and dilute it accordingly)

Ethanol absolute to be used as standard.

Multichannel pipette

Flat-bottom 96-well plate

### **Prepare culture samples**

Sample your culture and spin it 3 16000 rpm, 00:01:00

#### Prepare stocks

2

Prepare the following stocks using your tris buffer:

m protocols.io

2

Citation: Marcos Valenzuela-Ortega Ethanol Quantification assay <a href="https://dx.doi.org/10.17504/protocols.io.bnuemete">https://dx.doi.org/10.17504/protocols.io.bnuemete</a>

Stock	mg/ml	mM	Mw
Yeast	5		
NAD+	13.3	20	685.41
PMS	6.1	20	306.34
INTV	25.3	50	505.00
Ethano	46 (36.	1000	46

INTV needs to be dissolved in buffer-DMSO (50:50)

Keep reagents On ice (no need for ethanol)

Keep the buffer at 8 Room temperature

Prepare your ethanol standard ladder by making different dilutions of your 1M ethanol with tris buffer. The standard could be, for example, 400-200-100-50-25-10-5-0 mM ethanol.

# **Prepare reaction mastermix**

4 Prepare reaction mix for all your samples.

Calculate the volume to be prepared should be **200 \muL \* well \* 1.1** (so you have) When counting "wells", consider replicates for samples and standards.

That is the volum of Tris you need, to which you have to add:

Yeast ADH: 5 µL per mL
NAD+: 1µL per mL
PMS: 1 µL per mL
INTV: 20 µL per mL

## **Reaction plate and incubation**

5 Plan how to load samples in plate, and annotate it.

Consider that the reaction mix should be loaded to the replicates at different moments (so their variability



3

Citation: Marcos Valenzuela-Ortega Ethanol Quantification assay <a href="https://dx.doi.org/10.17504/protocols.io.bnuemete">https://dx.doi.org/10.17504/protocols.io.bnuemete</a>

comprises the variability caused by loading the mix at the different times.) 5.1 Load 200 µL of buffer to the blank wells. (some replicates here are good too). These wells wont receive reaction mix. 5.2 Load 10 µL of sample to the well bottom. Be sure all the volume of liquid is in the well and not in your pipette 5.3 Load 200 µL of reaction mix to the wells with your multichannel pipette with as little difference of time between sample as possible. 6 Mix & incubate 6.1 Mix the plate (the plate reader can do this). 6.2 Incubate the plate for \$\ \mathbb{R}\$ Room temperature for 15 minutes, in the dark.

7 Measure Abs 500 nm in plate reader at Room temperature

Note

Changes in temperature can cause condensation in the plate lead, which will affect readings.

### **Analysis and troubleshooting**

8 Blank should be done with water, samples measurements will be reliable only in the linear region of the standard.

BMG plate readers do most of the analysis for you, if you tell the software which wells had each blank,

protocols.io

4

Citation: Marcos Valenzuela-Ortega Ethanol Quantification assay https://dx.doi.org/10.17504/protocols.io.bnuemete

standarda and sample. You can introduce this in the software before or after the readings.

You can have a positive control (negative control culture spiked with ethanol).

Sensitivity can be increased using 20  $\mu$ L of sample in 200  $\mu$ L of mastermix, or increasing incubation time to 30 minutes.

9

