



Aug 19, 2021

Sequential extraction and immunoblotting

Isabel Lam¹¹Khurana Lab, Brigham and Womens Hospital, Harvard Medical School

Isabel Lam: Protocol from VolpicelliDaley et al Nat Protocols 2014.

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bu2wnyfe

hendersa

ABSTRACT

This protocol examines the fraction of alpha-synuclein (as assessed by alpha-synuclein and/or PS129 western blot) that is present in the triton-soluble or SDS-soluble fraction. Addition of alpha-synuclein pre-formed fibrils (PFFs) to neuron cultures seeds the recruitment of endogenous or transgenic alpha-synuclein into aggregates characterized by detergent insolubility. Monomeric alpha-synuclein will be present in the triton-soluble fraction, whereas PFF-induced alpha-synuclein oligomers will be present in the triton-insoluble, SDS-soluble fraction. This protocol encompasses preparation of cell extracts with Triton X-100, followed by sequential extraction of the Triton-insoluble material with SDS. The protein in the different detergent fractions are quantified by BCA, followed by SDS-PAGE and WB for alpha-synuclein, PS129, TUJ1, and loading control such as GAPDH.

ATTACHMENTS

[dh35bqa7.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.bu2wnyfe

PROTOCOL CITATION

Isabel Lam 2021. Sequential extraction and immunoblotting. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bu2wnyfe>



KEYWORDS

Immunoblotting, Alpha-synuclein, Sequential extraction

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

May 15, 2021

LAST MODIFIED

Aug 19, 2021

OWNERSHIP HISTORY

May 15, 2021 Urmilas


Jun 23, 2021 hendersa

PROTOCOL INTEGER ID

49974

MATERIALS TEXT

Materials:

-  **50 mL** Falcon tube
- Bucket with ice
- Cell scraper
- Sonicator (I used Misonix Sonicator S-4000 from Wolfe Lab, next to ultracentrifuge)
- Ultracentrifuge (I used Optima MAX-TL Ultracentrifuge TLA-120.2 S/N 12U1811)
- Polyallomar tube

Reagents:

- 1X PBS
- 1% (vol/vol) TX-100/TBS with protease and phosphatase inhibitors
- 2% SDS/TBS with protease and phosphatase inhibitors
- 4X Laemmli buffer
- TBS/5% (wt/vol) milk
- 10X TBS
- BCA assay reagents
- SDS-PAGE reagents




Sequential extraction and immunoblotting

7h 50m

1 

Rinse the neurons twice with PBS.

2 




Place the dish  **On ice** . By working one well at a time, completely aspirate the PBS and add the following volumes of ice-cold 1% (vol/vol) TX-100/TBS with protease and phosphatase inhibitors:  **250 µl** per well for a six-well plate, and  **500 µl** for a 6-cm dish.

3 Use a cell scraper to thoroughly scrape all neurons from each well.

4 Place the neurons in a polyallomar tube for a table top ultracentrifuge; keep the tube  **On ice** .

5 

Sonicate the tube ten times at a 0.5-s pulse and at 10% power.

5.1 Use Misonix Sonicator S-4000, with Program Settings: Amplitude 10, Process Time  **00:00:10** ,^{12s}
Pulse-ON time  **00:00:01** , Pulse-OFF time  **00:00:01** .

5.2 Wipe sonication tip with 1% SDS after coming into contact with PFFs, then 70% ethanol.

6 

30m


Incubate the tube ⚡ **On ice** for ⌚ **00:30:00**.

7 

30m

Centrifuge the tube at 🌀 **100000 x g** at ⚡ **4 °C** for ⌚ **00:30:00**.

7.1 Optima MAX-TL Ultracentrifuge TLA-120.2 S/N 12U1811 🌀 **100000 x g** ⌚ **00:30:00** ⚡ **4 °C**.^{30m}

8 

Add 4× Laemmli buffer to ~ 📏 **150 µl** – 📏 **200 µl** of TX-100 supernatant. Save ~ 📏 **20 µl** of supernatant for protein assay. Retain the supernatant ⚡ **On ice** or in a ⚡ **-20 °C** freezer.

Some of the supernatant is usually lost during Step 3–8; we recommend measuring the remaining supernatant before determining how much to remove for Laemmli buffer.

9 

To the pellet, add the same volume of ice-cold 1% (vol/vol) TX-100/TBS with protease and phosphatase inhibitors: 📏 **250 µl** per well for a six-well plate, and 📏 **500 µl** for a 📏 **6 cm** dish.

10 Sonicate ten times at a 0.5-s pulse and at 10% power. Keep the tip of the probe toward the bottom of the tube to prevent frothing. Make sure that the pellet is completely dispersed.

a. Sonicator: Fisherbrand™ Model 120 Sonic Dismembrator (120V, 50/60Hz)
b. I used Misonix Sonicator S-4000 with Program Settings: Amplitude 10, Process Time ⌚ **00:00:10**, Pulse-ON time: ⌚ **00:00:01**, Pulse-OFF time: ⌚ **00:00:01**

11 

30m

Centrifuge the mixture at 🌀 **100000 x g** at ⚡ **4 °C** for ⌚ **00:30:00**.

12 Discard the supernatant.

13 

Add 2% (vol/vol) SDS/TBS to the pellet with protease and phosphatases inhibitors. To a six-well plate, add **125 µl** of 2% (wt/vol) SDS/TBS per well, and to a **6 cm** dish add **250 µl** of 2% (wt/vol) SDS/TBS.

14 Sonicate 15 times, at a 0.5-s pulse and at 10% power. Keep the tip of the probe toward the bottom of the tube. Make sure that the pellet is completely dispersed.

a. I used Misonix Sonicator S-4000 with Program Settings: Amplitude 10, Process Time 00:00:15, Pulse-ON time: 00:00:01, Pulse-OFF time: 00:00:01

15 Remove the supernatant and place it into a new microcentrifuge tube.

16 Perform a BCA/protein assay on TX-100 supernatant and SDS extract. Typically, a 1:5 dilution for the BCA assay is sufficient.

17 Dilute 2% (wt/vol) SDS extract from Step 15 into Laemmli buffer to 2× volume for the corresponding TX-100 fraction (regardless of the protein concentration of the SDS fraction).

For example, if you have **180 µl** of TX-100 extract (from Step 8) at **1 mg/ml** and **90 µl** of SDS extract, add **60 µl** of 4× Laemmli buffer to the TX-100 extract and **30 µl** of 4× Laemmli buffer to the SDS extract. Load **12.5 µl** of both the TX-100 extract (**10 g**) and SDS extract. We suggest using a 2× volume because it makes the insoluble a-syn species more abundant and thus easier to visualize and quantify by immunoblotting.

18 Load the samples on a 4–20% (wt/vol) gel and run them according to the manufacturer's directions.

We use 85 V constant voltage until the dye front runs off the gel. Be sure not to let the 10-kDa marker to run off the gel.

2h 30m

19



Transfer the proteins from the gel to a nitrocellulose membrane according to the manufacturer's instructions at 100 V for ⌚ **01:15:00** or ⌚ **Overnight** at 40 V.

20

Block the membrane for ⌚ **01:00:00** with TBS/5% (wt/vol) milk.

1h

21



Dilute the primary antibodies in TBS/5% (wt/vol) milk and incubate them ⌚ **Overnight** at 🌡 **4 °C** with shaking.

1h

22



Rinse the membrane three times with TBS/T, ⌚ **00:10:00** each rinse.

10m

23



Incubate the membrane with HRP-conjugated secondary antibodies for ⌚ **01:00:00** at 🌡 **Room temperature**.

1h

24



Rinse the membrane three times with TBS/T, ⌚ **00:10:00** each rinse.

10m

25

Develop with enhanced chemiluminescence.