



Dec 21, 2020

# Effective Identification of Protein-Protein Interaction using RIME-IP

George Laliotis<sup>1</sup><sup>1</sup>The Ohio State University

George Laliotis: Corresponding author

**1** Works for me [dx.doi.org/10.17504/protocols.io.bqx8mxrw](https://dx.doi.org/10.17504/protocols.io.bqx8mxrw) George Laliotis

DOI

[dx.doi.org/10.17504/protocols.io.bqx8mxrw](https://dx.doi.org/10.17504/protocols.io.bqx8mxrw)

PROTOCOL CITATION

George Laliotis 2020. Effective Identification of Protein-Protein Interaction using RIME-IP. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bqx8mxrw>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 21, 2020

LAST MODIFIED

Dec 21, 2020

PROTOCOL INTEGER ID

45792

DISCLAIMER:

## DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## Buffer Preparation

2h

- 1 This is a modification of the Mohammed et al., 2016 (Nature protocols, 11(2), p.316.), that it has been edited and optimized for the need of Tschlis Lab. The authors declare no conflicts of interest.

## Immunoprecipitation

2d

- 2 Prepare the following buffers (**usually prepare 25ml of each, store at 4 °C And make fresh every month**) at the final concentration listed :

•1M Tris pH=7.5

for 100mL add 12.114g base  
Adjust pH to 7.5

•1M NaCl

For 100mL add 5.844g

•50mM EDTA

For 100mL add 1.4612g

•5mM EGTA

For 100mL add 0.19g

•50mM DTT

For 1,5mL add 75uL from 1M stock

#### Lysis Buffer 1 (LB1)-will contain Cytoplasmic Proteins

For 25ml Volume

- 50mM Tris-HCL ph 7.5 - add 1.25mL from 1M stock
- 20mM NaCl- add 500uL from 1M stock
- 1mM EDTA - add 500uL from 50mM stock
- 0.5% NP-40 - add 125uL from stock (Sigma # 74385-1L)
- 0.25% Triton X-100 - add 62.5uL from Stock
- 10% Glycerol - add 2.5 ml from stock
- 1mM DTT - add 75uL from 50mM stock

#### Lysis Buffer 2 (LB2)- wash to remove LB1 and cytosolic protein

For 25ml volume

- 10mM Tris-HCL ph 7.5 - add 250uL from 1M stock
- 20mM NaCl- add 500uL from 1M stock
- 1mM EDTA - add 500uL from 50mM stock
- 0.5mM EGTA - add 500uL from 5mM stock
- 1mM DTT - add 75mL from 50mM stock

#### Lysis Buffer 3 (LB3)

- 10mM Tris-HCL ph 7.5 - add 250uL from 1M stock
- 100mM NaCl- add 2.5ml from 1M stock
- 1mM EDTA - add 500uL from 50mM stock
- 0.5mM EGTA - add 500uL from 5mM stock
- 0.1% (w/v) sodium deoxycholate
- 0.5% (v/v) N-laurylsarcosine - add 625uL from 20% Sigma stock
- 1x protease inhibitors-add 250uL from 100x stocks

### 3 Day 1

#### Antibody-Bead Conjugation Day One

1. Magnetic Beads (Pierce Cat No. # 88802), lightly vortex to resuspend. With a cut pipette tip, transfer 50-75uL (some will be lost during the washes) of beads to 1mL Eppendorf tube. Place magnetic stand and remove supernatant.
2. Add 500uL ice-cold PBS and lightly vortex to resuspend. Spin briefly with table-top microcentrifuge at 2000xg for 1'. (Tip : by vortexing in higher speed you are destroying the magnetic beads). Place tube in magnetic stand and remove PBS wash. Repeat 2x for a total of 3 washes.
3. Add 300uL of LB3 to beads. Repeat vortex, spin and wash removal as in step 2 for a total of 3 washes.

4. Add antibody for IP (according to company specifications-5ug) to top of the beads and incubate on ice for 10 min. Add 300uL LB3 to beads. Vortex tube lightly to resuspend beads. Rotate bead/antibody mixture at  $4^{\circ}\text{C}$  O/N.

## 4 Day 2

### Harvest cells

1. Add media from 15cm plate at 80-90% confluent (around  $2 \times 10^7$  cells) to 50mL Falcon. Wash the cells with 5mL ice-cold PBS and add to Falcon. Add 3mL ice-cold PBS and collect cells with scraper. Add to Falcon. Rinse plate with 5-10mL ice-cold PBS and add to Falcon.
2. Pellet cells by centrifugation at 400xg for 6 min. Aspirate media
3. Resuspend cell pellet in 1mL of ice-cold PBS and transfer to Eppendorf tube. Rinse Falcon with 500uL of ice-cold PBS and transfer to Eppendorf. Spin as in step 2.
4. Remove PBS. Wash cell pellet with 500uL PBS, spin as before and remove PBS wash.
5. Resuspend cell pellet in 750uL of LB1+freshcocktail inhibitor. Rotate the tube at  $4^{\circ}\text{C}$  for 10 min. Centrifuge at 14,000rpm at  $4^{\circ}\text{C}$  for 6 min. Remove supernatant (cytosolic proteins) and label tube if needed for further analysis.
6. Add 500uL LB2+inhibitors to top of nuclei. Centrifuge at 12,000rpm  $4^{\circ}\text{C}$  for 6 min. Remove supernatant
7. Resuspend nuclei in 300uL LB3+inhibitors. Sonicate cells on high 30s on/off for 5 min. Centrifuge to pellet debris max. Centrifuge to pellet debris max speed at  $4^{\circ}\text{C}$  for 10 min. Do not forget the 10% input.

### Antibody-Beads Conjugates washes

1. Place beads on magnetic rack and remove supernatant. You can save supernatant and run in gel long with your samples to check to saturation of beads-Ab.
2. Wash beads with 300uL of LB3, vortex lightly to resuspend, spin briefly and remove supernatant. Repeat wash 4x for a total of 5 washes.
3. Quantitate amount of protein with Bradford prior to IP to ensure equal loading. Add the appropriate LB3 protein volume to antibody-bead conjugate. Rotate  $4^{\circ}\text{C}$  O/N

Western Blot analysis

2d

## 5 Day 3

1. Place beads in magnetic stand and remove protein mixture. Add 500uL LB3 to beads to wash, lightly vortex, spin briefly, place in stand and remove wash. Repeat this step at least 4x for a total of 5 washes.
2. For WB analysis, add about 50uL 5x Sample Buffer and 2.5uL 50mM DTT. Boil in  $99^{\circ}\text{C}$  for 10'. Place on magnetic rack and remove supernatant. This is your sample for Western Blot analysis.