

VERSION 2

MAY 04, 2023

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.eq2lypqzelx9/v2

Protocol Citation: Bjørn Peare Bartholdy, a.g.henry 2023. Biofilm growth with starch treatment.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.eq2lypqzelx9/v2> Version created by Bjørn Peare Bartholdy

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
 We use this protocol and it's working

Created: Mar 29, 2023

Last Modified: May 04, 2023

PROTOCOL integer ID:
 79636


Biofilm growth with starch treatment V.2

In 1 collection

Bjørn Peare Bartholdy¹, a.g.henry¹

¹Leiden University

BYOC



Bjørn Peare Bartholdy

Delft University of Technology, Leiden University

ABSTRACT

This is the main protocol to grow a calcifying oral biofilm model with starch treatments. It uses a 24 deepwell plate with the accompanying lid as substratum (polypropylene). The protocol takes 25 days to run, with daily tasks.

Modified from protocols by Sissons et al. (1991) and Extercate et al. 2010.

IMAGE ATTRIBUTION

Image created by Bjørn Peare Bartholdy using BioRender.com

GUIDELINES


The preferred substratum for inoculation is glass or hydroxyapatite. Plastic substrata can be used but are less effective, so surface treatment of the plastic is recommended (e.g. heated HCl or acetone).

Substrata should be autoclaved prior to use, both before and during the experiment.


MATERIALS

Solutions

Artificial saliva

Protocol	
	<div>NAME</div> <div>Artificial saliva</div>
<div>CREATED BY</div> <div>Bjørn Peare Bartholdy</div>	<div>PREVIEW</div>

CPMU

Protocol	
	<div>NAME</div> <div>CPMU</div>
<div>CREATED BY</div> <div>Bjørn Peare Bartholdy</div>	<div>PREVIEW</div>

20% (v/v) sterile glycerol in dH₂O

5% (w/v) sucrose in dH₂O

0.25% (w/v) potato starch in dH₂O

0.25% (w/v) wheat starch in dH₂O

0.50% (w/v) equal parts wheat (0.25%) + potato (0.25%) in dH₂O

Equipment

24 deepwell polypropylene plates (w. lid containing pegs suspended from the lid)

Shaking incubator

Powder free nitrile gloves

BEFORE START INSTRUCTIONS

Surfaces in the lab should be cleaned three times. Once with warm water and starch-free detergent, then with 5% NaOH, followed by a final cleaning with distilled water.

Saliva donor criteria

- Must have no/limited history of dental caries
- Must not have used antibiotics in the past 6 months
- Abstain from oral hygiene 24 hours prior to donation
- Refrain from eating and drinking (except water) 2 hours before donation

For experiments involving starches, donors should avoid eating starch-containing foods on the day of donation. To make this more bearable, saliva donation should take place in the morning before breakfast.

Prerequisites

Required solutions should be prepared beforehand.

- 20% (v/v) sterile glycerol in dH₂O
- Artificial saliva (required for day 0)
- Sucrose solution (required for day 0)
- Starch solutions (required for day 9)
- CPMU (required for day 15)

Protocol



NAME

Artificial saliva

CREATED BY

Bjørn Peare Bartholdy

PREVIEW

Protocol



NAME

CPMU

CREATED BY

Bjørn Peare Bartholdy

PREVIEW

Saliva collection

- 1 Saliva donors rinse their mouth with water for 30 seconds.
- 2 Stimulate saliva production by chewing tasteless gum or parafilm.
- 3 Collect the saliva by spitting into 50 ml plastic centrifuge tubes.

Note



Make sure donors wear gloves to avoid contamination with non-oral bacteria.

- 4 Make a 2-fold dilution of saliva in **sterile** [M] 20 % (v/v) glycerol and vortex the solution.

Day 0: Inoculation and feeding

- 5 Before inoculation, vortex the saliva solution again.
- 6 Pipette the saliva solution into the wells, so approx. 1-2 cm of the substratum is submerged.
- 7 Place the plate in the incubator at 36°C for 4 hours for static inoculation 🕒 04:00:00

4h

- 8 After inoculation, transfer the samples to a new plate containing the artificial saliva, and place in a shaking incubator at  30 rpm, 36°C for  04:00:00

4h

Protocol






NAME

Artificial saliva









CREATED BY

Bjørn Pearce Bartholdy

PREVIEW









- 8.1 Add  300 mL distilled (or deionized) dH₂O to a  1000 mL beaker, with stirring and heat  60 °C .





- 8.2 Add:


-  2.5 g  Mucin from porcine stomach (Type III) Sigma Aldrich Catalog #M1778
-  5 g  Trypticase™ Peptone Thermo Fisher Catalog #211921
-  10 g  Oxoid™ Proteose Peptone Thermo Fisher Catalog #LP0085B
-  5 g  Bacto Yeast Extract Becton-Dickinson



Let the reagents completely dissolve before continuing to the next step

- 8.3 Add:

-  2.5 g  KCl Contributed by users
-  0.35 g  NaCl Contributed by users
-  0.2 g  CaCl₂ Contributed by users
-  0.74 g  Sodium phosphate dibasic Sigma Aldrich Catalog #7558-79-4

-  0.54 g  NaHCO3 Contributed by users
-  2.5 mg  Hemin Contributed by users

8.4 Add the remaining  700 mL distilled (or deionized) dH₂O

8.5 Adjust to  7 with  NaOH Contributed by users and stirring

8.6 Transfer to two 1000 ml bottles, so half of each bottle is filled.

15m







Autoclave at  121 °C ,  1 Bar for  00:15:00 minutes

Safety information


Do NOT screw bottle caps on tightly.




Loosely screw the caps on the bottles or cover the tops with foil

8.7 Once the solution has cooled, add:

-  1 mg  Menadione Contributed by users
-  0.3 g  Urea Contributed by users
-  0.17 g  L-Arginine Contributed by users Catalog #A5006

8.8 Store in fridge at ca.  4 °C

Occasionally test the pH to ensure it stays around  7

9 Transfer the samples to a plate containing a 5% (m/v)  Sucrose Contributed by users solution for  00:06:00 , then transfer back to the artificial saliva and leave  Overnight .







6m



Plates with lids (substratum) in the incubator. Sucrose treatment plates covered with sterilised foil at the back.

Day 1-2: Feeding

8h 12m

- 10 First thing in the morning, transfer the samples to a new plate containing a 5% (m/v)  **Sucrose Contributed by users** solution for  00:06:00 . While in the sucrose solution, add more artificial saliva to the wells on the original plate that have been partially depleted overnight. 6m
- 11 After the 6 mins. return the samples to the plate with **artificial saliva**, and cover up the sucrose plate and leave for  08:00:00 . 8h
- 12 After 8 hours, transfer the samples back to the plate with 5% (m/v)  **Sucrose Contributed by users** solution for  00:06:00 . Transfer back to the artificial saliva and leave  Overnight . Dispose of the sucrose. 6m

Day 3: Inoculation and feeding

8h

- 13 Repeat steps from saliva collection and Day 0: Inoculation and feeding.

8h

Expected result

A layer of clear plaque should be visible on the substrata

- 14 Prepare a new plate with artificial saliva. Transfer the samples from the inoculation plate to the **artificial saliva**.

Day 4: Feeding

8h

- 15 Repeat steps 10 through 12

8h



Note

Prepare a new plate of artificial saliva every third day throughout the experiment. Every other morning, top up the wells with artificial saliva (so ca. 1-2 mm of the substratum is submerged).

Day 5: Inoculation

- 16 Repeat steps 1 through 9



Day 6-8: Feeding

- 17 Repeat steps 10 through 12




Day 9-14: Starch treatment

8h 6m

- 18 Transfer the samples to a plate with the starch treatment(s) for 00:06:00 at 60 rpm, 36°C

6m

Example setup:



Safety information



If you are using multiple starch treatments within a single plate, take care to avoid cross-contamination.

- 19 Transfer the samples back to the artificial saliva plate for 08:00:00 at 30 rpm, 36°C

8h



- 20 Transfer the samples to a plate with the starch treatment(s) for 00:06:00 at 60 rpm, 36°C

6m

- 21 Transfer the samples back to the artificial saliva plate at  30 rpm, 36°C and leave  Overnight

Day 15-24: Mineralisation

8h 30m

- 22 Transfer the samples to a plate containing the calcium phosphate monofluorophosphate urea (CPMU) solution for  00:06:00 at  60 rpm, 36°C

6m

Protocol






NAME

CPMU











CREATED BY

Bjørn Peare Bartholdy

PREVIEW

- 22.1 Add  300 mL distilled (or deionized) dH₂O to a  1000 mL beaker, with stirring and heat  60 °C

- 22.2 Add:

-  1.55 g  Calcium Chloride Sigma Aldrich
-  1.44 g  Sodium Phosphate monobasic Sigma Aldrich
-  0.72 g  Sodium Fluorophosphate Sigma Aldrich Catalog #344443
-  0.08 g  Magnesium Chloride Sigma Aldrich Catalog #AC223210010
-  30 g  Urea Sigma Aldrich

22.3 Add the remaining **700 mL** and keep stirring until precipitate has completely dissolved
Store in fridge at **4 °C**


23 Transfer the samples back to the **artificial saliva** for **02:00:00** at **30 rpm, 36°C** 2h

Put a lid on the plate with CPMU (or cover with foil) to prevent evaporation.

Repeat step 22 and 23, four more times every two hours.

24 Transfer the samples to a plate with the **starch treatment(s)** for **00:06:00** at **60 rpm, 36°C** 6m


25 Transfer the samples back to the **artificial saliva plate** at **30 rpm, 36°C** and leave **Overnight**



Analysis

26 Samples should be dried before sampling.

Transfer to a new plate with no liquid and leave in the incubator.



27 Once dried, samples are processed like archaeological samples.