

•





© DNA Barcoded Hydrogel Bead Synthesis V.1

Chris Cheshire¹

¹The Francis Crick Institute

In Development This protocol is published without a DOI.



Chris Cheshire The Francis Crick Institute

SUBMIT TO PLOS ONE

PROTOCOL CITATION

Chris Cheshire 2021. DNA Barcoded Hydrogel Bead Synthesis. protocols.io https://protocols.io/view/dna-barcoded-hydrogel-bead-synthesis-bjn2kmge

LICENSE

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 13, 2020

LAST MODIFIED

Feb 10, 2021

PROTOCOL INTEGER ID

40378

GUIDELINES

• Notes on <u>Acrylamide Polymerization</u>

MATERIALS TEXT

MATERIALS

Aldrich Catalog #W4502-1L

⊠Triton X-100 **Sigma**

Aldrich Catalog #T8787

Aldrich Catalog #P9541

Scientific Catalog #15568025

Aldrich Catalog #370533

Aldrich Catalog #P9416

Sodium chloride Sigma

Aldrich Catalog #S3014

⊠10 N NaOH **Sigma**

Aldrich Catalog #72068

Aldrich Catalog #A3678 **⊠** UltraPure™ 0.5 M EDTA pH 8.0 **Thermo Fisher** Scientific Catalog #15575020 ⊠ HFE-7500 3M Novec Engineered fluid Fluorochem Catalog #051243 Aldrich Catalog #227064 Span 80 Sigma Aldrich Catalog #S6760 Aldrich Catalog #A4058 Aldrich Catalog #A9926 ▼Tetramethylethylenediamine Sigma Aldrich Catalog #411019 Mineral oil Sigma Aldrich Catalog #M5310 **⊠** Brij™-35 30% Solution **Thermo Fisher** Scientific Catalog #20150 **⊠** 10x Isothermal Amplification Buffer Pack **New England** Biolabs Catalog #B0537S ■ Deoxynucleotide (dNTP) Solution Mix New England Biolabs Catalog #N0447L **⊠** Bst 2.0 DNA Polymerase **New England** Biolabs Catalog #M0537L ⊠ Exonuclease I (E. coli) New England Biolabs Catalog #M0293L ⊠ Exonuclease I Reaction Buffer New England Biolabs Catalog #B0293S Scientific Catalog #FD0804

□ FastDigest Buffer (10X) Thermo Fisher

Scientific Catalog #B64

Name	Recommended purchase
	amount
NaCl	500g-1Kg
KCI	500g-1Kg
Triton X-100	100 mL
Tween-20	100 mL
Perfluoro-1-octanol	25 g
HFE-7500	100 g
Tris-HCl	1 L
EDTA	4 x 100 mL
Hexane	1 L
Span 80	250 mL
NaOH	100 mL
Ammonium Persulphate	25 g
Acrylamide solution 40%	100 mL
Acrylamide/Bis-acrylamide 40% solution	100 mL
TEMED	100 mL
Mineral oil	100 mL
Brij-35	950 mL

SAFETY WARNINGS

COSHH sheet for this protocol can be found here:

https://docs.google.com/document/d/1uR8fdgWffDvPY2sY5ufCDKfpoQho4jAek9GcYnvBzvY/edit?usp=sharing

DISCLAIMER:

DISCLAIMER - FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

BEFORE STARTING

Glossary:

- Nuclease free water (NFW)
- Hydrogel Bead (HB)
- Tetramethylethylenediamine (TEMED)
- Microfluidic (MF)

NaCl

1 Prepare 5 M NaCl

Dissolve **■291.65** g of NaCl in **■700** mL approx. of NFW

Adjust final volume to 11 L

Store at & Room temperature for at least 1 year

KCl

2 Prepare 1 M KCI

Dissolve **□74.55** g of NaCl in **□900** mL approx. of NFW

Adjust final volume to 11 L

Store at & Room temperature for at least 1 year

Perfluorooctanol

3



Prepare 20% (vol/vol) Perfluorooctanol in HFE-7500 oil (20% (vol/vol)) - 20% PFO



Perform this step in a **fume-hood** with appropriate personal protective equipment and avoid contact with skin and eyes when handling this reagent.

Store the resulting mixture in a fume cupboard.

Combine **■8 mL** of **HFE-7500 oil** with **■2 mL** of **Perfluorooctanol**

Store at & Room temperature for at least 1 year

10% (vol/vol) Tween-20

4 Prepare 10% (vol/vol) Tween-20

Combine **■90 mL** of NFW with **■10 mL** of Tween-20

Store at & Room temperature for at least 6 months

10% (vol/vol) Triton X-100

5 Prepare 10% (vol/vol) Triton X-100

Combine ■90 mL of NFW with ■10 mL of Triton X-100

Store at A Room temperature for at least 6 months

DNA elution buffer

6 Prepare **■50 mL approx**. of **DNA elution buffer**

Reagent	Starting conc.	Reaction conc.	Volume to add
Tris-HCl	1 M	10 mM	500 μΙ
EDTA	0.5 M	100 μΜ	10 μΙ
NFW			49.5 mL
Total			50 mL

Make into □1 mL aliquots and store them at 8 -20 °C for up to 1 year.

Tris-EDTA-Tween (TET) buffer

7 Prepare **500 mL** of **TET buffer**

Reagent	Starting	Reaction	Volume
	conc.	conc.	to add
Tris-HCl	1 M	10 mM	5 mL
EDTA	0.5 M	10 mM	10 mL
Tween-20	10%	0.1 %	5 mL
NFW			480 mL
Total			500 mL

Filter solution through a → **0.2 µm** membrane.

Store at § Room temperature for up to 1 year

Tris-buffered saline-EDTA-Triton (TBSET) buffer

8 Prepare 11 L of TBSET buffer

Reagent	Starting	Reaction	Volume to
	conc.	conc.	add
Tris-HCl	1 M	10 mM	10 mL
EDTA	0.5 M	10 mM	20 mL
NaCl	5 M	137 mM	27.4 mL
KCI	1 M	2.7 mM	2.7 mL
Triton X-100	10 %	0.1 %	10 mL
NFW			822 mL
Total			1 L

Filter solution through a → 1-0.2 µm membrane.

Store at 8 Room temperature at least 6 months

9



Prepare 100 mL of 1% (vol/vol) Span-80 in hexane



Perform this step in a **fume-hood** with appropriate personal protective equipment and avoid contact with skin when handling this reagent.

Store the resulting mixture in a fume cupboard.

Combine **■99 mL** of hexane with **■1 mL** of Span-80

Store at & Room temperature in a fume cupboard for at least 6 months

Hydrogel bead wash buffer

10 Prepare 11 L of HB wash buffer

Reagent	Starting	Reaction	Volume
	conc.	conc.	to add
Tris-HCl	1 M	10 mM	10 mL
EDTA	0.5 M	100 μΜ	200 μΙ
Tween-20	10 %	0.05 %	5 mL
NFW			980 mL
Total			~ 1 L

Filter solution through a → 0.2 µm membrane.

Store at § Room temperature at least 6 months

STOP-25 buffer

11 Prepare 11 L of STOP-25 buffers

Reagent	Starting conc.	Reaction conc.	Volume to add
Tris-HCl	1 M	10 mM	10 mL
EDTA	0.5 M	25 mM	50 mL
Tween-20	10 %	0.1 %	10 mL
KCI	2 M	100 mM	50 mL
NFW			880 mL
Total			1 L

Filter solution through a → 0.2 µm membrane.

Store at & Room temperature at least 6 months

Prechill to § 4 °C before use.

STOP-10 buffer

12 Prepare 11 L of STOP-10 buffer

Reagent	Starting	Reaction	Volume
	conc.	conc.	to add
Tris-HCl	1 M	10 mM	10 mL
EDTA	0.5 M	10 mM	20 mL
Tween-20	10 %	0.1 %	10 mL
KCI	2 M	100 mM	50 mL
NFW			910 mL
Total			1 L

Filter solution through a → • 0.2 µm membrane.

Store at 8 Room temperature at least 6 months

Hybridization buffer

13 Prepare **500 mL** of **Hybridization buffer**

Reagent	Starting	Reaction	Volume
	conc.	conc.	to add
Tris-HCl	1 M	10 mM	5 mL
EDTA	0.5 M	100 μΜ	100 μΙ
Tween-20	10 %	0.1 %	5 mL
KCI	2 M	330 mM	82.5 mL
NFW			407 mL
Total			500 mL

Filter solution through a → **0.2 µm** membrane.

Store at § Room temperature at least 6 months

QC Buffer

14

Α	В	С	D
Reagent	Starting	Reaction	Volume to
	conc.	conc.	add
Tris-HCl	1 M	?	250 μΙ
EDTA	0.5 M	?	500 μΙ
Tween-20	10% (vol/vol)	?	250 μΙ
KCI	1 M	1 M	49 mL
NFW			0 mL
Total			50 mL

4× Acrylamide/bis-acrylamide (AB) solution

15





Acrylamide and bis-acrylamide are potent neurotoxins.

Perform this step in a **fume-hood** with appropriate personal protective equipment and avoid contact with skin and eyes when handling this reagent.

Reagent	Starting	Reaction	Volume
	conc.	conc.	to add
Acrylamide solution	40 %	-	2.58 mL
Acrylamide/bis-acrylamide solution	40 %	-	3.6 mL
NFW			3.82 mL
Total			10 mL

Filter solution through a → 0.2 µm membrane.

Store at 8 4 °C at least 6 months

10% (wt/vol) ammonium persulfate (APS)

16



Prepare _5 mL of 10% (wt/vol) ammonium persulfate (APS)



Perform this step in a **fume-hood** with appropriate personal protective equipment and avoid contact with skin when handling this reagent.

Combine **□500 mg** of **ammonium persulfate** with **□4.5 mL** of **NFW**

Adjust final volume to 35 mL

Filter solution through a → 0.2 µm membrane

Aliquots of □100 µl can be stored at 8-20 °C for up to 6 months

Acrydite-modified primer

17 5'-/acrydite/photocleavable spacer/CGATGACGTAATACGACTCACTATAG GGATACCACCATGGCTCTTTCCCTACA CGACGCTCTTC-3'

Resuspend and dilute to [M]250 Micromolar (µM) in DNA Elution Buffer

Acrydite-modified DNA primers are light-sensitive, as they contain a photocleavable (1-(2-nitrophenyl)-ethyl-based) moiety.

At all times protect the DNA primers from exposure to UV light (<400 nm). Noticeable primer deactivation may occur at prolonged ($\ge 20 \text{ min}$) exposure to ambient white light illumination.

Perform in dim illumination and away from UV light sources.

Carrier Oil (1% surfactant)

18 Prepare **2.5 mL** of carrier oil mix

Combine ■1.25 mL of HFE-7500 oil with ■1.25 mL of dSurf (2%)

Vortex well

Store at & Room temperature and prepare freshly for the experiment run

Carrier oil / TEMED mix

19



Prepare **■2.5 mL** of **TEMED mix**



Perform this step in a **fume-hood** with appropriate personal protective equipment and avoid contact with skin when handling this reagent.

Combine **□2.5 mL** of **carrier oil** with **□10 µl** of **TEMED**

Vortex well

Store at § Room temperature and prepare freshly for the experiment run

Neutralization buffer

20 Prepare **500 mL** of **Neutralization buffer**

Reagent	Starting	Reaction	Volume
	conc.	conc.	to add
Tris-HCl	1 M	100 mM	50 mL
EDTA	0.5 M	10 mM	10 mL
Tween-20	10 %	0.1 %	5 mL
NaCl	1 M	100 mM	50 mL
NFW			385 mL
Total			500 mL

Filter solution through a → • 0.2 µm membrane.

Store at & Room temperature at least 6 months

Prepare MF Equipment

- 21 Cut tubing for the following lengths using the template guides found in the template bag:
 - 2 x lengths of each size of 1/32
 - 3 x lengths of the 1/16 tube

Collect MF equipment:

- 2 x p-valve head block
- 2 x long 1/32 tubing
- 2 x p-valve attachment kit

Clean heads and o-rings and bring to fumehood

Acrylamide-primer mix

23



Prepare $\blacksquare 1$ mL of Acrylamide-primer mix in a $\blacksquare 2$ mL tube



Acrylamide and bis-acrylamide are potent neurotoxins.

Perform this step in a **fume-hood** with appropriate personal protective equipment and avoid contact with skin and eyes when handling reagents in this list.

Acrydite-modified DNA primers are light-sensitive, as they contain a photocleavable (1-(2-nitrophenyl)-ethyl-based) moiety.

At all times protect the DNA primers from exposure to UV light (<400 nm). Noticeable primer deactivation may occur at prolonged ($\ge 20 \text{ min}$) exposure to ambient white light illumination.

Perform in dim illumination and away from UV light sources.

Reagent	Starting	Reaction	Volume
	conc.	conc.	to add
TBSET buffer	?	?	100 μΙ
AB Solution	?	?	250 μΙ
Acrydite-modified DNA primer	250 μΜ	50 μΜ	200 μΙ
10% (wt/vol) APS	10 %	0.3 %	30 μΙ
NFW			420 µl
Total			1000
			μΙ

24



Vortex solution to mix and centrifuge - 31000 x g, Room temperature , 00:00:05

25 Place 500μ of mineral oil on top of the solution

Prepare MF run

- 26 Place Carrier oil/TEMED mix and Acrylamide-primer mix into p-valve system.
 - Connect the sample tube to p-valve head
 - Cover tube with foil

- Slide green spacer over 1/32 tube to center
- Thread green space and tube through the cream seal and place the blue plug in reverse over the green spacer facing the cream seal thread
- Place 1/32 tubing through the aluminium seal in top until it hits the bottom of tube
- Slide the green spacer and seals down until they hit the aluminium thread
- Screw seal thread until the 1/32 tubing does not move anymore
- Place both into tube stand
- $27 \qquad \text{Move tube stand from fume-hood to the MF bench} \\$
- 28 Attach flow units and connect conversion adapter to 1/32 tubing coming from each p-valve
 - Slide flow-unit seal over the 1/32 tube coming from p-valve and connect to LHS of flow unit. Tighten until the tube does not move anymore
 - Attach the short 1/32 tubing to the RHS of the flow unit using the same method
 - Attach red conversion adapter to the 1/16 tubing
 - Attach 1/32 from RHS flow-unit to the other side of the conversion adaptor
- 29 Place MF bead device onto microscope mount upside down and secure
- 30 Connecting tubing to MF device including output
- 31 Attach blue pressure tubing to p-valves
- 32 Place p-valve head unit over collection tube and place into tube rack.

Slide output tubing into top of p-valve but do not seal.

Wrap output tube in foil.

33 Turn on power

Turn on microscope lighting

Connect the camera to a laptop and get a live picture

34 Adjust lighting, slide position and focus until droplet junction is in ROI and focused/lit properly

MF run start-up

35 TODO

MF run shutdown

36 TODO

Polymerize beads

37



Transfer HB tubes to an incubator overnight at § 65 °C protected from light.

MF Clean

38 TODO

HB Clean-up

39 Remove the top (mineral oil) and bottom (carrier oil) phases from the 2-ml collection tubes, leaving **300 μl approx**. of polymerized HBs per tube.

TODO: PICTURE

40 Add **3500 μl** of **TBSET buffer** on top of the creamy-looking emulsion.

TODO: PICTURE

41



Release HBs from the emulsion by adding 1 mL of 20% PFO

Vortex well then centrifuge at \$\&\pmo 5000 x g, Room temperature , 00:00:30

HBs should appear as a milky (or semitransparent) phase.

TODO: PICTURE

42 Remove the bottom 20% PFO phase.

Repeat **o** go to step #41 until the milky phase is transformed into a solid, well-packed, translucent mass on top of the 20% PFO phase.

Remove the final bottom 20% PFO phase.

43



Add 1 mL of 1% (vol/vol) Span-80 in hexane

Vortex well then centrifuge at \$\&\pmo 5000 x g, Room temperature , 00:00:30

This time, the HBs will appear as a solid mass at the bottom of the tube.

Remove the hexane layer.

Repeat this step.

44 Add 11 mL of TBSET buffer and resuspend HBs by vortexing.

45 Transfer HBs to **□15 mL** tubes.





Fill tubes with TBSET buffer.

Vortex then centrifuge at 3000 x g, Room temperature, 00:03:00

The HBs will sediment whereas traces of hexane will appear as a milky layer on the top. Carefully aspirate the supernatant, keeping the tip of the serological pipette at the liquid-air interface.

If the HBs do not sediment, extend the centrifuge time.

47 Repeat **ogo to step #46** until the top milky layer disappears.

48 TODO: QC the bead size

Put 20 microlitre on slide and check in microscope, take picture and use imageJ to check measurement - using python script for size of particles etc

get a histogram



49 🐧

Pass the HBs through a → 70 µm cell strainer on top of a □50 mL tube.

Use ice-cold **TBSET buffer** to facilitate the passage of the beads through the mesh.

50



Centrifuge at 3000 x g, Room temperature, 00:10:00

Remove supernatant.

Transfer batches of HBs into 15 mL tubes and collect remaining into 50 mL tube. Wash the sides of the tube with 11 mL of TBSET buffer to release HBs stuck to the side for the 50 mL tube.

51



Store at § 4 °C for up to 6 months protected from light.

Bead QC

52 Prepare FAM-PE1 probe 5'-/6-FAM/AGATCGGAAGAGCGTCGTGTAGGGAAAGAG-3' (standard desalting)

Prepare a [M]100 Micromolar (µM) solution in DNA elution buffer.

This solution is stable for at least 1 year at 8 -20 °C

53 Spin down hbs

add □10 µl to □1.5 mL tubes and 1.4 mL of QC buffer

Vortex

Cent 1 min @ 1,000g

remove supermatent

adjust final QC buffer vol to 36 µl

Add **□4 µI** of probe

mix by pipetting up and down

incubate room temp for 20 mins

Wash x 3 - add 1.4ml of QC buffer, vortex, centrifuge 1 min @ 1000g - remove supernatant

should have ~40 left in each tube

Split into 2 tubes

Place on slide and image

expose other tube to UV

place on slide and image

Denaturation solution

Prepare **■250 mL** of denaturation solution.



NaOH is highly corrosive; wear appropriate personal protective equipment when handling this reagent.

Perform this step in a **fume-hood** with appropriate personal protective equipment and avoid contact with skin and eyes when handling reagents in this section.

Reagent	Starting	Reaction	Volume
	conc.	conc.	to add
NaOH	10 M	150 M	3.75 mL
Brij-35	30 %	0.5 %	4.2 mL
NFW			242 mL
Total			250 mL

Filter solution through a → • 0.2 µm membrane.

Store at & Room temperature

PREPARE FRESH