



VERSION 2

APR 05, 2024

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.3byl4qo5zv05/v2

Protocol Citation: daniel.dautan
daniel, Per Svenningsson 2024.
Rotarod (40 RPM and 20 RPM
with docking) . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.3byl4qo5zv05/v2>Version
created by [Eileen Ruth Torres](#)

License: This is an open
access protocol distributed under
the terms of the [Creative
Commons Attribution License](#),
which permits unrestricted use,
distribution, and reproduction in
any medium, provided the original
author and source are credited

Protocol status: Working
We use this protocol and it's
working

🌐 Rotarod (40 RPM and 20 RPM with docking) V.2

daniel.dautan daniel¹, Per Svenningsson¹

¹Karolinska Institute Stockholm

ASAP Collaborative Research Network

Kapliitt Protocols



Eileen Ruth Torres

Weill Cornell Medicine

ABSTRACT

This protocol describes motor function testing in mice using the rotarod by the Svenningsson lab.

Two options for the protocol are included: 1) 40 RPM constant speed after habituation and 2) 20 RPM with docking.

MATERIALS

Mouse Rotarod (Ugo Basile 47650)
50% ethanol

Created: Feb 22, 2024

Last Modified: Apr 05, 2024

PROTOCOL integer ID: 95631

Keywords: ASAPCRN

Funders Acknowledgement:

Aligning Science Across
Parkinson's
Grant ID: 020608

- 1 Transfer the cage into the experimental room at least 1 hr before testing.

- 2 To avoid stress, always put the cage into a ventilated cage rack.

- 3 Using a paper tube, transfer each mouse of the same cage onto the rotarod (Ugo Basile, 47650). Always maintain order, mice one into line one, mice 2 into line 2 etc...

- 4 Wait 2-3 minutes for the mice to habituate to the apparatus.

- 5 Lift the collecting platform and set up the protocol: either max speed of constant 40 RPM or 20 RPM docking.

STEP CASE

40 RPM 4 steps

- 6** Following habituation, start the rotarod program to reach a speed of 40 RPM within the first 5s of the protocol.
- 7** The time to fall is defined as the moment the mice fall in the receptacle or the moment the mice turn around the rotor more than 5 consecutive times. The maximum trial length is 2 min.
- 8** Repeat the experiment 3 times with at least 6 hr between each trial. If an animal does not fall or turn more than 5 times within 120s of the trial, stop the trial and record “120” as the individual time.
- 9** Analyze for each mouse, the time to fall for each session, the average and and max time.