



Nov 30, 2021

MPRA library preparation

Gerald Raffl¹, Boyan Bonev¹

¹Helmholtz Pioneer Campus





protocol



Protocol to generate both DNA and RNA MPRA libraries.

Gerald Raffl, Boyan Bonev 2021. MPRA library preparation. **protocols.io** https://protocols.io/view/mpra-library-preparation-bxdtpi6n

_____ protocol,

Aug 16, 2021

Nov 30, 2021

52371

Devices:

Centrifuge

2 ThermoMixers (Eppendorf)

Thermal Cycler

Magnet for magnetic bead purifications

Qubit fluorometer (ThermoFisher Scientific)

Real-Time Thermal Cycler for qPCR

Agilent 2100 Bioananlyzer (Agilent)

Reagents:

Nuclease free water (Ambion Invitrogen, Cat. N: AM9937)

2X digestion buffer (Zymo Research, Cat. N: D3050-1)

Proteinase K (lyophilized) & Storage Buffer (Zymo Research, Cat. N: D3001-2)

Quick-DNA/RNA Microprep Plus Kit (Zymo Research, Cat. N: D7005)

TURBO DNA-free™ Kit (ThermoFisher Scientific, Cat. N: AM1907)

Qubit™ Assay Tubes (ThermoFisher Scientific, Cat. N: Q32856)

Qubit™ dsDNA HS Assay Kit (ThermoFisher Scientific, Cat. N: Q32854)

Qubit™ RNA HS Assay Kit (ThermoFisher Scientific, Cat. N: Q32852)

RNA 6000 Pico Kit (Agilent, Cat. N: 5067-1513)

Maxima H Minus Reverse Transcriptase (200 U/μL) (ThermoFisher Scientific, Cat. N: EP0753)

Oligo(dT)18 Primer (ThermoFisher Scientific, Cat. N: SO132)

dNTP Solution Mix (NEB, Cat. N: N0447)

RNasin® Plus Ribonuclease Inhibitor (Promega, Cat. N: N2611)

NEBNext® Ultra™ II Q5® Master Mix (NEB, Cat. N: M0544)

Eppendorf PCR tubes, 0.2 mL (Eppendorf, Cat. N: 30124359)

DNA LoBind® Tubes, 1.5 mL (Eppendorf, Cat. N: 30108051)

AMPure XP (Beckman Coulter, Cat. N: A63881)



1

Primers:

name	sequence
RV_univ_MPRA	CGACGCTCTTCCGATCT
FWD_mScar_Tn7_10UMI_3	GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGNNNNNNNNNN
FWD_CRS_Tn7	GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGCCAGGACCGGATCAACT
P5NEXTPT5	AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCG*A*T*C*T
Ad2.1_TAAGGCGA	CAAGCAGAAGACGCCATACGAGATTCGCCTTAGTCTCGTGGGCTCGGAGATGT
Ad2.2_CGTACTAG	CAAGCAGAAGACGGCATACGAGATCTAGTACGGTCTCGTGGGCTCGGAGATGT
Ad2.3_AGGCAGAA	CAAGCAGAAGACGCCATACGAGATTTCTGCCTGTCTCGTGGGCTCGGAGATGT
Ad2.4_TCCTGAGC	CAAGCAGAAGACGGCATACGAGATGCTCAGGAGTCTCGTGGGCTCGGAGATGT
Ad2.6_TAGGCATG	CAAGCAGAAGACGGCATACGAGATCATGCCTAGTCTCGTGGGCTCGGAGATGT
Ad2.7_CTCTCTAC	CAAGCAGAAGACGCCATACGAGATGTAGAGAGGTCTCGTGGGCTCGGAGATGT
Ad2.8_CAGAGAGG	CAAGCAGAAGACGGCATACGAGATCCTCTCTGGTCTCGTGGGCTCGGAGATGT
Ad2.9_GCTACGCT	CAAGCAGAAGACGCATACGAGATAGCGTAGCGTCTCGTGGGCTCGGAGATGT
Ad2.10_CGAGGCTG	CAAGCAGAAGACGGCATACGAGATCAGCCTCGGTCTCGTGGGCTCGGAGATGT
Ad2.11_AAGAGGCA	CAAGCAGAAGACGGCATACGAGATTGCCTCTTGTCTCGTGGGCTCGGAGATGT
Ad2.12_GTAGAGGA	CAAGCAGAAGACGGCATACGAGATTCCTCTACGTCTCGTGGGCTCGGAGATGT
Ad2.13_GTCGTGAT	CAAGCAGAAGACGGCATACGAGATATCACGACGTCTCGTGGGCTCGGAGATGT
Ad2.14_ACCACTGT	CAAGCAGAAGACGGCATACGAGATACAGTGGTGTCTCGTGGGCTCGGAGATGT
Ad2.15_TGGATCTG	CAAGCAGAAGACGGCATACGAGATCAGATCCAGTCTCGTGGGCTCGGAGATGT
Ad2.16_CCGTTTGT	CAAGCAGAAGACGGCATACGAGATACAAACGGGTCTCGTGGGCTCGGAGATGT
Ad2.17_TGCTGGGT	CAAGCAGAAGACGGCATACGAGATACCCAGCAGTCTCGTGGGCTCGGAGATGT
Ad2.18_GAGGGGTT	CAAGCAGAAGACGCCATACGAGATAACCCCTCGTCTCGT

5m

5m

DNA & RNA extraction from fixed nuclei

- 1 Preheat a ThermoMixer to 55°C.
- Pellet sorted nuclei in a 1.5 mL tube at 2500 x g, 4°C for 3min.

 Resuspend pellet in 1X ProtK digestion buffer (95 μL nuclease free water, 95 μL 2X digestion buffer, 10 μL Proteinase K)

3 Incubate at 55°C, 600 rpm for 1h. (Proteinase K digestion)
In the meantime, preheat another ThermoMixer to 65°C.

- 4 Incubate at 65°C, 600 rpm for 15min. (Reverse Crosslinking)
- 5 Centrifuge at maximum speed for 3min. (Get rid of cell debris) In the meantime, set the 2 ThermoMixers to 37°C and 26°C, respectively.

protocols.io

2

- Transfer supernatant (190 μL) to a new tube, add 190 μL RNA/DNA lysis buffer (from Quick-DNA/RNA Microprep Plus Kit) and mix thoroughly.
- 7 Continue with the Quick-DNA/RNA Microprep Plus Kit following the manufacturer's protocol. Perform in-column DNase digest for the RNA samples.

30m

- 8 For **DNA** samples:
 - Elute in 21 μL nuclease free water.
 - Measure DNA concentration with the Qubit[™] dsDNA HS Assay Kit, using 1 µL DNA.
 - Freeze at -20°C for later use.
- 9 For RNA samples:
 - Elute in 15 µL nuclease free water.
 - Continue with TURBO DNase digestion and Reverse Transcription.

TURBO DNase digest and Reverse Transcription

1h 40m

- 10 Add to each RNA sample:
 - 1.5 µL TURBO DNase buffer (10X)
 - 1 μL TURBO DNase

Incubate at 37°C, 600rpm for 20min

11 Add 2 μ L inactivation reagent to each sample. Incubate at 26°C, 1600rpm for 5min.

5m

20m

12 /

Transfer each sample to $0.5\,\text{mL}$ tubes (critical, allows better aspiration of supernatant). Centrifuge at $10'000\,\text{x}$ g for 90s.

Aspirate the supernatant (15 μ L). Do not take up any inactivation reagent! Measure RNA concentration with the Qubit[™] RNA HS Assay Kit, using 1 μ L RNA.

14 🛠

Optional: Determine the RIN value by analyzing the RNA samples with the RNA 6000 Pico Kit, following the manufacturer's protocol.

- 15 Use 12 µL RNA for Reverse Transcription. For that, set up the following reaction in 0.2 mL PCR tube:
 - 12 μL RNA
 - 1.5 µL dNTP mix
 - 1.5 µL Oligo(dT)18 Primer
- 16 Incubate in a thermocycler with heated lid:

5m

protocols.io

- 5min at 65°C
- 17 Place reaction on ice for 1min.

Add 15 μ L 2X RT mix to each sample (6 μ L 5X RT buffer, 7.5 μ L nuclease free water, 1 μ L Maxima H- RT, 0.5 μ L RNasin plus).

18 Incubate in a thermocycler with heated lid:

40m

- 30min at 50°C
- 5min at 85°C

19 (1)

30m

Transfer each sample to a DNA low binding 1.5 mL tube.

Clean up the cDNA with 1.5X vol. (= $45 \,\mu$ L) of AMPure XP magnetic beads following the manufacturer's protocol. Elute the cleaned cDNA in 21 μ L nuclease free water.

Freeze cDNA at -20°C for later use.

Attaching UMIs

In the two subsequent sections, DNA sequencing libraries are generated from both DNA and cDNA (RNA) samples. This is done in the exact same way for both sample types.

In this first PCR, UMIs are attached by 3 PCR cycles. Primer binding sites are in the C-terminal part of mScarlet-I (FWD_mScar_Tn7_10UMI_3) and immediately downstream of the 12bp barcode (RV_univ_MPRA). The fragment created in this PCR has a length of 196bp.

Set up the following PCR reaction:

- 17 µL DNA/cDNA from step 🕁 go to step #8 and step 🕁 go to step #19 , respectively
- 3 uL nuclease free water
- 2.5 µL primer FWD_mScar_Tn7_10UMI_3 (10µM stock)
- 2.5 μL primer RV_univ_MPRA (10μM stock)
- 25 µL NEBNext Ultra II Q5 Master Mix

Don't use all 20 µL of cDNA: Leave 3 µL of cDNA for qPCR on marker genes.

- 21 Incubate in a PCR cycler with heated lid:
 - 30s 98°C
 - 3 cycles (10s 98°C; 30s 65°C; 1min 72°C)
 - 5min 72°C
 - hold 4°C
- Transfer each sample to a DNA low binding 1.5 mL tube.

Clean up each sample with 1.2X vol. (= $60 \, \mu L$) of AMPure XP magnetic beads following the manufacturer's protocol.

Elute in 20 µL nuclease free water.

Library amplification

73 The library amplification step is split up in 2 separate PCRs. In the first PCR, the P7 illumina adapter is added.



Then, in the remaining cycles, also the P5 illumina adapter is added. By splitting up the library amplification in two PCRs (with a total of 20-30 cycles, depending on input), the likelihood of primer dimer formation is drastically reduced in low-input samples.

Set up the following PCR reaction:

- 20 µL sample from step **⑤ go to step #22**
- 2.5 μL primer Ad2.X_index (2μM stock)
- 2.5 μL primer RV_univ_MPRA (2μM stock)
- 25 µL NEBNext Ultra II Q5 Master Mix
- 24 Incubate in a PCR cycler with heated lid:
 - 30s 98°C
 - X cycles (10s 98°C; 90s 65°C)
 - 5min 65°C
 - hold 4°C

X cycles: 10 for DNA samples, 12 for cDNA samples

25 Transfer each sample to a DNA low binding 1.5 mL tube.

Clean up each sample with 1.2X vol. (= $60 \, \mu L$) of AMPure XP magnetic beads following the manufacturer's protocol.

Elute in 10 µL nuclease free water.

26 Use 1.5 μL sample from step **go to step #25** for qPCR to estimate the remaining PCR cycles.

Set up the following qPCR reaction:

- 1.5 µL sample
- 0.5 μL primer Ad2.X_index (2μM stock)
- 0.5 µL primer P5NEXTPT5 (2µM stock)
- 2.5 μL nuclease free water
- 5 µL 2X Luna Universal qPCR Master Mix
- 27 Incubate in a Real-Time thermocycler:
 - 5mins 95°C
 - 45 cycles (10s 98°C; 10s 60°C; 20s 72°C)
 - Cooling
- 28 /

Determine the cycle number, where the amplification curve was at its half maximum. Subtract 1 from this number. This is the required remaining cycle number.

- 29 Set up the following PCR reaction:
 - 8 µL sample from step **⑤ go to step #25**
 - 12 µL nuclease free water
 - 2.5 μL primer Ad2.X_index (2μM stock)
 - 2.5 μL primer P5NEXTPT5 (2μM stock)
 - 25 µL NEBNext Ultra II Q5 Master Mix
- 30 Incubate in a PCR cycler with heated lid:
 - 30s 98°C
 - X cycles (10s 98°C; 90s 65°C)



- 5min 65°C
- hold 4°C

X cycles: number determined in step ogo to step #28

Transfer each sample to a DNA low binding 1.5 mL tube. Clean up each sample with 0.8X vol. (= $40 \,\mu$ L) of AMPure XP magnetic beads following the manufacturer's protocol.

The lower AMPure XP beads ratio here gets rid of potential primer dimers.

Elute in 10 µL nuclease free water. This is the final library.

32 Quantify by measuring 1 µL with the Qubit™ dsDNA HS Assay Kit.

Analyze the final library with the High Sensitivity DNA Kit on the Agilent 2100 Bioanalyzer.



The final library should have a size of 269bp.

CRE-barcode association

In this section, the DNA sequencing library for the CRE-barcode association is created. For that, the **MPRA pre- pool** (without min.Promoter_mScarlet-I) is used as input.

Set up the following PCR reaction:

- 5 ng MPRA pre-pool
- to 20 µL: nuclease free water
- 2.5 μL primer FWD_CRS_Tn7 (10μM stock)
- 2.5 μL primer RV_univ_MPRA (10μM stock)
- 25 µL NEBNext Ultra II Q5 Master Mix
- 34 Incubate in a PCR cycler with heated lid:
 - 30s-98°C
 - 3 cycles (10s 98°C; 30s 65°C; 3min 72°C)
 - 5min 72°C
 - hold 4°C

 $The \ longer \ extension \ time \ prevents \ chimeric \ DNA \ annealing \ caused \ by \ incomplete \ elongation.$

35 Transfer the sample to a DNA low binding 1.5 mL tube.

Clean up the sample with 1.2X vol. (= $60 \, \mu L$) of AMPure XP magnetic beads following the manufacturer's protocol.

Elute in 20 µL nuclease free water.

- 36 Set up the following PCR reaction:
 - 20 μL PCR product
 - 2.5 μL primer Ad2.X_index (2μM stock)

- 2.5 μL primer P5NEXTPT5 (2μM stock)
- 25 µL NEBNext Ultra II Q5 Master Mix
- 37 Transfer the sample to a DNA low binding 1.5 mL tube. Clean up the sample by double size selection with 0.5X vol. and 0.8X vol. of AMPure XP magnetic beads following the manufacturer's protocol. Elute in 10 μ L nuclease free water.
- 38 Quantify by measuring 1 μL with the Qubit™ dsDNA HS Assay Kit.

 Analyze the final library with the High Sensitivity DNA Kit on the Agilent 2100 Bioanalyzer.



The final library should have a size of 455bp.