



Dec 08, 2020

Protein Transfer

Steven J Burgess¹¹University of Illinois at Urbana-Champaign*In Development*

This protocol is published without a DOI.

Ag SynBio Lab UIUC



Steven Burgess

University of Illinois at Urbana-Champaign

ABSTRACT

This protocol is based on the assumption that TGX pre-cast gels have been used, the procedure is very similar for self-made gels, just that the running time is longer. Using the transfer packs is economically similar to wet-transfer using MeOH.

Note: Clean the apparatus immediately after use with dH₂O and dry. This prevents the build-up of salts and rust which impair performance.

Literature:

Full manual: <https://www.bio-rad.com/webroot/web/pdf/lsr/literature/10020688.pdf>Quick-start: <https://www.bio-rad.com/webroot/web/pdf/lsr/literature/10016505D.pdf>Transfer pack: <https://www.bio-rad.com/webroot/web/pdf/lsr/literature/10019593D.pdf>

PROTOCOL CITATION

Steven J Burgess 2020. Protein Transfer. **protocols.io**
<https://protocols.io/view/protein-transfer-bqhvm66>



LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 07, 2020

LAST MODIFIED

Dec 08, 2020

PROTOCOL INTEGER ID

45333

PARENT PROTOCOLS

In steps of

[SDS-PAGE gel electrophoresis](#)

MATERIALS TEXT

- TransBlot^R-TurboTM Transfer System (Bio-Rad Laboratories; [1704150](#))
- Trans-Blot Turbo Mini 0.2 µm Nitrocellulose Transfer Packs (Bio-Rad Laboratories; [1704158](#))

ABSTRACT

This protocol is based on the assumption that TGX pre-cast gels have been used, the procedure is very similar for self-made gels, just that the running time is longer. Using the transfer packs is economically similar to wet-transfer using MeOH.

Note: Clean the apparatus immediately after use with dH₂O and dry. This prevents the build-up of salts and rust which impair performance.

Literature:

Full manual: <https://www.bio-rad.com/webroot/web/pdf/lsr/literature/10020688.pdf>

Quick-start: <https://www.bio-rad.com/webroot/web/pdf/lsr/literature/10016505D.pdf>

Transfer pack: <https://www.bio-rad.com/webroot/web/pdf/lsr/literature/10019593D.pdf>

- 1 Place the membrane and bottom stack in the middle of the cassette base
- 2 Place pre-cast gel in the middle of the membrane and roll to remove air-bubbles
- 3 Place the top stack on top of the gel, gently roll
- 4 Close the cassette lid, taking care not to disturb the gel
- 5 Place the cassette in the Trans-Blot Turbo and follow the instructions on the machine (Fast protocol TGX gel).

Trans-Blot® Turbo™
Protein transfer apparatus

Bio-rad Laboratories 1704150EDU [↗](#)

- 6 Either dry membrane and store at **4 °C** for later use or proceed immediately to fluorescent western protocol