



Jun 30, 2021

T cell isolation from mouse tissues

Kristin Anderson¹¹University of Washington and Fred Hutchinson Cancer Research Center

1 Works for me



Share

dx.doi.org/10.17504/protocols.io.sqnedve

Kristin Anderson

ABSTRACT

This protocol outlines the steps for isolating murine T cells from mouse blood, ascites, spleen, and ID8 tumors.

In brief: Spleens were disrupted by mechanical force through a 70mm filter (BD Falcon, ref: 352350) and red blood cells lysed using ACK lysis buffer (Gibco, cat: A1049201). Ovarian tumors were excised and normal adjacent tissue removed. Tumors were minced using a GentleMACS dissociator (Miltenyi), incubated for 10 minutes on a rotator (Miltenyi) at 37°C (medium spin setting), disrupted through a 70mm filter, washed with RPMI containing 5% FBS, and lymphocytes were purified on a 44/67% Percoll gradient (800 $\times g$ at 23°C for 20 minutes; GE Healthcare, cat: 17-0891-01).

DOI

dx.doi.org/10.17504/protocols.io.sqnedve

PROTOCOL CITATION

Kristin Anderson 2021. T cell isolation from mouse tissues. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.sqnedve>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 17, 2018

LAST MODIFIED

Jun 30, 2021

PROTOCOL INTEGER ID

14830

MATERIALS TEXT

MATERIALS

Percoll Sigma

Aldrich Catalog #17-0891-01

ACK Lysing Buffer Gibco - Thermo

Fisher Catalog #A10492-01

Falcon cell strainer, 70um,

[nylon Corning Catalog #352350](#)

Spleen digestion

1. Harvest cells by smashing spleen through a 70um filter.
2. Spin down cells at 369 $\times g$ for 4-5 minutes at 4C.

3. Remove supernatant and lyse red blood cells with 2 mL ACK lysis buffer for 2 minutes at room temperature.
4. Quench lysis using 10mL T cell media.
5. Spin down cells at 369*xg* for 4-5 minutes at 4C.
6. Remove supernatant and resuspend cells in T cell media. Store on ice until use.

Blood and Ascites digestion

- 2
1. Spin down cells at 369*xg* for 4-5 minutes at 4C.
2. Remove supernatant and lyse red blood cells with 2 mL ACK lysis buffer for 2 minutes at room temperature.
3. Quench lysis using 2mL T cell media. Spin down cells at 369*xg* for 4-5 minutes at 4C.
4. Remove supernatant and lyse red blood cells a second time with 2 mL ACK lysis buffer for 2 minutes at room temperature.
5. Quench lysis using 2mL T cell media. Spin down cells at 369*xg* for 4-5 minutes at 4C.
6. Remove supernatant and resuspend cells in T cell media. Store on ice until use.

ID8 ovary tumor digestion

- 3
1. Excise tumor from animal and remove normal tissue. Be sure to remove any lymph nodes from tumor! Store in RPMI containing 5% FBS on ice until processing.
2. Transfer tumor tissue into 3mL RPMI+5%FBS in a GentleMACs C tube.
3. Mince tumors using a GentleMACS dissociator (Miltenyi). Use protocol 'm_imptumor_01'.
4. Transfer tissue to a 50 mL conical tube. Add 35 mL of collagenase (pre-warmed to 37C).
5. Put tubes on a MACS Mix Rotator and place in the 37C incubator rotating on medium speed for 10 minutes.
6. Spin down tubes at 369*xg* for 4-5 minutes at 4C.
7. Remove supernatant and resuspend tissue in 5mL RPMI+5%FBS. Smash sample through a 70um filter.
8. Spin down tubes at 369*xg* for 4-5 minutes at 4C.
9. Resuspend pellet in 5mL of 44% Percoll.
10. Underlay with 3-5mL of 67% Percoll.
11. Spin samples for 20 minutes at 809*xg* at room temperature with NO BRAKE.
12. Remove and discard fat layer floating on top. Gently remove the cells at the interface and transfer to a fresh tube.
13. Wash cells with 15mL RPMI+5% FBS. Spin cells at 518*xg* for 10 minutes at 4C.
14. Remove supernatant. Resuspend cells in either RPMI+5% FBS and store on ice until use, or proceed directly to cell staining for FACS.