



FEB 08, 2024

MSD V-PLEX Proinflammatory panel

patricia.garcia¹

¹The Francis Crick Institute



patricia.garcia

ABSTRACT

The V-PLEX panel offers analytically validated multiplex assay kits. Developed under rigorous design control, V-PLEX kits provide accurate and reproducible results with consistency from lot to lot. The V-PLEX Proinflammatory Panel 1 (human) Kit includes ten biomarkers that are associated with the inflammatory response and immune system regulation.

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.n92ldmd28l5b/v1

Protocol Citation: patricia.garcia
2024. MSD V-PLEX
Proinflammatory panel.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.n92ldmd28l5b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Feb 06, 2024

MSD V-PLEX Proinflammatory Panel

- 1** Prepare calibrator dilutions
 - a. Prepare the highest calibrator (Calibrator 1) by adding 1000mL of Diluent 2 to the lyophilized calibrator vial. After reconstituting, invert at least 3 times. Let it equilibrate for 15-30 minutes at room temperature and then vortex briefly using short pulses.
 - b. Prepare the next calibrator by transferring 100mL to the highest calibrator to 300mL of Diluent 2.
 - c. Use Diluent 2 as the zero calibrator.
- 2** Prepare Detection Antibody Solution

For one plate, combine 2.400mL of Diluent 3 with:

 - 60mL of SULFO-TAG Anti-hu IFN- γ Antibody
 - 60mL of SULFO-TAG Anti-hu IL-1b Antibody
 - 60mL of SULFO-TAG Anti-hu IL-2 Antibody
 - 60mL of SULFO-TAG Anti-hu IL-4 Antibody
 - 60mL of SULFO-TAG Anti-hu IL-6 Antibody
 - 60mL of SULFO-TAG Anti-hu IL-8 Antibody
 - 60mL of SULFO-TAG Anti-hu IL-10 Antibody
 - 60mL of SULFO-TAG Anti-hu IL-12p70 Antibody
 - 60mL of SULFO-TAG Anti-hu IL-13 Antibody
 - 60mL of SULFO-TAG Anti-hu TNF- α Antibody
- 3** Prepare Wash Buffer

For one plate combine:

 - 15mL of MSD Wash Buffer (20X)
 - 285mL of deionized water
- 4** Prepare Read Buffer T

For one plate combine:

 - 10mL of Read Buffer T (4X)
 - 10mL of deionized water
- 5** Assay Protocol

5.1 Wash and Add Sample

- a. Wash the plate 3 times with at least 150mL/well of Wash Buffer.
- b. Add 50uL of prepared samples, calibrators, or controls per well.
- c. Seal the plate with an adhesive plate seal.
- d. Incubate at room temperature with shaking for 2h.

5.2 Wash and Add Detection Antibody Solution

- a. Wash the plate 3 times with at least 150mL/well of Wash Buffer.
- b. Add 25uL of detection antibody solution to each well.
- c. Seal the plate with an adhesive plate seal.
- d. Incubate at room temperature with shaking for 2h.

5.3 Wash and Read

- a. Wash the plate 3 times with at least 150mL/well of Wash Buffer.
- b. Add 150uL of 2X Read Buffer T to each well.
- c. Seal the plate with an adhesive plate seal.
- d. Analyze the plate on an MSD instrument.