



Aug 18, 2022

Coating superfrost microscope slides with gelatin-chromium potassium sulfate

Benjamin Trist¹, Alejandra Rangel¹, Rain Kwan¹¹The University of Sydney, Sydney, Australia

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.n2bvj65owlk5/v1 Benjamin Trist

ABSTRACT

This protocol describes how to coat microscope slides with gelatin-chromium potassium sulfate (gelatin-chrom alum) in preparation for histology or immunohistochemical analysis of thin tissue sections. Slides coated with gelatin-chrom alum exhibit much better retention of tissues mounted after free-floating immunohistochemical staining, especially during alcohol and xylene dehydration steps immediately prior to mounting media embedding and cover slipping.

ATTACHMENTS

[it74bj7ap.docx](#)

DOI

dx.doi.org/10.17504/protocols.io.n2bvj65owlk5/v1

PROTOCOL CITATION

Benjamin Trist, Alejandra Rangel, Rain Kwan 2022. Coating superfrost microscope slides with gelatin-chromium potassium sulfate. **protocols.io** <https://protocols.io/view/coating-superfrost-microscope-slides-with-gelatin-ce7tthnn>



FUNDERS ACKNOWLEDGEMENT

Michael J Fox Foundation
Grant ID: ASAP-000497

KEYWORDS

Microscope slide coating, Histology, Immunohistochemistry, Tissue sections, Mounting

LICENSE

_____ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 12, 2022

LAST MODIFIED

Aug 18, 2022

OWNERSHIP HISTORY

Aug 12, 2022  madhavi.d

Aug 18, 2022  Benjamin Trist

PROTOCOL INTEGER ID

68563

MATERIALS TEXT



Equipment:

- Oven
- Heated magnetic stirrer
- Thermometer
- Chemical spatulas
- Slide racks
- Heat resistant beaker compatible with slide racks

Consumables:

- Magnetic stirrer bars
- Slide storage box
- Superfrost microscope slides

Key reagents:

-  [Gelatin Sigma](#) –
- [Aldrich Catalog #G2500](#)
-  [Chromium potassium sulfate Merck](#)
- [Millipore Catalog #101036](#)

Experimental Outline

2d 0h 0m 50s

Place heated magnetic stirrer in fume hood.

1

2 Pre-heat oven to \uparrow 42 °C .

3 Heat \square 1 L dH₂O to \uparrow 50 °C - \uparrow 60 °C and completely dissolve \square 10 g gelatin with aid of magnetic stirrer.

4



Add \square 1 g chromium potassium sulfate - solution should turn a pale green/blue and be completely clear.

5 Once dissolved reduce temperature to \uparrow 40 °C - \uparrow 50 °C .

6 Place slides for coating into designated slide coating rack(s).

7 Dip rack of slides into warm (\uparrow 40 °C - \uparrow 50 °C) gel mixture for approximately ⌚ 00:00:20^{50s} - ⌚ 00:00:30 .

8 Shake the excess liquid from the rack.

9 Repeat steps 7 and 8.

10 Place slide rack(s) to dry in the oven at \uparrow 42 °C for ⌚ 48:00:00 .

2d

11 Store slides in dust-free slide storage box.



12 Discard Gelatin-Chrom Alum solution into appropriate waste container.