



Oct 11, 2022

Purification of Plasmid DNA by Miniprep

Georgina Elisa Diego Macías¹, Ana Belem García González¹,
Irán Alessandra Chaparro Rodríguez¹, Jair Alexis Gardea Sáenz¹

¹Tecnologico de Monterrey Campus Chihuahua

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.n92ldp97ol5b/v1

G Georgina Elisa Diego Macías

ABSTRACT

Extraction of Plasmid DNA by Miniprep PureLink T M Quick Plasmid Miniprep Kits

DOI

dx.doi.org/10.17504/protocols.io.n92ldp97ol5b/v1

PROTOCOL CITATION

Georgina Elisa Diego Macías, Ana Belem García González, Irán Alessandra Chaparro Rodríguez, Jair Alexis Gardea Sáenz 2022. Purification of Plasmid DNA by Miniprep. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ldp97ol5b/v1>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 11, 2022

LAST MODIFIED












Oct 11, 2022



PROTOCOL INTEGER ID



71135

MATERIALS TEXT

- PureLink T M Quick Plasmid Miniprep Kits
- Resuspension Buffer (R3)
- RNase A
- Lysis Buffer (L7)
- Precipitation Buffer (N4)
- Wash Buffer (W9)
- Wash Buffer (W10)
- Buffer TE
- Centrifugation columns in recollection tubes
- Overnight LB-culture
- Eppendorf Tubes of 1.6 mL
- Etanol 96-100%

- 1 Centrifuge  **1-5 mL** of the overnight LB-culture. Remove all medium.
- 2 Add  **250 µL** Resuspension Buffer (R3) with RNase A to the cell pellet and resuspend the pellet until it is homogeneous
- 3 Add  **250 µL** Lysis Buffer (L7). Mix gently by inverting the capped tube until the mixture is ^{5m} homogeneous. Do not vortex. Incubate the tube at room temperature for  **00:05:00** .
- 4 Add  **350 µL** Precipitation Buffer (N4). Mix immediately by inverting the tube, or for large ^{10m} pellets, vigorously shaking the tube, until the mixture is homogeneous. Do not vortex. Centrifuge the lysate at $>12,000 \times g$ for  **00:10:00**
- 5 Load the supernatant from step 4 onto a spin column in a 2-mL wash tube. Centrifuge the ^{1m} column at $12,000 \times g$ for  **00:01:00** . Discard the flowthrough and place the column back into the wash tube.
- 6 Add  **500 µL** Wash Buffer (W10) with ethanol to the column. Incubate the column for ^{2m}  **00:01:00** at room temperature. Centrifuge the column at $12,000 \times g$ for  **00:01:00** . Discard the flowthrough and place column back into the wash tube
- 7 Add  **700 µL** Wash Buffer (W9) with ethanol to the column. Centrifuge the column at $12,000$ ^{2m}

× g for  00:01:00 . Discard the flowthrough and place the column into the wash tube.
Centrifuge the column at 12,000 × g for  00:01:00 . Discard the wash tube with the flowthrough

- 8 Place the Spin Column in a clean 1.5-mL elution tube. Add  75 µL of preheated TE Buffer^{1m} (TE) to the center of the column. Incubate the column for  00:01:00 at room temperature.
- 9 Centrifuge the column at 12,000 × g for 2 minutes. The elution tube contains the purified plasmid DNA. Discard the column. Store plasmid DNA at 4°C (short term) or store the DNA in aliquots at -20°C (long term).