




Jun 03, 2020

Coverslipping: Use of the Leica CV5030 Coverslipper

Allen Institute for Brain Science¹¹Allen Institute**1** Works for me dx.doi.org/10.17504/protocols.io.bg4xyxn

BICCN Allen Institute for Brain Science

 Dillan Brown

ABSTRACT

This protocol describes the coverslipping of glass slides with the Leica CV5030.

Note: Research reported in this publication was supported by the National Institute Of Mental Health of the National Institutes of Health under Award Number U19MH114830. The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

ATTACHMENTS

[EQ0003_Coverslipping-Use_of_Leica_CV5030_Coverslipper.docx](#) [Leica_CV5030_Instruction_Manual.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.bg4xyxn

PROTOCOL CITATION

Allen Institute for Brain Science 2020. Coverslipping: Use of the Leica CV5030 Coverslipper. **protocols.io**
dx.doi.org/10.17504/protocols.io.bg4xyxn

KEYWORDS

EQ0003, coverslipped, slides, slide

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 03, 2020

LAST MODIFIED

Jun 03, 2020

PROTOCOL INTEGER ID

37751