



Aug 23, 2022

Single-Cell Isolation of Human Knee Meniscus

Forked from [Single-Cell Isolation of Human Articular Cartilage](#)

hswahn¹, Martin Lotz¹

¹Scripps Research

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.n2bvj6k7blk5/v1

Human BioMolecular Atlas Program (HuBMAP) Method Development Community
Tech. support email: Jeff.spraggins@vanderbilt.edu

molmer

ABSTRACT

This is a protocol that describes the process of isolating single cells from human knee meniscus for scRNA-seq.

DOI

dx.doi.org/10.17504/protocols.io.n2bvj6k7blk5/v1

PROTOCOL CITATION

hswahn, Martin Lotz 2022. Single-Cell Isolation of Human Knee Meniscus. **protocols.io**
<https://protocols.io/view/single-cell-isolation-of-human-knee-meniscus-ce5ftg3n>

FORK NOTE

FORK FROM

Forked from [Single-Cell Isolation of Human Articular Cartilage](#), hswahn

KEYWORDS

single cell isolation, meniscus

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 10, 2022

LAST MODIFIED

Aug 23, 2022

PROTOCOL INTEGER ID

68487

1 ~1g of tissue from healthy donor knees (Grades 0-1) is collected from the meniscus. For details regarding

the tissue harvesting procedure please see
[dx.doi.org/10.17504/protocols.io.6qpvr614zvmk/v1](https://doi.org/10.17504/protocols.io.6qpvr614zvmk/v1)

- 2 Meniscal tissue is washed with **Room temperature**
Dulbeccos Phosphate Buffered Solution [DPBS] Gibco - Thermo
Fisher Catalog #14190-144
Bovine Calf Serum [CS] VWR international
supplemented with 10% **Ltd Catalog #10158-358**, 1%
Antibiotic-Antimycotic 9100x0 [Anti-Anti] Thermo Fisher
Scientific Catalog #15240062 and 1%
Penicillin-Streptomycin-Glutamine [PSG]
Corning Catalog #30-009-CI
- 3 The tissue is then finely minced with a #21 Feather disposable scalpel, and digested in 20mL 30m
DMEM with L-Glutamine 4.5g/L Glucose and Sodium Pyruvate Fisher
Scientific Catalog #MT-10-013-CV
supplemented with 1% Anti-Anti and 2%
Collagenase Type II Worthington Biochemical
Corporation Catalog #L5004177 using 100 rpm
shaking at **37 °C** for **00:30:00**.
- 4 Cells are gently passed through a 100 µm filter into a 50 mL centrifuge tube followed by gentle passage through a 40 µm filter into a fresh 50 mL centrifuge tube.
- 5 Filtered cells are spun down at 1200 rpm for **00:05:00** at **Room temperature** 5m
- 6 Carefully remove and collect the collagenase supernatant; the collagenase is reused in subsequent steps.
- 7 Cells are then resuspended in DMEM supplemented with 10% CS, 1% Anti-Anti and 1% PSG and stored at **37 °C**.
- 8 Repeat the digestion, spin down and filtration steps for **01:00:00** and then **02:00:00**, for a total of 6h 30m
03:30:00 of digestion.
- 9 Combine all collected cells and spin down at 1200 rpm for **00:05:00** at **Room temperature** 5m

- 10 The supernatant is carefully removed, and the remaining cell pellet is delicately resuspended in 10mL of **⚡ Room temperature** DPBS supplemented with 5% CS and 5 mM

[⌘ Ethylenediamine tetraacetic acid \[EDTA\] Fisher](#)

Scientific Catalog #BP120-500

- 11 Single cells are spun down at 1200 rpm for ⌚ **00:05:00** at **⚡ Room temperature** .

5m

- 12 The supernatant is carefully removed, and the remaining cell pellet is delicately resuspended in 10 mL of DPBS supplemented with 0.04%

[⌘ Bovine Serum Albumin \[BSA\] Fisher](#)

Scientific Catalog #9048-46-8

- 13 The Invitrogen Countess II FL automated cell counter is used to quantify single cells and determine cell viability. Live cells are determined by trypan blue staining. If >70% cell viability is confirmed, the single cell suspension is diluted to a concentration of 1×10^6 cells/mL for single cell RNA-seq library preparation.