



Sep 09, 2020

TPGY BROTH (TRYPTONE PEPTONE GLUCOSE YEAST EXTRACT BROTH)

Roey Angel¹, Ana Lara-Rodriguez¹, Eva Petrova¹¹Soil and Water Research Infrastructure**1** *Works for me* This protocol is published without a DOI.

SoWa RI Anaerobic and Molecular Microbiology (public)
Tech. support email: eva.petrova@bc.cas.cz

Ana Lara

ABSTRACT

For the enrichment of *Clostridium* sp. from soil samples, sufficient for 100 cultivation tubes.

EXTERNAL LINK

<http://www.himedialabs.com/TD/M969.pdf>

PROTOCOL CITATION

Roey Angel, Ana Lara-Rodriguez, Eva Petrova 2020. TPGY BROTH (TRYPTONE PEPTONE GLUCOSE YEAST EXTRACT BROTH). **protocols.io**
<https://protocols.io/view/tpgy-broth-tryptone-peptone-glucose-yeast-extract-3e5gg6>



EXTERNAL LINK

<http://www.himedialabs.com/TD/M969.pdf>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

May 28, 2019

LAST MODIFIED

Sep 09, 2020

PROTOCOL INTEGER ID

23741

GUIDELINES

taken from: <http://www.himedialabs.com/TD/M969.pdf>

MATERIALS

NAME	CATALOG #	VENDOR
Peptone A FROM MEAT	G213.SIZE.500g	Bio Basic Inc.
Sodium thioglycolate	S0265.SIZE.25g	Bio Basic Inc.
BD Bacto™ Yeast Extract	212750	BD Biosciences
Enzymatic Digest Casein Hydrolysate	AAJ12855Q1	
Dextrose	D9434	Sigma – Aldrich

SAFETY WARNINGS

Use heat reasistant gloves to take media out of the autoclave and to put and take out the media from the water bath.

BEFORE STARTING

1. Make sure that the glass tubes are sterile.
2. Pour media only in the flow box.

Medium preparation

1

1. Dilute in 1000 ml distilled water:

 **50 g Casein Hydrolysate**

 **5 g Peptic Digest from Animal Tissue (Peptone A)**

 **20 g Yeast Extract**


 **4 g Dextrose**

 **1 g Sodium Thioglycolate**

Final pH (at 25°C) 7.0±0.2




Heat if necessary, to dissolve the medium completely.





- 2 Sterilize by autoclaving at 15 lbs pressure  **121 °C** for  **00:15:00** .

- 3 Refrigerate  **4 °C** the sterile medium until use.

Before inoculation of media

4

Pour  **10 mL** of media in sterile tubes.

- 5 Heat media in water bath  **100 °C**  **00:10:00** **OR** alternative do a gass exchange using the gassing station for 3 minutes to each tube **OR** put tubes in to the anaerobic tent for  **24:00:00** to  **48:00:00** before inoculating.