

MAR 20, 2024

OPEN ACCESS



Protocol Citation: Callen Hyland, Rob Steele 2024. Blackwater medium. protocols.io https://protocols.io/view/blackwat er-medium-daj72crn

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development We are still developing and optimizing this protocol

Created: Mar 13, 2024

Last Modified: Mar 20, 2024

PROTOCOL integer ID: 96607

Keywords: Hydra, invertebrates,

animal care

Blackwater medium

Callen Hyland¹, Rob Steele²

¹University of San Diego; ²UC Irvine



Callen Hyland University of San Diego

ABSTRACT

Blackwater medium is a laboratory growth medium for Hydra that replicates the chemical composition of blackwater swamps. This medium has lower low ionic strength than traditional Hydra medium and contains humic acid. Blackwater medium is ideal for Hydra species isolated from blackwater swamps and improves the health of some transgenic Hydra lines that are prone to infection.

IMAGE ATTRIBUTION

Callen Hyland

GUIDELINES

Please see Low cost methods for Hydra care for general guidelines for working with Hydra.

MATERIALS

- Sodium Chloride Fisher Scientific Catalog #S271
- Humic Acid Sodium Salt Merck MilliporeSigma (Sigma-Aldrich) Catalog #H16752
- Lab scoops
- Weigh boats or weigh papers
- Electronic balance
- DI water
- 50 mL conical tubes
- Syringe filter or vacuum filter
- Adjustable volume micropipette and tips
- Liter bottles or plastic carboy

Stock solutions

1 Make stock solutions in 50 mL conical tubes with DI water.

A	В	С
Component	Stock conc.	Weight for 50 mL
Calcium chloride dihydrate (CaCl2*2H2O)	2.5 M	18.4 g
Magnesium chloride hexahydrate (MgCl2*6H2O)	1.0 M	10.2 g
Sodium chloride (NaCl)	5.0 M	14.6 g
Potassium chloride (KCI)	1.0 M	3.73 g
Humic acid sodium salt	20 g/L	1 g

Table 1. Stock solutions for blackwater medium.

2 Not all of the humic acid sodium salt will go into solution. Filter the humic acid solution with a syringe filter or vacuum filter to remove undissolved particles.

Working solution

3 Prepare blackwater medium by diluting stock solutions in DI water. Use a 1 L bottle, or scale up volumes for a large plastic carboy. Mix well before using.

А	В	С
Component	Final conc.	Volume of stock for 1 L
Calcium chloride dihydrate (CaCl2*2H2O)	0.015 mM	6.0 µL
Magnesium chloride hexahydrate (MgCl2*6H2O)	0.0185 mM	18.4 µL
Sodium chloride (NaCl)	0.125 mM	25 μL
Potassium chloride (KCI)	0.0075 mM	75 μL
Humic acid sodium salt	0.01 g/L	0.5 mL

Figure 2. Chemical composition of Blackwater medium.